medi® Soft OA light

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Knies-Softorthese zur Entlastung

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medi® Soft OA light

Donning instructions

Soft knee orthosis for pressure relief

Note: It is best if you sit on the edge of a chair or stool when putting on the brace. After every donning step check that the brace is still fitted correctly to your leg and has not become twisted.

1 Bend the affected leg to an approx. 70 degree angle and slip the brace over it (image 1).

2 Pull the brace upwards until the fabric touches the lower edge of the kneecap, and the kneecap is positioned in the centre of the opening. Make sure that the brace is not twisted when fitted to the leg (image 2).

3 Now close the upper hook and loop closures (image 3).

4 Close straps 1 and 2. Adjust the straps to the values agreed with your orthopaedic technician.

Note: The initial steps are intended as a basis for orientation. According to your pain levels, you can individually adjust the strap and raise or reduce the pressure (image 4).

5 Adjust the straps in such a way that you feel a pleasant sensation of pressure and pain relief. Walk a few steps and make sure that the brace is fitted correctly to the leg (image 5).

medi tip

Use the QR code or the short URL medi.biz/soft-oa-light-aza to access the video with step-by-step instructions.