medi GmbH & Co. KG
Medicusstraße 1
D-95448 Bayreuth
Germany
T +49 921 912-0
F +49 921 912-780
ortho@medi.de
www.medi.de

medi Australia Pty Ltd
Unit 4/13-21 Hallmark Street
Pendle Hill NSW 2145
Australia
T +61 2 9890 8696
F +61 2 9890 8439
sales@medi.au
www.mediaustralia.com.au

medi Austria GmbH
Adamgasse 16/7
6020 Innsbruck
Austria
T +43 512 57 95 15
F +43 512 57 95 15 45
vertrieb@medi-austria.at
www.medi-austria.at

medi Belgium NV
Staatsbaan 77/0099
3945 Ham
Belgium
T: + 32-11 24 25 60
F: +32-11 24 25 64
info@medibelgium.be
www.medibelgium.be

medi Canada Inc / médi Canada Inc
104-1375 Lionel-Boulet,
Varennes, Québec,
QC Canada J3X 1P7
T +1 450-583-3317 / +1 800 361 3153
F +1 888-583-6827
service@medicanada.ca
www.medicanada.ca

MAXIS a.s.,
medi group company
Slezská 2127/13
120 00 Prague 2
Czech Republic
T: +420 571 633 510
F: +420 571 616 271
info@maxis-medica.com
www.maxis-medica.com

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Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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M.4[®] AGR

Intended purpose

M.4 AGR is a brace with a frame that is intended to reinforce and stabilise the knee joint using limitation of extension and flexion.

Indications

All indications that require physiological stabilisation and reinforcement of the knee joint and/or stress relief of the joint in at least two planes, such as:

- For genu recurvatum

Contraindications

None have been noted to date.

Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in conjunction with sweating) or due to the composition of the material. The immobilisation which is achievable by the device presents a risk for development of thrombosis in the lower veins of the legs, especially in the context of surgical interventions or trauma.

Intended users and patient target groups


Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Important Note

Only change or modify the extension limits on the instruction of your healthcare provider or orthopaedic technician. To avoid hyperextension, always insert extension inserts at 0°. Furthermore, ensure the correct positioning of the brace: the central parts of the hinge should be in line with the top of the kneecap; the hinge should be positioned behind the medial lateral line of the leg.

Instructions for the orthopaedic technician

Adjusting the extension:

Extension limitation inserts (grey; extension limits):
0°, 10°, 20°, 30°.

- The medial adjustment (inner) and lateral adjustment (outer) of the degrees should always be equal.

Extension inserts at 0° are already integrated into the brace.

- To change the inserts, loosen the cover screws, remove the cover and unscrew the insert screws. To easily remove the inserts, gently bend the knee device.
- The extension inserts are marked R-E (RIGHT) or L-E (LEFT) and indicate the respective degrees. Position the insert in such a way that the letter L or R is the same as the hinge plate. Push the new extension insert between the two plates and fasten it with the screws.
- Bend and flex the hinge to check that the inserts have been correctly positioned. Lastly, reposition the hinge cover.

Replacing the rubber washers

The rubber washer inside the hinge can be replaced if it becomes worn or in order to increase or reduce the resistance by selecting the correct density of the mechanical device (increase or decrease, respectively, the hardness of the rubber washers). Before replacement, select the rubber washer with the most suitable hardness (shore - hardness of the rubber washer: 80 shore = min, 90 shore = medium, 95 shore = max), cross-referencing the weight and recurvatum degree as shown in Table 1.

NB: The rubber washer positioned in the medial part (inner) and the rubber washer located in the lateral part (outer) should always have the same hardness and degree of wear.

Degree of recurvatum



Table 1

Instructions for use

- Remove the Velcro tab located on the internal part of the hinge.
- Unscrew the six screws and then remove the plate.
- Carefully extract the rubber washer located in the anterior part of the hinge and insert the new rubber washer as shown in the photograph.
- Reposition the hinge plate making sure to correctly fitted and fasten it using the four 4x6 screws (the shorter screws) for the upper part and the two 4x8 screws for the lower part.
- Lastly, affix the Velcro oval (provided in the packaging) to the cover.

Fitting instructions

- (fig. 1) The tibia support is positioned with the wider part towards the top and the concave part cushioning the inside. Remove the paper and affix the

two adhesive hook and loop tabs to the outer part of the tibia support so that they correspond with the two tibial straps (#1 and #4).

- (fig. 2) Open all of the hook and loop closures. Sit on a chair and position the brace around the affected leg; flex the leg so that it forms a 45° angle. Position the hinges in such a way that the centres of the hinge are positioned in line with the top part of the kneecap.
- (fig. 3+4) Keep the knee in a flexed position and support the two hinges at the back of the leg. Fasten the upper rear strap (#2) and then the upper front strap (#6).
- (fig. 5+6) Keep the knee in a flexed position and support the two hinges at the back of the leg. Fasten the upper rear strap (#2) and then the upper front strap (#6).
- (fig. 7+8) Fasten the lower rear strap (#5); and then the upper rear strap (#3).
- For knee devices worn without a tibia support, follow the numbering of the straps: 1, 2, 3, 4, 5 and 6.
To prevent the brace from slipping, gently flex the knee and check that the Velcro straps have been correctly fastened.

Care instructions

Hook and loop fasteners should be closed for washing. Soap residues, lotions and ointments can cause skin irritation and material wear.

- Wash the paddings and covers by hand, preferably using medi clean washing agent.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

Check regularly that the joint screws are secure, tighten them if necessary.

Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

Material composition

Avional, velours, polyester, PU-foam, Elastomer

Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

Disposal

The product can be disposed of in the domestic waste.

Your medi Team

Wishes you a quick return to full fitness.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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