medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth
Germany
T +49 921 912-0
F +49 921 912-780
ortho@medi.de
www.medi.de

medi Australia Pty Ltd
83 Fennell Street
North Parramatta NSW 2151
Australia
T +61-2 9890 8696
F +61-2 9890 8439
sales@mediaustralia.com.au
www.mediaustralia.com.au

medi Austria GmbH
Adamgasse 16/7
6020 Innsbruck
Austria
T +43 512 57 95 15
F +43 512 57 95 15 45
vertrieb@medi-austria.at
www.medi-austria.at

medi Belgium NV
Staatsbaan 77/0099
3945 Ham
Belgium
T: + 32-11 24 25 60
F: +32-11 24 25 64
info@medibelgium.be
www.medibelgium.be

medi Canada Inc / médi Canada Inc
104-1375 Lionel-Boulet,
Varenes, Québec,
QC Canada J3X 1P7
T +1 450-583-3317 / +1 800 361 3153
F +1 888-583-6827
service@medicanada.ca
www.medicanada.ca

MAXIS a.s.,
medi group company
Slezská 2127/13
120 00 Prague 2
Czech Republic
T: +420 571 633 510
F: +420 571 616 271
info@maxis-medica.com
www.maxis-medica.com

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
E006500/10.2024

medi 3 C

Orthese zur Entlastung der LWS / BWS in Sagittalebene ·
Brace for relieving load on the lumbar spine/thoracic spine
in the sagittal plane

Gebrauchsanweisung. Instructions for use. Mode d'emploi. Instrucciones de uso.
 Instruções para aplicação. Istruzioni per l'uso. Gebruiksaanwijzing.
 Brugsvejledning. Bruksanvisning. Návod k použití. Uputa za uporabu.
 Инструкция по использованию. Kullanna kilavuzu. Instrukcja zakładowania.
 Navodila za uporabo. Інструкція з використання. دليل الاستخدام. Navodila za
 uporabo. Instrucțiunile de utilizare. שימוש הוראות.

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Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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medi 3 C

Intended purpose

The medi 3C is a brace for relieving load  on the lumbar spine/thoracic spine (sagittal limitation of range of motion)

Indications

All indications requiring a partial relief by correction of the lumbar/thoracic spine and movement restriction in sagittal plane, such as:

- Stable vertebral compression fractures of the middle and lower thoracic and lumbar spine
- Hyperkyphosis in Scheuermann's disease
- Temporary postoperative orthotic care

Contraindications

Unstable fractures of the lumbar/thoracic spine with or without neurologic deficit.

Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in conjunction with sweating) or due to the composition of the material.

Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Fitting instructions

- Select the correct basic model according to the waist circumference. To determine the length of the sternal stick, measure the distance between the top of the pubic crest and the lower edge of the sternoclavicular joints.
- Loosen the screws in the base plate to adjust the length of the brace. Fully extend the sternal pad and adjust the length between the upper and lower edges. Make sure that the waistband lies above the iliac crest.
- Now adjust the symphysis and sternal rods. The reclination effect can be set to suit the individual patient.
- Secure the dorsal pad in place and fit the waistband.
- The sternal rod and the lateral extensions of the base plate can be shaped by hand to the patient's body by bending slowly over your knee or the back of a chair. Round bending irons may be used. Do not bend the sternal rod in the radius of movement of the sternal pad. In the interest of safety remove the sternal pad from the sternal rod before bending.
- Slide the sleeve nut into the opening to

adjust the sternal pad and fasten with the screws. We recommend securing the screws in place with Loctite.

- Check the position of the pads: in the hyperextended position the upper edge of the sternal pad is level with the lower edge of the sternoclavicular joints. The sternal contours can be shaped by hand. The lower edge of the symphysis pad is level with the pelvic crest. Bend the lower part of the pad outwards.
- Please verify that all screws have been sufficiently tightened.
- We recommend checking the function of the medi 3 C hyperextension orthosis the day after fitting and adjusting it if necessary.

Closing the ratchet fastener

Attention: do not open the fastener covers beyond 180°.

Fold the covers of the two fasteners back and push the two ends of the belt forwards through the openings in the fasteners until the tips of the belt are visible at the front. The belt is tightened by opening and closing the fastener covers, and the patient can easily bring himself into the hyperextended position until the correct degree of hyperextension is reached. Both covers must then be closed. The ends of the belt can be trimmed to the desired length with scissors. When removing the brace it is helpful to open one of the fastener covers completely with the thumb, then the belt can be pulled out.

Caution: Do not close the fastener by force or beyond the pain trigger point. Normal breathing should be possible at all times!

Care instructions

The brace can be worn in both fresh and

salt water. For this purpose, please remove the cotton covers of the pads. Soap residues, lotions and ointments can cause skin irritation and material wear. Rinse the brace with clean water afterwards.

- Wash the paddings and covers by hand, preferably using medi clean washing agent.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

Material composition

Aluminum, polyethylene (foam), cotton (covers)

Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

Disposal

The product can be disposed of in the domestic waste.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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