**medi GmbH & Co. KG**  
**Medicusstraße 1**  
**95448 Bayreuth**  
**Germany**  
**T +49 921 912-0**  
**F +49 921 912-780**  
**ortho@medi.de**  
**www.medi.de**

medi Australia Pty Ltd  
 83 Fennell Street  
 North Parramatta NSW 2151  
 Australia  
 T +61-2 9890 8696  
 F +61-2 9890 8439  
 sales@mediaustralia.com.au  
 www.mediaustralia.com.au

medi Austria GmbH  
 Adamgasse 16/7  
 6020 Innsbruck  
 Austria  
 T +43 512 57 95 15  
 F +43 512 57 95 15 45  
 vertrieb@medi-austria.at  
 www.medi-austria.at

medi Belgium NV  
 Staatsbaan 77/0099  
 3945 Ham  
 Belgium  
 T: + 32-11 24 25 60  
 F: +32-11 24 25 64  
 info@medibelgium.be  
 www.medibelgium.be

medi Canada Inc / médi Canada Inc  
 104-1375 Lionel-Boulet,  
 Varennes, Québec,  
 QC Canada J3X 1P7  
 T +1 450-583-3317 / +1 800 361 3153  
 F +1 888-583-6827  
 service@medicanada.ca  
 www.medicanada.ca

MAXIS a.s.,  
 medi group company  
 Slezská 2127/13  
 120 00 Prague 2  
 Czech Republic  
 T: +420 571 633 510  
 F: +420 571 616 271  
 info@maxis-medica.com  
 www.maxis-medica.com

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
medi Trading (Shanghai) Co. Ltd.  
 Room 1703, T1, CES WEST BUND CENTER,  
 No. 277, Long Lan Road  
 Xuhui District  
 200232 Shanghai  
 People's Republic of China  
 T: +86-21 50582319  
 F: +86-21 50582319

UK Responsible Person  
 medi UK Ltd.  
 Plough Lane  
 Hereford HR4 0EL  
 Great Britain  
 T +44-1432 37 35 00  
 F +44-1432 37 35 10  
 enquiries@mediuk.co.uk  
 www.mediuk.co.uk

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#### Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If undue pain or an unpleasant sensation develops while you are wearing the product, please consult your doctor or orthotist immediately. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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# medi 4C flex

## Intended purpose

The medi 4C flex is a brace for relieving load on the lumbar spine/thoracic spine (sagittal and frontal limitation of range of motion)

## Indications

All indications that require partial load relief via correction of the lumbar/thoracic vertebrae and sagittal and frontal limitation of the range of motion, such as:

- Stable vertebral body compression fractures of the middle and lower lumbar and thoracic spine
- Temporary postoperative immobilisation
- Conservative therapy for inoperable tumours and metastases of the spine, or as a post-operative aid should they be operable
- Axial stabilisation for osteoporotic sintering of the thoracic and lumbar spine

## Contraindications

Unstable fractures of the lumbar/thoracic vertebrae with or without neurologic deficit.

## Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related

to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in conjunction with sweating) or due to the composition of the material.

## Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

## Fitting instructions

Please wear the product on thin outerwear.

### • Adjusting the width

First loosen the screws on the chassis and adjust the orthosis to the width of the patient's body. Please always check that the functioning of the hinged fastener is not restricted.

Now loosen the two screws on the side of the chassis and the screws on the fastener. Adjust the orthosis to suit the patient's length – distance from the symphysis pubis to the subclavicular contact area. Then retighten all the screws.

### • Adjusting the strap length

Using the hook and loop fastener and the reversing loop adjust the length of the strap so that the pad comes to lie in the middle of the back. Remember that a few centimetres of the strap will have to be pulled through the fastener.

- **Opening and closing the fastener**

Slide the grey bar sideways to open the fastener. To close the orthosis flip the cover over until it audibly clicks into place.

---

#### Care instructions

The orthosis can be cleaned with a damp cloth, the back pad with a brush. Soap residues, lotions and ointments can cause skin irritation and material wear.

- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

---

#### Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

---

#### Material composition

Aluminium, PVC

#### Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

#### Disposal

The product can be disposed of in the domestic waste.

#### Your medi Team

Wishes you a quick return to full fitness.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities

in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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