medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth
Germany
T +49 921 912-0
F +49 921 912-780
ortho@medi.de
www.medi.de

medi Australia Pty Ltd
 83 Fennell Street
 North Parramatta NSW 2151
 Australia
 T +61-2 9890 8696
 F +61-2 9890 8439
 sales@mediaustralia.com.au
 www.mediaustralia.com.au

medi Austria GmbH
 Adamgasse 16/7
 6020 Innsbruck
 Austria
 T +43 512 57 95 15
 F +43 512 57 95 15 45
 vertrieb@medi-austria.at
 www.medi-austria.at

medi Belgium NV
 Staatsbaan 77/0099
 3945 Ham
 Belgium
 T: +32-11 24 25 60
 F: +32-11 24 25 64
 info@medibelgium.be
 www.medibelgium.be

medi Canada Inc / médi Canada Inc
 104-1375 Lionel-Boulet,
 Varennes, Québec,
 QC Canada J3X 1P7
 T +1 450-583-3317 / +1 800 361 3153
 F +1 888-583-6827
 service@medicanada.ca
 www.medicana.ca

MAXIS a.s.,
 medi group company
 Slezská 2127/13
 120 00 Prague 2
 Czech Republic
 T: +420 571 633 510
 F: +420 571 616 271
 info@maxis-medica.com
 www.maxis-medica.com

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medi Trading (Shanghai) Co. Ltd.
 Room 1703, T1, CES WEST BUND CENTER,
 No. 277, Long Lan Road
 Xuhui District
 200232 Shanghai
 People's Republic of China
 T: +86-21 50582319
 F: +86-21 50582319

UK Responsible Person
 medi UK Ltd.
 Plough Lane
 Hereford HR4 0EL
 Great Britain
 T +44-1432 37 35 00
 F +44-1432 37 35 10
 enquiries@mediuk.co.uk
 www.mediuk.co.uk

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Important notes

This medical device **MD** is made only for multiple uses on a single patient . If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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medi SAS[®] 45

Intended purpose

The medi SAS 45 is a shoulder joint support for immobilisation  and stress relief .

Indications

For all indications which require immobilisation and stress relief of the shoulder joint in abduction, such as:

- After rotator cuff repair
- Humeral head fractures (nonsurgical, postoperative)
- After insertion of shoulder prosthesis
- After stabilization of acromioclavicular joint dislocations
- After biceps tenodesis
- Impingement syndrome (nonsurgical, postoperative)

Contraindications

All indications which expressly require internal rotation or abduction of more than 45° or less than 30°.

Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical

irritation of the skin (particularly in conjunction with sweating) or due to the composition of the material.

Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Important notes

Only change the angle of abduction on the instructions of your doctor.

Fitting instructions

the medi SAS 45 is supplied in a single universal size and has two options for use, providing either 45° or 30° of abduction.

- Loosen all the straps including the hand section.
- Place the cushion on the hip of the affected side and lay the shoulder strap over the opposite shoulder. Fasten the straps to the cushion with the strip, first the back and then the front. Attach one end of the hip belt to the front of the cushion.
- Now secure the hand rest to the top of the cushion with the hook and loop fastener.
- Attach the two forearm straps to the hand rest, and secure the upper arm strap to the cushion slightly obliquely. The upper arm strap prevents the arm slipping backwards.

- Slip your arm into the sling and fasten all the strips. Make sure that the straps are not fastened too tightly. Simply rotate the cushion to select 45° or 30° abduction.

Note: Make sure that the head of the humerus is not pushed upwards when the upper arm is resting on the cushion.

Care instructions

The product should not be worn in water. Please close hook and loop fasteners prior to washing. Traces of soap, creams or ointments can cause skin irritations and material wear.

- Pillow case, bag and straps can be washed by hand, preferably using mediclean washing agent, or in the washing machine with the delicate cycle at 30°C using mild detergent without fabric softener.
- Do not bleach.
- Air dry.
- Do not iron.
- Do not dry clean.

Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

Material composition

PU-foam, Polyamid, aluminum

Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

Disposal

The product can be disposed of in the domestic waste.

Your medi team

wishes you a quick return to full fitness.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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