**medi GmbH & Co. KG**  
**Medicusstraße 1**  
**D-95448 Bayreuth**  
**Germany**  
**T +49 921 912-0**  
**F +49 921 912-780**  
**ortho@medi.de**  
**www.medi.de**

medi Australia Pty Ltd  
 83 Fennell Street  
 North Parramatta NSW 2151  
 Australia  
 T +61-2 9890 8696  
 F +61-2 9890 8439  
 sales@mediaustralia.com.au  
 www.mediaustralia.com.au

medi Austria GmbH  
 Adamgasse 16/7  
 6020 Innsbruck  
 Austria  
 T +43 512 57 95 15  
 F +43 512 57 95 15 45  
 vertrieb@medi-austria.at  
 www.medi-austria.at

medi Belgium NV  
 Staatsbaan 77/0099  
 3945 Ham  
 Belgium  
 T: +32-11 24 25 60  
 F: +32-11 24 25 64  
 info@medibelgium.be  
 www.medibelgium.be

medi Canada Inc / médi Canada Inc  
 104-1375 Lionel-Boulet,  
 Varennes, Québec,  
 QC Canada J3X 1P7  
 T +1 450-583-3317 / +1 800 361 3153  
 F +1 888-583-6827  
 service@medicanada.ca  
 www.medicanada.ca

MAXIS a.s.,  
 medi group company  
 Slezská 2127/13  
 120 00 Prague 2  
 Czech Republic  
 T: +420 571 633 510  
 F: +420 571 616 271  
 info@maxis-medica.com  
 www.maxis-medica.com

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medi Trading (Shanghai) Co. Ltd.  
 Room 1703, T1, CES WEST BUND CENTER,  
 No. 277, Long Lan Road  
 Xuhui District  
 200232 Shanghai  
 People's Republic of China  
 T: +86-21 50582319  
 F: +86-21 50582319

UK Responsible Person  
 medi UK Ltd.  
 Plough Lane  
 Hereford HR4 0EL  
 Great Britain  
 T +44-1432 37 35 00  
 F +44-1432 37 35 10  
 enquiries@mediuk.co.uk  
 www.mediuk.co.uk

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#### Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If undue pain or an unpleasant sensation develops while you are wearing the product, please consult your doctor or orthotist immediately. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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# medi ROM Walker

## Intended purpose

medi ROM Walker is a lower-leg foot brace for mobilisation with an adjustable range of motion.

## Indications

All indications which require early functional mobilisation with limitation of the range of motion in joints of the lower-leg foot region, such as:

- Following fibula fractures (post-operative and conservative)
- Following fracture-dislocations of the upper ankle
- Following injuries of the forefoot and midfoot and the tarsus (conservative)
- After operations on ligament, soft tissue and tendon injuries
- Secondary treatment of talus fractures, calcaneal fractures and medial malleolus fractures
- For early functional therapy of Achilles tendon rupture (conservative/post-operative)

## Contraindications

Unstable fractures or fractures of the proximal tibia or fibula

## Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)

- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in conjunction with sweating) or due to the composition of the material.

## Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

## Adjusting plantar- and dorsiflexion

To adjust the range of movement turn the opening in the plastic cover to the current position of the pins. Press the pin in and slide to the desired degree figure.

Plantarflexion: 0°, 10°, 20°, 30°, 40°

Dorsiflexion: 20°, 10°, 0°

The medial (inside) and the lateral (outside) wedge settings must always be the same.

## Information for technicians

Where the joint is fixed at an angle other than 0° (e.g. immobilised at 30° plantarflexion), the patient will require additional walking supports. If the patient is intended to be able to put weight on the foot a plantarflexion

can also be achieved by elevating the heel using the heel wedge set which is available separately, whereby the joint is then fixed at 0°.

### Fitting instructions

- First open all the hook and loop fasteners and the leg/foot padding (A). Place your foot in the padding so that the heel stands firmly in the pre-formed heel bed and close the padding (B). Patients with very slim calves can trim the pad to the right length using scissors if there is a large overlap of the padding.
- Place your foot in the rigid shell of the Walker. The two guide supports should lie exactly along the middle of the inside and outside of the calf. Then close both the hook and loop fasteners of the foot section and the back heel strap (C).
- Finally fasten the three leg straps from the bottom up (D). Make sure that the straps lie flat against the pad and are tightly closed at the front. However, be careful not to impair the blood circulation.

For additional stabilisation of the lower leg the enclosed plastic shells can be fixed to the pad. The broader one is intended for the calf area. The padding can be used for additional cushioning (e.g. across the instep, at the heel strap, between the ankle bone and joint).

### Important notes

Only change the extension and flexion limits on instruction from your doctor or orthotist to do so.

### Care instructions

Hook and loop fasteners should be closed for washing. Soap residues, lotions and ointments can cause skin irritation and material wear.

- Wash the paddings by hand, preferably using medi clean washing agent.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

### Storage instructions

Keep the brace in a dry place and do not expose to direct sunlight.

### Material Composition

Aluminum, PU, cotton

### Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

### Disposal

The product can be disposed of in the domestic waste.

### Your medi team

Wishes you a speedy recovery!

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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