**medi GmbH & Co. KG**  
**Medicusstraße 1**  
**D-95448 Bayreuth**  
**Germany**  
**T +49 921 912-0**  
**F +49 921 912-783**  
**ortho@medi.de**  
**www.medi.de**

medi Australia Pty Ltd  
 83 Fennell Street  
 North Parramatta NSW 2151  
 Australia  
 T +61-2 9890 8696  
 F +61-2 9890 8439  
 sales@mediaustralia.com.au  
 www.mediaustralia.com.au

medi Austria GmbH  
 Adamgasse 16/7  
 6020 Innsbruck  
 Austria  
 T +43 512 57 95 15  
 F +43 512 57 95 15 45  
 vertrieb@medi-austria.at  
 www.medi-austria.at

medi Belgium NV  
 Staatsbaan 77/0099  
 3945 Ham  
 Belgium  
 T: + 32-11 24 25 60  
 F: +32-11 24 25 64  
 info@medibelgium.be  
 www.medibelgium.be

medi Canada Inc / médi Canada Inc  
 104-1375 Lionel-Boulet,  
 Varennes, Québec,  
 QC Canada J3X 1P7  
 T +1 450-583-3317 / +1 800 361  
 3153  
 F +1 888-583-6827  
 service@medicanada.ca  
 www.medicanada.ca

MAXIS a.s.,  
 medi group company  
 Slezská 2127/13  
 120 00 Prague 2  
 Czech Republic  
 T: +420 571 633 510  
 F: +420 571 616 271  
 info@maxis-medica.com  
 www.maxis-medica.com

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medi Trading (Shanghai) Co. Ltd.  
 Room 1703, T1, CES WEST BUND CENTER,  
 No. 277, Long Lan Road  
 Xuhui District  
 200232 Shanghai  
 People's Republic of China  
 T: +86-21 50582319  
 F: +86-21 50582319

UK Responsible Person  
 medi UK Ltd.  
 Plough Lane  
 Hereford HR4 0EL  
 Great Britain  
 T +44-1432 37 35 00  
 F +44-1432 37 35 10  
 enquiries@mediuk.co.uk  
 www.mediuk.co.uk

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#### Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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# medi SAS<sup>®</sup> multi

## Intended Purpose

medi SAS multi is a shoulder joint support for immobilisation and stress relief.

## Indications

All indications requiring immobilization of the shoulder joint in abduction and optional in external rotation such as:

- After rotator cuff repair
- After shoulder stabilization
- After anterior shoulder dislocation (nonsurgical)
- Humeral head fractures (nonsurgical, postoperative)
- After insertion of shoulder prosthesis
- After stabilization of acromioclavicular joint dislocations
- After biceps tenodesis
- Impingement syndrome (nonsurgical, postoperative)

## Contraindications

All indications that specifically require an internal rotation or an abduction of more than 15°.

## Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the

area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in conjunction with sweating) or due to the composition of the material.

## Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

## Fitting instructions

The orthosis is normally fitted for the first time in the operating theatre. It should be put on and taken off on a table wherever possible. Please pay attention to the following tips:

The snap fasteners on the front of the cushion should be opened. The shoulder strap already lies on the table.

1. Place the arm that has been operated on in the open pocket (Fig. 1).
2. Now close the pocket (Fig. 2).

Also use the enclosed Velcro straps to secure it more strongly (Fig. 3).

3. In the case of the given indication, also connect the outside rotary wedge (Fig. 4). Beforehand, however, it is necessary under all circumstances to ensure that the indication also enables a positioning in external rotation.

4. Pass the waist strap around your back and close the snap fastener. (Fig.5)
5. Now pass the shoulder strap over your uninjured shoulder to the front. Grasp the snap fastener end with the hand on the injured side (Fig. 6). Now climb into the opening between the shoulder and underarm straps (Fig. 7).
6. Finally close the snap fastener. (Fig. 8). The shoulder strap can be put on more easily if you have a home help. Have the person remove the shoulder strap completely (open snap fasteners). They can now help you into the shoulder strap as if you were putting on a jacket. (Fig. 7). Then proceed in the manner described above. The shoulder strap is then simply clicked into the cushion. If you have the feeling during the night that the arm or shoulder is not securely positioned, you can put a pillow under the brace or you can tighten the strap around the waist.

In principle, this product must be used under medical supervision. When it is applied for the first time, this must be under the instructions of trained and specialised personnel, such as an orthopaedic technician. If the support is individually fitted by a qualified person, this will ensure that it is effective and comfortable.

If it is to be used with other products, this must be discussed with the prescribing doctor.

#### Care instructions

The product should not be worn in water. Soap residues, lotions and ointments can cause skin irritation and material wear.

- Pillow case, bag and straps can be washed by hand, preferably using medi clean washing agent, or in the washing machine with the delicate

cycle at 30°C using mild detergent without fabric softener.

- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

#### Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

#### Material composition

PU foam, polyamide, aluminum

#### Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

#### Disposal

The product can be disposed of in the domestic waste.

#### Your medi team

wishes you a quick return to full fitness.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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