**medi GmbH & Co. KG**  
**Medicusstraße 1**  
**95448 Bayreuth**  
**Germany**  
**T +49 921 912-0**  
**F +49 921 912-780**  
**ortho@medi.de**  
**www.medi.de**

medi Australia Pty Ltd  
 Unit 4/13-21 Hallmark Street  
 Pendle Hill NSW 2145  
 Australia  
 T +61 2 9890 8696  
 F +61 2 9890 8439  
 sales@medi.au  
 www.medi.australia.com.au

medi Austria GmbH  
 Adamgasse 16/7  
 6020 Innsbruck  
 Austria  
 T +43 512 57 95 15  
 F +43 512 57 95 15 45  
 vertrieb@medi-austria.at  
 www.medi-austria.at

medi Belgium NV  
 Staatsbaan 77/0099  
 3945 Ham  
 Belgium  
 T: +32-11 24 25 60  
 F: +32-11 24 25 64  
 info@medibelgium.be  
 www.medibelgium.be

medi Canada Inc / médi Canada Inc  
 104-1375 Lionel-Boulet,  
 Varennes, Québec,  
 QC Canada J3X 1P7  
 T +1 450-583-3317 / +1 800 361 3153  
 F +1 888-583-6827  
 service@medicanada.ca  
 www.medicanada.ca

MAXIS a.s.,  
 medi group company  
 Slezská 2127/13  
 120 00 Prague 2  
 Czech Republic  
 T: +420 571 633 510  
 F: +420 571 616 271  
 info@maxis-medica.com  
 www.maxis-medica.com

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# M.4 X-lock®

## Intended purpose

M.4 X-lock is a brace with a frame that is intended to reinforce and stabilise  the knee joint using limitation of extension and flexion.

## Indications

All indications requiring a physiological guidance and stabilization of the knee joint and / or relief of the joint in two directions such as:

- Meniscus intervention
- Cartilage surgery
- Meniscus refixation
- Postoperative after patellar or quadriceps tendon ruptures
- Postoperative after fractures of patella
- Restraint operations of the patella (e.g. MPFL reconstruction)
- Nonsurgical therapy of patellar dislocations
- Cruciate ligament ruptures (conservative / postoperative)
- Collateral ligament ruptures (conservative / postoperative)
- Knee instability

## Contraindications

None known at present.

## Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in conjunction with sweating) or

due to the composition of the material. The immobilisation which is achievable by the device presents a risk for development of thrombosis in the lower veins of the legs, especially in the context of surgical interventions or trauma.

## Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

## Note regarding use of the knee brace

Please heed the specifications of your physician when using the knee brace, especially with respect to stressing the leg and to knee flexibility. In general, stress the operated leg only in the extended position - the knee brace joints are locked in the extended position (red area in the joint window). This applies until your physician instructs otherwise. When the brace is not stressed (e.g. when you are sitting or lying down), you can unlock the joints (green area in the joint window) and bend your leg to the extent permitted by your physician. BEFORE GETTING UP, however, make sure that the joints are locked in the extended setting and then stand up with your leg stretched (Fig. A-D). One click and the additional display (red pin is visible on the front edge of the joint) confirms that the mechanism has been snapped in correctly.

## Setting the flexion / extension limits

Extension limit wedge, silver:

0°, 10°, 20°, 30°, 45°

Flexion limit wedge, black:

0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°

Can be immobilised at:

0°, 10°, 20°, 30°, 45°

For immobilisation, first insert the extension

wedge and then the flexion wedge. For removal please proceed in the reverse order.

1. Your knee brace is provided with pre-installed 0° extension wedges.
2. In order to change the wedges, press the joint cover on the rear side (see profile) outwards.
3. You can change the wedges by removing the tweezers from the joint cover and pulling the wedges forwards (extension) or backwards (flexion) out of the joint using the tweezers. The wedges are accordingly labelled with „FLEX“ or „EXT“ and the corresponding degree. The wedges may also be distinguished by the colour. Extension wedges are silver, flexion wedges are black. Push the new wedge from the front side (extension) or from the back (flexion) into the joint, so the markings „UP“ on the joint and the wedge are matched accordingly. Upon proper locking you will hear a click sound.
4. Move the joints of the orthosis to check the position of the wedges. Put the tweezers back into the mounting and attach the joint covers by hooking the two noses of the joint cover into the front side of the joint and then pushing the cover onto the joint. The cover is firmly attached if you hear a click sound.

### Important notes

Only change the extension and flexion limits on instruction from your doctor or orthotist. In order to prevent hyperextension, always insert extension wedges or 0° wedges. The wedges define the exact excursion of the orthosis. In difficult soft-tissue conditions, e.g. where the soft-tissues are under compression, it may be necessary to insert the larger wedge to limit the joint movement to the desired excursion. In addition, make sure that the orthosis joints are correctly positioned: the middles of the joints should be level with the upper border of the patella. The joint must lie behind the midline of the leg (Fig.12).

### Shaping the orthosis

When shaping the orthosis, please ensure that the joints stay parallel. This guarantees that the orthosis will function well and run easily and that there is no excessive wear.

### Locking mechanism

**WARNING!** The locking mechanism only works if a 0° extension wedge is inserted. With other wedge settings the lock mechanism does not work. The mechanism's slider can be pulled off. When reattaching the slider, make sure that it is pushed in again from behind with the green side first.

### Fastening the belt eyes

The belt eyes have been adequately fastened and are, to some extent, deliberately mobile. Please do not tighten the screws; they are tight enough. If you apply too much force, you can break off the screws.

### Fitting instructions

- Whenever possible, put on the knee brace in the extended position. To do so, set the knee brace to the extended position and lock both joints by moving the slider forward (you can see a red area in the windows on the joints); see Figures 1+2.
- Open all the straps. Sit on the edge of a chair and put the brace on. Position the brace so that the centre of the hinges (silver wedges) are on a level with the top edge of the kneecap (Fig. 9). Make sure that the brace is not twisted on your leg.
- First fasten the lower front strap #1 (Fig. 4), then the upper rear thigh strap #2 (Fig. 5). Position the two joints close to your joint towards the front of the leg (Fig. 6).
- Now fasten the lower rear thigh strap #3 (Fig. 7). Pull the joints slightly forward again. Make sure that you do not pull the joints in front of the midline of the leg.
- Now fasten the upper front calf strap #4 (Fig. 8), then the upper back calf strap #5 (Fig. 9 + 12), and finally the front thigh strap #6 (Fig. 10).

### Care instructions

The brace cannot be worn in saltwater. If the brace is worn in chlorinated water, it should be rinsed out thoroughly with clean water afterwards. The pads are water-repellent; nevertheless, they should be wiped off and dried briefly after you take a shower – the same

applies to the straps. Soap residues, lotions and ointments can cause skin irritation and material wear.

- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

With normal use, the hinge mechanism should not require any maintenance, but it can be lubricated with PTFE-spray. Check regularly that the joint screws are secure, tighten them if necessary.

### Storage instructions

Keep the brace in a dry place and do not expose to direct sunlight.

### Important notes

This medical device (MD) is made only for multiple uses on a single patient (P). If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

### Material composition

Avional, Velours, Polyester, PU foam

### Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

### Disposal

Dispose of properly after use.

### Your medi Team

Wishes you a quick return to full fitness.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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