medi GmbH & Co. KG
Medicusstraße 1
D-95448 Bayreuth
Germany
T +49 921 912-0
F +49 921 912-780
ortho@medi.de
www.medi.de

medi Australia Pty Ltd
 Unit 4/13-21 Hallmark Street
 Pendle Hill NSW 2145
 Australia
 T +61 2 9890 8696
 F +61 2 9890 8439
 sales@medi.au
 www.medi.australia.com.au

medi Austria GmbH
 Adamgasse 16/7
 6020 Innsbruck
 Austria
 T +43 512 57 95 15
 F +43 512 57 95 15 45
 vertrieb@medi-austria.at
 www.medi-austria.at

medi Belgium NV
 Staatsbaan 77/0099
 3945 Ham
 Belgium
 T: +32-11 24 25 60
 F: +32-11 24 25 64
 info@medibelgium.be
 www.medibelgium.be

medi Canada Inc / médi Canada Inc
 104-1375 Lionel-Boulet,
 Varennes, Québec,
 QC Canada J3X 1P7
 T +1 450-583-3317 / +1 800 361 3153
 F +1 888-583-6827
 service@medicanada.ca
 www.medicanada.ca

MAXIS a.s.,
 medi group company
 Slezská 2127/13
 120 00 Prague 2
 Czech Republic
 T: +420 571 633 510
 F: +420 571 616 271
 info@maxis-medica.com
 www.maxis-medica.com

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medi Trading (Shanghai) Co. Ltd.
 Room 1703, T1, CES WEST BUND CENTER,
 No. 277, Long Lan Road
 Xuhui District
 200232 Shanghai
 People's Republic of China
 T: +86-21 50582319
 F: +86-21 50582319

UK Responsible Person
 medi UK Ltd.
 Plough Lane
 Hereford HR4 OEL
 Great Britain
 T +44-1432 37 35 00
 F +44-1432 37 35 10
 enquiries@mediuk.co.uk
 www.mediuk.co.uk

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Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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

Diagnosis and choice of the appropriate orthosis

How to select the correct protect.4 OA for your patient

1. Anatomical shape of the legs	Varus-position	Valgus-position
	 <p>The diagram shows two legs from a front view. The right leg is on the left and the left leg is on the right. In both, the knee joint is shifted towards the inner side (medial) relative to the ankle, illustrating a varus alignment.</p> <p>right leg left leg</p>	 <p>The diagram shows two legs from a front view. The right leg is on the left and the left leg is on the right. In both, the knee joint is shifted towards the outer side (lateral) relative to the ankle, illustrating a valgus alignment.</p> <p>right leg left leg</p>
2. Indication	medial osteoarthritis, bow legged medial fracture of the tibial head medial meniscus refixation medial cartilage surgery	lateral osteoarthritis, knock knee lateral fracture of the tibial head lateral meniscus refixation lateral cartilage surgery
3. Therapy	off-loading the medial joint compartment	off-loading the lateral joint compartment
4. Position of the off-loading hinge	on the lateral side of the knee (condyle)	on the medial side of the knee (condyle)
5. Push principle applied	valgus force applied	varus force applied
6. Orthosis	protect.4 OA VARUS right/left	protect.4 OA VALGUS right/left

protect.4 OA

Intended purpose

protect.4 OA is a rigid brace for relief  and stabilisation  of the knee joint.

Indications

All indications requiring unicompartimental relief by 3 point principle and additional physiological guidance and stabilization by 4 point principle and if necessary limitations of Range of motion such as:

- Severe medial or lateral gonarthrosis (at least level III) combined with ligament instability
- Unicompartimental relief at serve gonarthrosis before HTO (brace test)
- Cartilage surgery or chondroplasty
- Meniscus intervention

Contraindications

Combined lateral and medial osteoarthritis

Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in

conjunction with sweating) or due to the composition of the material. The immobilisation which is achievable by the device presents a risk for development of thrombosis in the lower veins of the legs, especially in the context of surgical interventions or trauma.

Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Setting the flexion / extension limits

Extension limit wedge (limits extension):
0°, 10°, 20°, 30°, 45°

Flexion limit wedge (limits flexion):
0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°

Can be immobilised at:

0°, 10°, 20°, 30°, 45°

The medial (inside) and the lateral (outside) wedge settings must always be the same.

1. 10° extension wedges are pre-installed in your knee brace as standard. They are found in the front aspect of the brace's hinge joint.
2. To replace the wedges first loosen the screws in the joint cover, then take the cover off.
3. To change the silver extension wedge loosen the screw and pull the wedge

- out from the front of the joint. You will find the replacement wedge in the wedge set. It is marked with the corresponding wedge setting. Slide the new extension wedge into the front of the joint. Now replace the screw in the threaded hole and screw it tight.
4. The flexion wedges are installed from the back of the joint. Select the suitable wedge from the wedge set. Insert the wedge and secure it with the appropriate screw (enclosed).
 5. Bend and flex the brace 'joint' a couple of times and check that the wedges are positioned correctly. Finally, secure the joint cover.

Important notes

Only change the extension and flexion limits on instruction from your doctor or orthotist to do so.

In order to prevent hyperextension, always insert extension wedges or 0° wedges. The wedges define the exact excursion of the orthosis. In difficult soft-tissue conditions, e.g. where the soft-tissues are under compression, it may be necessary to insert the larger wedge to limit the joint movement to the desired excursion. In addition, make sure that the orthosis joints are correctly positioned: the middles of the joints should be level with the upper border of the patella. The joint must lie behind the midline of the leg.

Fitting instructions

- Open all the straps. Sit on the edge of a chair, pull the brace up the affected leg and flex your knee to 45°. Position the joints so that the centres of the hinges (extension wedges) are about 2,5 cm above the kneecap and slightly past midline of the leg. Make sure that the brace is not twisted on your leg.

- First fasten the lower front calf strap (strap #1), then the upper back thigh strap (#2). Keep your knee in a flexed position and push both joints towards the back of the leg.
- Now fasten the lower back thigh strap (#3). Make sure that you do not have the hinge joints positioned in front of the midline of the leg.
- Now fasten the upper front calf strap (#4). Then the upper back calf strap (#5), and finally the front thigh strap (#6).
- Bend your knee slightly a few times and refasten the straps again in the order described above to prevent the brace slipping.

Adjusting Varus / Valgus load

Varus / valgus relief should be undertaken by the orthotist or by the patient - only if the patient has been given detailed instructions from a qualified person.

First sit down. Use the Allen key to set the varus / valgus adjustment hinges to the neutral position. The neutral position is reached when the cushion (at the knee joint of the brace) lies flat against the knee.

Now stand up. Set the loading to one quarter of a turn (on both adjustment screws). This setting is a good starting point when wearing the brace for the first time. You can also use the various hinge cushions supplied to adjust the support to the width of your knee if necessary.

Make sure that the brace is fitted to your knee properly. In conservative therapy (without operation) you should not wear the orthosis longer than 2 hours per day. Then increase the wearing time gradually so that the brace remains comfortable. Above all you should wear the brace for activities that would normally cause you knee joint pain. In postoperative therapy the orthosis must

always be worn when putting weight on the leg. Please ask the doctor who is treating you for his instructions. Caution: "Extreme" settings have no therapeutic value, they only cause additional discomfort.

Care instructions

The brace cannot be worn in saltwater. If the brace is worn in chlorinated water, it should be rinsed out thoroughly with clean water afterwards. Soap residues, lotions and ointments can cause skin irritation and material wear.

- Wash the paddings and covers by hand, preferably using medi clean washing agent.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

With normal use, the hinge mechanism should not require any maintenance, but it can be lubricated with PTFE spray.

Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

Material composition

Aluminum, polyamide, PU

Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

Disposal

The product can be disposed of in the domestic waste.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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