**medi GmbH & Co. KG**  
**Medicusstraße 1**  
**95448 Bayreuth**  
**Germany**  
**T +49 921 912-0**  
**F +49 921 912-780**  
**ortho@medi.de**  
**www.medi.de**

medi Australia Pty Ltd  
 83 Fennell Street  
 North Parramatta NSW 2151  
 Australia  
 T +61-2 9890 8696  
 F +61-2 9890 8439  
 sales@mediaustralia.com.au  
 www.mediaustralia.com.au

medi Austria GmbH  
 Adamgasse 16/7  
 6020 Innsbruck  
 Austria  
 T +43 512 57 95 15  
 F +43 512 57 95 15 45  
 vertrieb@medi-austria.at  
 www.medi-austria.at

medi Belgium NV  
 Staatsbaan 77/0099  
 3945 Ham  
 Belgium  
 T: +32-11 24 25 60  
 F: +32-11 24 25 64  
 info@medibelgium.be  
 www.medibelgium.be

medi Canada Inc / médi Canada Inc  
 104-1375 Lionel-Boulet,  
 Varennes, Québec,  
 QC Canada J3X 1P7  
 T +1 450-583-3317 / +1 800 361 3153  
 F +1 888-583-6827  
 service@medicanada.ca  
 www.medicanada.ca

MAXIS a.s.,  
 medi group company  
 Slezská 2127/13  
 120 00 Prague 2  
 Czech Republic  
 T: +420 571 633 510  
 F: +420 571 616 271  
 info@maxis-medica.com  
 www.maxis-medica.com

4

0515261540453

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medi Trading (Shanghai) Co. Ltd.  
 Room 105, Building No. 4,  
 Lane 500 ShengXia Road  
 Pudong New District  
 201210 Shanghai  
 People's Republic of China  
 T: +86-21 61761988  
 F: +86-21 61769607  
 infocn@medi.cn

Medi Turk Ortopedi Medikal Ic ve Dis Tic. AS  
 Mustafa Kemal Mah. 2156 Sokak No: 16 /4-5  
 06510 Çankaya Ankara  
 Turkey  
 T: +90 312 435 20 26  
 F: +90 312 434 17 67  
 info@medi-turk.com  
 www.medi-turk.com

medi Ukraine LLC  
 Evhena Sverstiuka str, 11.  
 Kiev 02002  
 Ukraine  
 T: +380 44 591 11 63  
 F: +380 44 392 73 73  
 info@medi.ua  
 http://medi.ua

medi USA L.P.  
 6481 Franz Warner Parkway  
 Whitsett, N.C. 27377-3000  
 USA  
 T +1-336 4 49 44 40  
 F +1-888 5 70 45 54  
 info@mediusa.com  
 www.mediusa.com

UK Responsible Person  
 medi UK Ltd.  
 Plough Lane  
 Hereford HR4 0EL  
 Great Britain  
 T +44-1432 37 35 00  
 F +44-1432 37 35 10  
 enquiries@mediuk.co.uk  
 www.mediuk.co.uk

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### Important notes

This medical device **MD** is made only for multiple uses on a single patient . If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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# protect.ROM telescopic

## Intended purpose

The protect.ROM telescopic is a universal knee splint for mobilisation.

## Indications

All indications which require early functional mobilisation and limitation of the range of motion of the knee joint, such as:

- Following ligament injury (post-traumatic / post-operative)
- Following meniscus injury (post-traumatic / post-operative)
- Following tendon injury (post-traumatic / post-operative)

## Contraindications

None known at present.

## Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in conjunction with sweating) or due to the composition of the material.

The immobilisation achieved by the device presents a risk for development of thrombosis in the lower veins of the legs, especially in the context of surgical interventions or trauma.

## Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

## Adjusting the flexion and extension stops

To adjust the range of movement first turn the opening in the plastic to expose one of the pins. Now press and hold the red button, at the same time turn to set the desired degree figure. Repeat this for the second pin.

Extension limitation: 0°, 10°, 20°, 30°

Flexion limitation: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°, 105°, 120°

The medial (inside) and the lateral (outside) wedge settings must always be the same.

## Setting the length

Bring the yellow slider into the 'open position'. Pull the orthosis to the desired length and close the lock mechanism again until you hear it snap into place. The medial and lateral lengths must always be equal.

## Fitting instructions

- Open all the straps and take the joint braces out of the pads.
- Wrap the calf and thigh pads tightly around the leg and fasten at the front.
- First attach the lateral (outside) brace to the pad. Make sure that the brace joint is level with the patella and the midline of the knee.
- Now attach the medial (inside) brace in the same way. Make sure that the two joints are at the same level.
- First fasten the upper thigh strap, then the lower thigh strap. Now fasten the other straps, starting with the bottom one and working upwards.

- To take the orthosis off just open the straps and pads and pull the orthosis of backwards.

---

### Important notes

Only change the extension and flexion limits on instruction from your doctor or orthotist to do so.

---

### Care instructions

Hook and loop fasteners should be closed for washing. Soap residues, lotions and ointments can cause skin irritation and material wear.

- Wash the paddings by hand, preferably using medi clean washing agent.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

---

### Storage instructions

Keep the brace in a dry place and do not expose to direct sunlight.

---

### Material composition

Aluminum, PU-foam, cotton

### Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

### Disposal

The product can be disposed of in the domestic waste.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the

EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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