**medi GmbH & Co. KG**  
 Medicusstraße 1  
 95448 Bayreuth  
 Germany  
 T +49 921 912-0  
 F +49 921 912-780  
 ortho@medi.de  
 www.medi.de

medi Australia Pty Ltd  
 Unit 4/13-21 Hallmark Street  
 Pendle Hill NSW 2145  
 Australia  
 T +61 2 9890 8696  
 F +61 2 9890 8439  
 sales@medi.au  
 www.medi.australia.com.au

medi Austria GmbH  
 Adamgasse 16/7  
 6020 Innsbruck  
 Austria  
 T +43 512 57 95 15  
 F +43 512 57 95 15 45  
 vertrieb@medi-austria.at  
 www.medi-austria.at

medi Belgium NV  
 Staatsbaan 77/0099  
 3945 Ham  
 Belgium  
 T: +32-11 24 25 60  
 F: +32-11 24 25 64  
 info@medibelgium.be  
 www.medibelgium.be

medi Canada Inc / médi Canada Inc  
 104-1375 Lionel-Boulet,  
 Varennes, Québec,  
 QC Canada J3X 1P7  
 T +1 450-583-3317 / +1 800 361 3153  
 F +1 888-583-6827  
 service@medicanada.ca  
 www.medicanada.ca

MAXIS a.s.,  
 medi group company  
 Slezská 2127/13  
 120 00 Prague 2  
 Czech Republic  
 T: +420 571 633 510  
 F: +420 571 616 271  
 info@maxis-medica.com  
 www.maxis-medica.com

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medi Trading (Shanghai) Co. Ltd.  
 Room 613, HnR Group Building  
 No.1178, Tian Yao Qiao Road  
 Xuhui District, Shanghai, 200232  
 People's Republic of China  
 T +86-21 50582318  
 info@medi.cn  
 www.medi.cn

UK Responsible Person  
 medi UK Ltd.  
 Plough Lane  
 Hereford HR4 0EL  
 Great Britain  
 T +44-1432 37 35 00  
 F +44-1432 37 35 10  
 enquiries@mediuk.co.uk  
 www.mediuk.co.uk

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### Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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# medi PT control wrap

## Intended purpose

medi PT control is a brace with adjustable hinges for correction and control of patellar tracking .

## Indications

All indications requiring correction and / or relief and protection of patellatracking in defined Range of Motion, such as:

- Patella luxation (conservative/operative)
- Patellofemoral pain syndrom
- Lateral malalignment with knee instability
- Retropatellar cartilage surgery

## Contraindications

None known to date.

## Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical

irritation of the skin (particularly in conjunction with sweating) or due to the composition of the material. The immobilisation which is achievable by the device presents a risk for development of thrombosis in the lower veins of the legs, especially in the context of surgical interventions or trauma.

## Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

## Setting the flexion / extension limits Only to be carried out by an orthopedic technician.

Extension limit wedge, EXT (limits extension): 0°, 10°, 20°, 30°, 45°

Flexion limit wedge, FLEX (limits flexion): 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°

Can be immobilised at:

0°, 10°, 20°, 30°, 45°

For immobilising the orthosis first insert the extension wedge and then the corresponding flexion wedge.

The medial (inside) and the lateral (outside) wedge settings must always be the same.

1. Your knee brace is provided with pre-installed 0° extension wedges.
2. In order to change the wedges, push the joint cover from the back toward the outside (see profile).

3. To change the wedges, remove the tweezers from the joint cover and use them to simply pull the wedges toward the front (extension) or toward the back (flexion) out of the joint. Each wedge is labelled with „FLEX“ or „EXT“ and with the respective degree setting. Another distinguishing feature is the different colour of the wedges. Extension wedges are silver and flexion wedges are black. Insert the new wedge from the front (extension) or from the back (flexion) into the joint so that the „UP“ marking on the joint and on the wedge correspond. A clicking sound confirms that the wedge has interlocked.
4. Move the brace joint and test the position of the wedges. Now put back the tweezers into the holder and attach the joint cover again by hooking the two catches of the joint cover onto the front of the joint and then pushing the cover onto the joint. A clicking sound confirms that the joint cover is fastened securely.

#### Important notes

Only change the extension and flexion limits on instruction from your doctor or orthotist. In order to prevent hyperextension, always insert extension wedges or 0° wedges. The wedges define the exact excursion of the orthosis. In difficult soft-tissue conditions, e.g. where the soft-tissues are under compression, it may be necessary to insert the larger wedge to limit the joint movement to the desired excursion. In addition, make sure that the orthosis hinges are correctly positioned: the middles of the hinges should be level with the upper border of the patella. The hinge must lie behind the midline of the leg.

#### Fitting instructions

- Open all the straps so that they will not

interfere with you putting on the brace.

- Now put your leg into a 30-degree angle flexion position. Wrap the brace tightly around your leg from the back and close it with the two Velcro strips (Image 1),
- Make sure that the kneecap is positioned centrally in the space for the patella. Position the brace so that the middles of the hinges are level with the upper border of the kneecap and just behind the lateral centre line of the leg (Image 2).
- Now close the straps in the following order:
  - Lower calf strap in the back (Image 3).
  - Lower thigh strap in the back (Image 4).
  - Patellar straps so that the kneecap is guided from the outside to the inside
    - please ensure that the patellar air cushion is positioned to the side of the kneecap and that it does not push onto the kneecap from the front (Image 5 + 6).
  - Upper thigh strap in the back (Image 7).
  - Thigh strap in the front (Image 8).
  - Calf strap in the front (Image 9).
- If possible, walk a few steps in order to check that the brace is in the correct position. Please adjust the position as required (Image 10).

---

#### Care instructions

Please close the hook and loop fasteners before washing and remove the orthoses. Soap residues, lotions and ointments can cause skin irritation and material wear.

- Wash the product by hand, preferably using medi clean washing agent.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.

- Do not dry clean.

### Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

### Material composition

Avional, polyamide, elasthane, PU

### Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

### Disposal

Dispose of properly after use.

### Your medi team

wishes you a quick return to full fitness.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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