medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth
Germany
T +49 921 912-0
F +49 921 912-780
ortho@medi.de
www.medi.de

medi Australia Pty Ltd
 83 Fennell Street
 North Parramatta NSW 2151
 Australia
 T +61-2 9890 8696
 F +61-2 9890 8439
 sales@mediaustralia.com.au
 www.mediaustralia.com.au

medi Austria GmbH
 Adamgasse 16/7
 6020 Innsbruck
 Austria
 T +43 512 57 95 15
 F +43 512 57 95 15 45
 vertrieb@medi-austria.at
 www.medi-austria.at

medi Belgium NV
 Staatsbaan 77/0099
 3945 Ham
 Belgium
 T: +32-11 24 25 60
 F: +32-11 24 25 64
 info@medibelgium.be
 www.medibelgium.be

medi Canada Inc / médi Canada Inc
 104-1375 Lionel-Boulet,
 Varennes, Québec,
 QC Canada J3X 1P7
 T +1 450-583-3317 / +1 800 361 3153
 F +1 888-583-6827
 service@medicanada.ca
 www.medicanada.ca

MAXIS a.s.,
 medi group company
 Slezská 2127/13
 120 00 Prague 2
 Czech Republic
 T: +420 571 633 510
 F: +420 571 616 271
 info@maxis-medica.com
 www.maxis-medica.com

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Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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Polyamid, Elastan, Aluminium, PU, Polyester

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protect.Co

Intended purpose

protect.Co is a knee support brace which functions based on the 4-point principle and with limitation of extension / flexion.

Indications

All indications that require reinforcement of the physiological stabilisation of the knee joint and/or stress relief of the joint, such as:

- Following injury to the collateral ligament (post-traumatic/postoperative)
- Following injury to the cruciate ligament (posttraumatic/postoperative)
- For knee instability
- Following meniscal injury
- For knee osteoarthritis

Contraindications

None known at present.

Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in conjunction with sweating) or due to the composition of the material.

Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals. Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Setting the flexion / extension limits

Extension limit wedge (limits extension):

0°, 10°, 20°, 30°, 45°

Flexion limit wedge (limits flexion):

0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°

Can be immobilised at:

0°, 10°, 20°, 30°, 45°

The medial (inside) and the lateral (outside) wedge settings must always be the same.

1. 10° extension wedges are pre-installed in your knee brace as standard. They are found in the front aspect of the brace's hinge joint.
2. To change the silver extension wedge loosen the screw and pull the wedge out from the front of the joint. You will find the replacement wedge in the wedge set. It is marked with the corresponding wedge setting. Slide the new extension wedge into the front of the joint. Now replace the screw in the threaded hole and screw it tight.
3. The flexion wedges are installed from the back of the joint. Select the suitable wedge from the wedge set. Insert the wedge and secure it with the appropriate screw (enclosed).
4. Bend and flex the brace 'joint' a couple of times and check that the

wedges are positioned correctly.
Finally, secure the joint cover.

Important notes

Only change the extension and flexion limits on instruction from your doctor or orthotist to do so. In order to prevent hyperextension, always insert extension wedges or 0° wedges. The wedges define the exact excursion of the orthosis. In difficult soft-tissue conditions, e.g. where the soft-tissues are under compression, it may be necessary to insert the larger wedge to limit the joint movement to the desired excursion. In addition, make sure that the orthosis joints are correctly positioned: the middles of the joints should be level with the upper border of the patella. The joint must lie behind the midline of the leg.

Fitting instructions

- Undo all straps and the two hook and loop strips.
- Now put your leg into a 45-degree angle flexion position. Wrap the brace tightly around your leg from the back and close it with the two hook and loop strips (Image 1).
- Make sure that the kneecap is positioned centrally in the space for the patella. Position the brace so that the middles of the hinges are level with the upper border of the kneecap and just behind the lateral centre line of the leg (Image 2).
- First, fasten the lower calf strap (Image 3) and then the upper thigh strap (Image 4).
- Next, fasten the strap directly below the kneecap (Image 5).
- Finally, fasten the remaining lower thigh strap (Image 6).
- If possible, walk a few steps in order to check that the brace is in the correct position (Image 7).
Please adjust the position as required.

Care instructions

Please close the hook and loop fasteners before washing and remove the orthoses. Soap residues, lotions and ointments can cause skin irritation and material wear.

- Wash the product by hand, preferably using medi clean washing agent.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

Storage instructions

Keep the brace in a dry place and do not expose to direct sunlight.

Material composition

Polyamide, elastane, aluminum, PU, polyester

Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

Disposal

The product can be disposed of in the domestic waste.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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