medi GmbH & Co. KG  
 Medicusstraße 1  
 95448 Bayreuth  
 Germany  
 T +49 921 912-0  
 www.medi.de

medi Australia Pty Ltd  
 83 Fennell Street  
 North Parramatta NSW 2151  
 Australia  
 T +61-2 9890 8696  
 F +61-2 9890 8439  
 sales@mediaustralia.com.au  
 www.mediaustralia.com.au

medi Austria GmbH  
 Adamgasse 16/7  
 6020 Innsbruck  
 Austria  
 T +43 512 57 95 15  
 F +43 512 57 95 15 45  
 vertrieb@medi-austria.at  
 www.medi-austria.at

medi Belgium NV  
 Staatsbaan 77/0099  
 3945 Ham  
 Belgium  
 T: +32-11 24 25 60  
 F: +32-11 24 25 64  
 info@medibelgium.be  
 www.medibelgium.be

medi Canada Inc / médi Canada Inc  
 104-1375 Lionel-Boulet,  
 Varennes, Québec,  
 QC Canada J3X 1P7  
 T +1 450-583-3317 / +1 800 361 3153  
 F +1 888-583-6827  
 service@medicanada.ca  
 www.medicanada.ca

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








medi

## circaid® juxtacures® Inelastic compression system

Gebrauchsanweisung. Instructions for use. Instructions d'utilisation. Instrucciones de uso. Instruções para aplicação. Istruzioni per l'uso. Gebruiksaanwijzing. Brugsvejledning. Bruksanvisning. Návod k použití. Инструкция по использованию. Instrukcja zakłádania. Használati útmutató. دليل الاستخدام. Інструкція з використання. 使用说明. אהאנת שימוש. Bruksanvisning. Kasutusjuhend. Käyttöohje. Arahán penggunaan.  
 20–30, 30–40, or 40–50mmHg\*

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


### Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If undue pain or an unpleasant sensation develops while you are wearing the product, please consult your doctor or orthotist immediately. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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medi Norway AS  
Vestheimvegen 35  
4250 Kopervik  
Norway  
T: +47 52 84 45 00  
F: 47 52 84 45 09  
post@medinorway.no  
www.medinorway.no

medi Polska Sp. z o.o.  
ul. Łabędzka 22  
44-121 Gliwice  
Poland  
T: +48-32 230 60 21  
F: +48-32 202 87 56  
info@medi-polska.pl  
www.medi-polska.pl

medi Bayreuth Unipessoal, Lda  
Rua do Centro Cultural, no. 43  
1700-106 Lisbon  
Portugal  
T +351-21 843 71 60  
F +351-21 847 08 33  
medi.portugal@medibayreuth.pt  
www.medi.pt

medi RUS LLC  
Business Center NEO GEO  
Butlerova Street 17  
117342 Moscow  
Russia  
T +7-495 374 04 56  
F +7-495 374 04 56  
info@medirus.ru  
www.medirus.ru

medi Sweden AB  
Box 6034  
19206 Sollentuna  
Sweden  
T: +46-8-969798  
F: +46-8-6266870  
info@medi.se  
www.medi.se

medi Ukraine LLC  
Evhena Sverstiuka str, 11.  
Kiev 02002  
Ukraine  
T: +380 44 591 11 63  
F: +380 44 392 73 73  
info@medi.ua  
http://medi.ua

medi USA L.P.  
6481 Franz Warner Parkway  
Whitsett, N.C. 27377-3000  
USA  
T +1-336 4 49 44 40  
F +1-888 5 70 45 54  
info@mediusa.com  
www.mediusa.com

medi Trading (Shanghai) Co. Ltd.  
Room 105, Building No. 4,  
Lane 500 ShengXia Road  
Pudong New District  
201210 Shanghai  
People's Republic of China  
T: +86-21 61761988  
F: +86-21 61769607  
info@medi.cn

Medi Turk Ortopedi Medikal Iç ve Dis Tic. AS  
Mustafa Kemal Mah. 2156 Sokak No: 16 /4-5  
06510 Çankaya Ankara  
Turkey  
T: +90 312 435 20 26  
F: +90 312 434 17 67  
info@medi-turk.com  
www.medi-turk.com

MAXIS a.s.  
medi group company  
Slezská 2127/13  
120 00 Prague 2  
Czech Republic  
T: +420 571 633 510  
F: +420 571 616 271  
info@maxis-medica.com  
www.maxis-medica.com

Tachezy Sanit Srl.  
Via Druso 48/a  
39100 Bolzano  
Italy  
T: +39 0471 288 068  
F: +39 0471 282 773  
info@tachezysanit.com  
www.tachezysanit.com

medi Japan K.K.  
5-1 Nihonbashi Kabutocho Chuoku  
103-0026 Tokyo  
Japan  
T: +81 3 6778 2590  
F: +81 3 5847 7901  
info@medi-japan.jp  
www.medi-japan.co.jp

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# circaid® juxtacures®

### Intended purpose

The compression system is designed to provide compression to the lower leg for patients with venous and lymphatic disorders. The device consists of a series of juxtapositioned, inelastic bands extending from the central portion of the garment. The wrap contours the leg and is designed to cover from just above the ankle bone to just below the knee.

### Indications

- Chronic venous insufficiency
  - Varicose Veins
  - Lipodermatosclerosis
  - Ulcus cruris venosum/  
mixed Ulcus cruris
- Deep vein thrombosis/Thrombosis Prevention
- Acute DVT
- Post Thrombotic Syndrome
- Post Sclerotherapy
- Dependent oedema

### Contraindications

- Severe peripheral arterial disease
- Uncontrolled Congestive Heart Failure
- Septic phlebitis
- Phlegmasia cerula dolens
- Untreated infection in the leg
- Any circumstance where increased venous and lymphatic return is undesirable

### Risks / Side effects

- Mild to Moderate Peripheral Arterial Disease
- Impaired Sensation - must be able to sense applied pressure
- Intolerance to material(s)
- Pediatric & assisted use - must be able to sense & communicate applied pressure

### Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

### Note on wearing

*Note: Unless your treating physician orders or recommends otherwise, we recommend wearing the the legwear together with the respective footwear included in the set. If in doubt, please consult your treating physician.*

Always ensure that the correct prescribed pressure range is being applied. Slightly loosen the bands of the garment for night-time wear. Each vertical BPS line on the bands may be at a different location on the leg depending on the size and shape of the leg. They do not need to be aligned down the front center of the leg or in any other specific location on the leg. If experiencing any pain, or discomfort, reduce the compression (mmHg) applied or remove device. Inform your clinician of your action. When bands are placed correctly they should slightly overlap, leaving no space between bands. Do not tuck bands. Hook tabs should never touch the skin!

### Selecting the circaid juxtacures length

*circaid juxtacures is available in short length (28cm), standard length (33 cm) and long length (38cm).*

The appropriate length is determined by the patient's floor to knee crease measurement. The patient should be measured in standing position.

If the patients measures less than 44 cm, select the short legging. If the patient measures between 44 cm and 49 cm, select the standard legging, and if the patient measures more than 49 cm, select the long legging.

It is important that the length of the juxtacures covers the calf for reinforcement of the calf-muscle pump.

**Note:** *circaid juxtacures will fit a 19 cm - 42 cm ankle circumference and a 25 cm - 64 cm calf circumference. circaid juxtacures will not accommodate more than a 30 cm difference in circumference from*

the ankle to the calf measurement. If your patient exceeds these measurements, alternative compression garments are available for order.

The circaid compression ankle will fit a maximum arch circumference of 36 cm for the Standard size and 41 cm for the Large size, the maximum ankle circumference is 27 cm for the Standard size and 42 cm for the Large size. If your patient exceeds these measurements, further foot options are available for order.

### Customizing the circaid juxtacures

- Separate the body and spine from one another, unroll the bottom band and locate the number markings along the bottom and top edge of the body piece.
- Locate the patient's ankle circumference measurement amongst the printed numbers along the bottom edge. Then, locate the patient's calf circumference measurement amongst the printed numbers along the top edge of the garment.
- Take the spine and place the bottom hook tab on the appropriate number (2a). Take the top hook tab of the spine and place it on the appropriate number (2b). Adjust the center spine tab to create a smooth spine curve (2c).

**Note:** There is no one correct spine shape. The garment may not lay flat once the spine has been set (2d). This is to be expected when dealing with more extreme spine angles. If there is a large difference between the two circumferences, the spine may not reach both the bottom and top edges of the body piece. In this case, center the spine (2e). When dealing with more extreme spine angles, excess material sticking out beyond the top of the spine piece can be trimmed (2f).

- With the spine in place, fold over the smaller garment piece. Cut the excess material from the body piece leaving approximately 5 cm of overlap (3). Be sure not to cut any other part of the garment.

**Note:** If desired, attach the separate hook tab stays along the spine for additional adhesion (4).

**Tip:** You may also trace the spine to easily reattach it in its original position in the unlikely event of the spine becoming detached.

- Take the BPS-card and highlight the range that contains the patient's ankle circumference. If the patient's

measurement falls between two ranges on the card, 30.5 cm for example, choose the smaller range.

### Donning the circaid juxtacures

- Slide the circaid undersleeve onto the lower leg and over the knee. The cuffed edge of the undersleeve should cover the heel. Ensure that there are no wrinkles in the undersleeve (5a).
- Align the garment just above the ankle bone and below the knee crease. The black side of the material should face the skin (5b).  
Unroll the top two bands and loosely secure them to the outside of the garment to anchor it in place (5c). This enables you to easily raise or lower the garment for correct positioning.

**Note:** The bottom most band can be identified by the circaid-logo on the beige side of the material.

- Unroll the bottom two bands and secure the bottom most band to the outside of the legging to a firm and comfortable compression level (5d).
- While holding the second band, detach the next band (5e). Secure the second band and continue this process with the other bands, going up the leg (5f).
- The garment should lay flat and wrinkle-free against the leg. Adjust each band to a firm and comfortable compression level (5g). Ensure that the bands are placed in an alternating order.
- Once all bands are secured, inspect the garment for any gaps or creases. Adjust the bands as necessary.

**Note:** Once the correct pressure is set (5h), you may fold the top and bottom part of the circaid undersleeve over the top and bottom edge of the circaid juxtacures garment (5i).

### Donning the compression ankle

First fold the bottom part of the circaid undersleeve over the bottom edge of the circaid lower leg component (5i). Slide the circaid compression ankle onto the foot and ensure that it is pulled high enough so that it overlaps with the circaid garment (5j). Ensure that there are no wrinkles in the ankle.

### Dooffing the circaid juxtacures

To remove the garment, remove the circaid compression ankle first. Then unfold the circaid undersleeve. Detach all bands of the legging starting with the top band. Fold each band onto itself to keep the hook and loop tabs lint-free. Remove the garment and then remove the circaid undersleeve.

## Using the built-in pressure system guide card

- Locate the built-in pressure card in your packaging.
- Identify the appropriate edge of the card and pressure scale. The edge of the card is determined by ankle circumference. It is recommended to highlight the correct edge and scale.
- Starting with the bottom band, line up the black triangle on the card with one of the BPS lines on the bottom band.
- Note where the second BPS line lines up with the card's compression ranges (6a).
- If the BPS line on the garment either falls short (6b) or goes beyond the correct compression range (6c), readjust the band as necessary so that the second BPS line on the garment is aligned with the prescribed compression range (6a). Completely loosen the band before reapplying it if too much tension has been applied.
- Repeat the steps 3 - 5 for each band going up the leg. Be sure to adjust all bands to the same tension to ensure graduated compression. Bands may need adjusting throughout the day to maintain the prescribed compression range

## Follow-up visits

Remind the patient to bring the entire circaid juxtacures system to every follow-up appointment for proper adjustment. At every follow-up visit with the patient, the circaid juxtacures should be checked for necessary adjustments.

## Reduction of the leg circumferences

Should the patient reduce significantly, ankle and calf circumference measurements can be re-taken and the garment can be adjusted to the patient's new size by resetting the spine per customizing steps 1 - 4.

Once the device has been properly adjusted, re-apply the legging to the patient's limb. Be sure to also check whether or not the patient should be using a different edge of the BPS-card as the appropriate range on the card may have changed. Make the appropriate adjustments, if necessary.

**Tip** Once the patient's ulcer has healed, circaid juxtalite and mediaven round-knitted stockings are available for order to prevent ulcer recurrence.

## Care instructions

The garment can be machine washed using a gentle cycle and dried on a low

heat setting. We recommend using a laundry bag. To extend the life of the garment, hand wash and drip dry.

---

|   |                      |
|---|----------------------|
|  | Machine wash warm    |
|  | Do not bleach        |
|  | Tumble dry, low heat |
|  | Do not iron          |
|  | Do not dry clean     |
|  | Do not wring         |

---

## Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

## Lifespan

6 months - due to material wear-and-tear the medical efficacy can only be expected for a defined lifespan. This assumes correct handling of the garment (e.g. proper care, application, and removal).

## Material composition

55% Nylon, 39% Polyurethane, 5% Elastane, 1% other fibres

## Disposal

The product can be disposed of in the domestic waste.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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