medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth
Germany
T +49 921 912-0
www.medi.de

medi Australia Pty Ltd
 83 Fennell Street
 North Parramatta NSW 2151
 Australia
 T +61-2 9890 8696
 F +61-2 9890 8439
 sales@mediaustralia.com.au
 www.mediaustralia.com.au

medi Austria GmbH
 Adamgasse 16/7
 6020 Innsbruck
 Austria
 T +43 512 57 95 15
 F +43 512 57 95 15 45
 vertrieb@medi-austria.at
 www.medi-austria.at

medi Belgium NV
 Staatsbaan 77/0099
 3945 Ham
 Belgium
 T: +32-11 24 25 60
 F: +32-11 24 25 64
 info@medibelgium.be
 www.medibelgium.be

medi Canada Inc / médi Canada Inc
 104-1375 Lionel-Boulet,
 Varennes, Québec,
 QC Canada J3X 1P7
 T +1 450-583-3317 / +1 800 361 3153
 F +1 888-583-6827
 service@medicanada.ca
 www.medicanada.ca

MAXIS a.s.,
 medi group company
 Slezská 2127/13
 120 00 Prague 2
 Czech Republic
 T: +420 571 633 510
 F: +420 571 616 271
 info@maxis-medica.com
 www.maxis-medica.com

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medi Trading (Shanghai) Co. Ltd.
 Room 105, Building No. 4,
 Lane 500 ShengXia Road
 Pudong New District
 201210 Shanghai
 People's Republic of China
 T: +86-21 61761988
 F: +86-21 61769607
 infocn@medi.cn

Medi Turk Ortopedi Medikal İc ve Dis Tic. AS
 Mustafa Kemal Mah. 2156 Sokak No: 16 /4-5
 06510 Çankaya Ankara
 Turkey
 T: +90 312 435 20 26
 F: +90 312 434 17 67
 info@medi-turk.com
 www.medi-turk.com

medi Ukraine LLC
 Evhena Sverstiuka str, 11.
 Kiev 02002
 Ukraine
 T: +380 44 591 11 63
 F: +380 44 392 73 73
 info@medi.ua
 http://medi.ua

UK Responsible Person
 medi UK Ltd.
 Plough Lane
 Hereford HR4 OEL
 Great Britain
 T +44-1432 37 35 00
 F +44-1432 37 35 10
 enquiries@mediuk.co.uk
 www.mediuk.co.uk

medi

circaid® reduction kit™ arm inelastic compression system


Gebrauchsanweisung. Instructions for use. Mode d'emploi. Instrucciones de uso. Instruções para aplicação. Istruzioni per l'uso. Gebruiksaanwijzing. Bruksanvisning. Инструкция по использованию. Instrukcja zakładania. دليل الاستخدام. Használati útmutató. 使用说明. הוראות שימוש. Bruksanvisning. Kasutusjuhend. Käyttöohje.

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

Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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Polyurethan, Nylon, Elasthan, Polyethylen

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circaid® reduction kit arm

Intended purpose

The compression system is designed to provide compression to the arm for patients with venous and lymphatic disorders. The compression system is designed to provide compression to the hand for patients with venous and lymphatic disorders.

Indications

- Lymphedema
- Other forms of edema
- Lipedema

Contraindications

- Severe Peripheral Arterial Disease
- Uncontrolled Congestive Heart Failure
- Septic Phlebitis
- Phlegmasia Cerulea Dolens
- Untreated infection
- Any circumstance where increased venous and lymphatic return is undesirable
- Suspected or known untreated acute deep vein thrombosis

Risks / Side Effects

- Mild to Moderate Peripheral Arterial Disease
- Impaired Sensation - must be able to sense applied pressure
- Intolerance to material(s)
- Extensive weeping dermatoses
- Pediatric & assisted use - must be able to sense & communicate applied pressure

Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Fitting instructions

- Apply the appropriate finger compression, if needed (Fig. 1).
- The undersleeve is worn for added comfort and when exudate is expected (Fig. 2).
- To fit the circaid reduction kit arm, slide the sleeve onto the patient's arm using the two inner bands. The black side of the material should face the skin and the sleeve should extend from the wrist to the axilla.

Tip: If the sleeve is too long for the patient, it can be trimmed by 3 cm and 7 cm. Trim along one of the indicated cut lines on the black side of the material (Fig.3). Trim the corresponding hook tab as needed.

The hand wrap may be applied after the reduction kit arm.

- Adjust the two inner bands to fit snugly on the patient's arm. The inner bands should not provide compression.
- Trim the excess material from the inner bands and refasten the bands (Fig.5).
- Starting with the wrist band, pull each band across the arm and note the amount it overlaps with the opposing side of the sleeve. Mark the band where it overlaps about 3 cm with the sleeve and trim the band. Repeat this step with all of the bands (Fig.6).
- To tighten the circaid reduction kit arm, instruct the patient to pull the wrist band to a firm and comfortable compression level across the arm and to attach it to the outside of the sleeve using the provided hook tabs (Fig.7). Be sure the patient does not tuck the bands under the sleeve. The hook tab should never make contact with the skin or under-

sleeve. Repeat this process for the remaining bands. Instruct the patient to angle each band individually to achieve a snug and wrinkle-free fit while ensuring that there are no gaps between any of the bands of the garment.

Tip: Once the patient reduces significantly, the garments can be adjusted to the patient's new size by cutting the bands and resetting the hook tabs following steps 3-5.

- Bands should be pulled to a firm but comfortable tension to provide therapeutic compression within the tolerance level of the patient. If the patient and therapist desire, the BTS card can be used to identify and record the appropriate tension settings (Fig.8). It is not required to have all bands set to the same tension level. If more than one level of tension is appropriate for treatment, use the BTS card to identify and record the appropriate tension levels. To do so, align one of the vertical guidelines on the band being calibrated with the triangle on the BTS card. Then record where the second vertical guideline falls. This will allow the patient to safely repeat the settings on their own, away from the clinic. Record tension setting for each set of markings so the patient can easily repeat the process and ensure it is donned to the same or higher tension.

Reminders for the patient

- The arm component should feel firm and comfortable.
- The patient may re-adjust the bands as necessary throughout the day, but should ensure that all vertical guidelines are set to the appropriate level of tension.
- Each set of vertical guidelines may be at a different location on the arm depending on the size and shape of the arm.
- For improved comfort, the patient may slightly loosen the bands of the arm component for nighttime wear.
- When bands are placed in parallel they should slightly overlap, leaving no space between bands.

Care instructions

The garment can be machine washed using a gentle cycle and dried on a low heat setting. Please roll up the bands before washing, the hook-and-loop closure must not be open. We recommend using a laundry bag. To extend the life of the garment, hand wash and drip dry.

 Machine wash warm

 Do not bleach

 Tumble dry, low heat

 Do not iron

 Do not dry clean

 Do not wring

Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

Lifespan

2 months - due to material wear-and-tear the medical efficacy can only be expected for a defined lifespan. This

assumes correct handling of the garment (e.g. proper care, application, and removal).

Material composition

Polyurethane, Nylon, Elastane, Polyethylene

Disposal

The product can be disposed of in the domestic waste.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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