medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth
Germany
T +49 921 912-0
www.medi.de

medi Australia Pty Ltd
83 Fennell Street
North Parramatta NSW 2151
Australia
T +61-2 9890 8696
F +61-2 9890 8439
sales@mediaustralia.com.au
www.mediaustralia.com.au

medi Austria GmbH
Adamgasse 16/7
6020 Innsbruck
Austria
T +43 512 57 95 15
F +43 512 57 95 15 45
vertrieb@medi-austria.at
www.medi-austria.at

medi Belgium NV
Staatsbaan 77/0099
3945 Ham
Belgium
T: + 32-11 24 25 60
F: +32-11 24 25 64
info@medibelgium.be
www.medibelgium.be

medi Canada Inc / médi Canada Inc
104-1375 Lionel-Boulet,
Varenes, Québec,
QC Canada J3X 1P7
T +1 450-583-3317 / +1 800 361 3153
F +1 888-583-6827
service@medicanada.ca
www.medicanada.ca

MAXIS a.s.,
medi group company
Slezská 2127/13
120 00 Prague 2
Czech Republic
T: +420 571 633 510
F: +420 571 616 271
info@maxis-medica.com
www.maxis-medica.com

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
medi

circaid® juxtalite® lower leg inelastic compression system

Gebrauchsanweisung. Instructions for use. Mode d'emploi.
Instrucciones de uso. Instruções para aplicação. Istruzioni per
l'uso. Gebruiksaanwijzing. Bruksveiledning. Brugsanvisning.
Návod k použití. Uputa za uporabu. Инструкция по
использованию. Instrukcja zakładania. Használati útmutató.
إرشادات الاستخدام. Інструкція з використання. 使用说明. 使用說明.
שימוש. Bruksanvisning. Kasutusjuhend. Käyttöohje.

20 – 30, 30 – 40, 40 – 50 mmHg*

circaid juxtalite lower leg

circaid juxtalite ankle foot
wrap

circaid juxtalite with ankle
foot wrap

circaid juxtalite lower leg with
circaid undersleeve lower leg

circaid juxtalite with
undersleeve ankle

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Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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


medi Trading (Shanghai) Co. Ltd.
Room 1703, T1, CES WEST BUND CENTER,
No. 277, Long Lan Road
Xuhui District
200232 Shanghai
People's Republic of China
T: +86-21 50582319
F: +86-21 50582319

Medi Turk Ortopedi Medikal Iç ve Dis Tic. AS
Mustafa Kemal Mah. 2156 Sokak No: 16/4-5
06510 Çankaya Ankara
Turkey
T: +90 312 435 20 26
F: +90 312 434 17 67
info@medi-turk.com
www.medi-turk.com

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circaid® juxtalite® lower leg

Intended purpose

The compression system is designed to provide compression to the lower leg for patients with venous disorders. The legging consists of a series of juxtapositioned, inelastic bands extending from the central portion of the wrap. The wrap contours the leg and is designed to cover from just above the ankle bone to just below the knee.

Indications

- Chronic venous insufficiency
 - Varicose veins
 - Lipodermatosclerosis
 - Venous leg ulcer/Mixed leg ulcer
- Acute DVT
- Deep vein thrombosis/Thrombosis Prevention
- Post Thrombotic Syndrome
- Post Sclerotherapy
- Dependent oedema

Contraindications

- Severe Peripheral Arterial Disease
- Uncontrolled Congestive Heart Failure
- Septic Phlebitis
- Phlegmasia Cerulea Dolens
- Untreated infection
- Any circumstance where increased venous and lymphatic return is undesirable

Risks / Side effects

- Mild to Moderate Peripheral Arterial Disease
- Impaired Sensation - must be able to sense applied pressure
- Intolerance to material(s)
- Pediatric & assisted use - must be able to sense & communicate applied pressure

Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Note on wearing

Average compression levels are at the

ankle. Due to the inelastic nature of the material, compressions levels will vary with position. Always ensure that the correct prescribed pressure range is being applied. Slightly loosen the bands of the garment for night-time wear. Each vertical BPS line on the bands may be at a different location on the leg depending on the size and shape of the leg. They do not need to be aligned down the front center of the leg or in any other specific location on the leg. If experiencing any pain, or discomfort, reduce the compression (mmHg) applied or remove device. Inform your clinician of your action. When bands are placed correctly they should slightly overlap, leaving no space between bands. Do not tuck bands. The Velcro tabs should never touch the skin!

Donning instructions

- Slide the circaid undersleeve onto the lower leg and over the knee. The cuffed edge of the undersleeve should cover the heel. Ensure that there are no wrinkles in the undersleeve.
- Align the garment just above the ankle bone and below the knee crease. The black side of the material should face the skin. Unroll the top two bands and loosely secure them to the outside of the garment to anchor the legging in place. This enables you to easily raise or lower the garment for correct positioning.
- Unroll the bottom two bands and secure the bottom most band to the outside of the legging to a firm and comfortable compression level.
- While holding the second band, detach the next band. Secure the second band and continue this process with the other bands, going up the leg.
- The garment should lay flat and wrinkle free against the leg. Adjust each band to a firm and comfortable compression level. Ensure that the bands are placed in an alternating order.
- Once all bands are secured, inspect the garment for any gaps or creases. Adjust the bands as necessary.

Note: Once the correct pressure is set using the BPS card, you may fold the top and bottom part of the circaid undersleeve over the top and bottom edge of the circaid juxtalite garment.

Donning the circaid compression anklet.

First fold the bottom part of the circaid undersleeve over the bottom edge of the lower leg component (5i). Slide the circaid compression anklet onto the foot and ensure that it is pulled high enough so that it overlaps with the circaid juxtalite garment. Ensure that there are no

wrinkles in the anklet.

Doffing instructions

To remove the garment, remove the circaid compression anklet first. Then unfold the circaid undersleeve. Detach all bands of the legging starting with the top band. Fold each band onto itself to keep the Velcro tabs lint-free. Remove the garment and then remove the circaid undersleeve. This garment is designed to provide gradient compression based on these average ankle compression levels.

Using the Built-In-Pressure card

This garment is designed to provide gradient compression based on these average ankle compression levels.

- Locate the Built-In-Pressure card in your packaging.
- Identify your prescribed pressure scale using the color-coded system.
- Starting with the bottom band, line up the vertical line on the card with one of the BPS lines on the bottom band.
- Note where the second BPS line lines up with the card's compression ranges (6a).
- If the BPS line on the garment either falls short (6b) or goes beyond the correct compression range (6c), readjust the band as necessary so that the second BPS line on the garment is aligned with the prescribed compression range (6a). Completely loosen the band before reapplying it if too much tension has been applied.
- Repeat the steps 3 - 5 for each band going up the leg. Be sure to adjust all bands to the same compression level. Bands may need adjusting throughout the day to maintain the prescribed compression range.

Care instructions

The garment can be machine washed using a gentle cycle and dried on a low heat setting. We recommend using a laundry bag. To extend the life of the garment, hand wash and drip dry.

-
-  Machine wash warm
 -  Do not bleach
 -  Tumble dry, low heat
 -  Do not iron
 -  Do not dry clean
 -  Do not wring
-

Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

Lifespan

6 months - due to material wear-and-tear the medical efficacy can only be expected for a defined lifespan. This assumes correct handling of the garment (e.g. proper care, application, and removal).

Material composition

Nylon, Polyurethane, Elastane

Disposal

The product can be disposed of in the domestic waste.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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