medi GmbH & Co. KG
 Medicusstraße 1
 95448 Bayreuth
 Germany
 T +49 921 912-0
 www.medi.de

medi Australia Pty Ltd
 Unit 4/13-21 Hallmark Street
 Pendle Hill NSW 2145
 Australia
 T +61 2 9890 8696
 F +61 2 9890 8439
 sales@medi.au
 www.mediaustralia.com.au

medi Austria GmbH
 Adamgasse 16/7
 6020 Innsbruck
 Austria
 T +43 512 57 95 15
 F +43 512 57 95 15 45
 vertrieb@medi-austria.at
 www.medi-austria.at

medi Belgium NV
 Staatsbaan 77/0099
 3945 Ham
 Belgium
 T: +32-11 24 25 60
 F: +32-11 24 25 64
 info@medibelgium.be
 www.medibelgium.be

medi Canada Inc / médi Canada Inc
 104-1375 Lionel-Boulet,
 Varennes, Québec,
 QC Canada J3X 1P7
 T +1 450-583-3317 / +1 800 361 3153
 F +1 888-583-6827
 service@medicanada.ca
 www.medicanada.ca

MAXIS a.s.,
 medi group company
 Slezská 2127/13
 120 00 Prague 2
 Czech Republic
 T: +420 571 633 510
 F: +420 571 616 271
 info@maxis-medica.com
 www.maxis-medica.com

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medi. I feel better.

circaid® juxtafit®
hand wrap

circaid® juxtafit®
essentials glove

circaid® juxtafit®
essentials glove with
dorsum strap

circaid® juxtafit®
essentials open palm
glove

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Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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Step 5


Step 6


circaid® juxtafit® essentials glove


Step 1

Step 2

Step 3

Step 4

circaid® juxtafit® essentials golve with dorsum strap

Step 1

Step 2

Step 3

Step 4

circaid® juxtafit® essentials glove removal strap

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circaid juxtafit essentials hand wrap

Intended purpose

This garment is designed to provide compression to the hand of patients with venous and lymphatic disorders.

Indications

- Lymphedema

Contraindications

- Severe Peripheral Arterial Disease
- Uncontrolled Congestive Heart Failure
- Septic Phlebitis
- Phlegmasia Cerulea Dolens
- Untreated infection
- Any circumstance where increased venous and lymphatic return is

undesirable

- Suspected or known untreated acute DVT

Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Risks / Side Effects

- Mild to Moderate Peripheral Arterial Disease
- Impaired Sensation - must be able to sense applied pressure
- Intolerance to material(s)
- Pediatric & assisted use - must be able to sense & communicate applied pressure

Donning instructions

- Slide the thumb through the hand wrap's thumb hole. The seam should run from the base of the thumb to the wrist. The curved portion of the wrap should be in the palm of the hand.
- Unroll the band closest to the fingers and secure it with the velcro tab to the

hand wrap with firm and comfortable compression. The top edge of the hand wrap should line up with the base of the fingers.

- Wrap the remaining band around the wrist and secure the Velcro tab. Be sure not to over tighten the wrist band so that fluid can freely move out of the hand. The circaid juxtafit essentials hand wrap should overlap with the circaid juxtafit essentials arm.

Note: Once all bands are secured inspect the garment for any gaps or creases. Adjust the bands as necessary.

If using multiple circaid juxtafit garments, only lightly apply bands overlapping the legging below to ensure gradient compression. Set the remaining bands to firm and comfortable compression or prescribed compression level.

Wearing recommendation

Slightly loosen the bands of the garment for night-time wear. If experiencing any pain, or discomfort, reduce the compression (mmHg) applied or remove device. Inform your clinician of your action. When bands are placed correctly they should slightly overlap, leaving no space between bands. Do not tuck bands. Hook tabs should never touch the skin!

Care instructions

The garment can be machine washed using a gentle cycle and dried on a low heat setting. If doing so, it is recommended that you place the garment in a pil-

low case to aid in keeping the hook and loop tap lint-free and to protect other clothing from getting damaged.

To extend the life of the garment, hand wash and drip dry.

 Machine wash warm

 Do not bleach

 Tumble dry, low heat

 Do not iron

 Do not dry clean

 Do not wring

Storage instructions

Store the product in a dry place and do not expose it to direct sunlight.

Lifespan

3 months – due to material wear and/or slackening, medical efficacy can only be guaranteed for a defined lifespan. This assumes correct handling of the product (e.g. during care and when putting on and taking off the product).

Material composition

Nylon, Polyurethane, Elastane, Other Fibres

Disposal

The product can be disposed
Dispose of properly after use.

In the event of any complaints regarding the product such as damage to the fabric

or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

FRANÇAIS / FRENCH

circaid juxtafit essentials hand wrap

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medi Trading (Shanghai) Co. Ltd.
medi Ukraine LLC