**medi GmbH & Co. KG**  
 Medicusstraße 1  
 95448 Bayreuth  
 Germany  
 T +49 921 912-0  
 www.medi.de

medi Australia Pty Ltd  
 83 Fennell Street  
 North Parramatta NSW 2151  
 Australia  
 T +61-2 9890 8696  
 F +61-2 9890 8439  
 sales@mediaustralia.com.au  
 www.mediaustralia.com.au

medi Austria GmbH  
 Adamgasse 16/7  
 6020 Innsbruck  
 Austria  
 T +43 512 57 95 15  
 F +43 512 57 95 15 45  
 vertrieb@medi-austria.at  
 www.medi-austria.at

medi Belgium NV  
 Staatsbaan 77/0099  
 3945 Ham  
 Belgium  
 T: +32-11 24 25 60  
 F: +32-11 24 25 64  
 info@medibelgium.be  
 www.medibelgium.be

medi Canada Inc / médi Canada Inc  
 104-1375 Lionel-Boulet,  
 Varennes, Québec,  
 QC Canada J3X 1P7  
 T +1 450-583-3317 / +1 800 361 3153  
 F +1 888-583-6827  
 service@medicanada.ca  
 www.medicanada.ca

MAXIS a.s.,  
 medi group company  
 Slezská 2127/13  
 120 00 Prague 2  
 Czech Republic  
 T: +420 571 633 510  
 F: +420 571 616 271  
 info@maxis-medica.com  
 www.maxis-medica.com

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medi. I feel better.

circaid® juxtafit®  
hand wrap

circaid® juxtafit®  
essentials glove

circaid® juxtafit®  
essentials glove with  
dorsum strap

circaid® juxtafit®  
essentials open palm  
glove

### Important notes

This medical device **MD** is made only for multiple uses on a single patient . If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If undue pain or an unpleasant sensation develops while you are wearing the product, please consult your doctor or orthotist immediately. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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Step 5


Step 6


## circaid® juxtafit® essentials glove


Step 1

Step 2

Step 3

Step 4

## circaid® juxtafit® essentials golve with dorsum strap

Step 1

Step 2

Step 3

Step 4

## circaid® juxtafit® essentials glove removal strap

Step 1

Step 2

Step 3

medi Bayreuth Espana SL  
C/Canigo 2-6 bajos  
Hospitalet de Llobregat  
08901 Barcelona  
Spain  
T +34-932 60 04 00  
F +34-932 60 23 14  
medi@mediespana.com  
www.mediespana.com

medi France  
Z.I. Charles de Gaulle  
25, rue Henri Farman  
93297 Tremblay en France Cedex  
France  
T +33-1 48 61 76 10  
F +33-1 49 63 33 05  
infos@medi-france.com  
www.medi-france.com

medi Hungary Kft.  
Bokor u. 21.  
1037 Budapest  
Hungary  
T +36 1 371-0090  
F +36 1 371-0091  
info@medi.hu  
www.medi.hu

medi Japan K.K.  
5-1 Nihonbashi Kabutocho Chuoku  
103-0026 Tokyo  
Japan  
T: +81 3 6778 2590  
F: +81 3 5847 7901  
info@medi-japan.jp  
www.medi-japan.co.jp

medi Nederland BV  
Heusing 5  
4817 ZB Breda  
The Netherlands  
T +31-76 57 22 555  
F +31-76 57 22 565  
info@medi.nl  
www.medi.nl

medi Polska Sp. z o.o.  
ul. Łąbiedzka 22  
44-121 Gliwice  
Poland  
T: +48-32 230 60 21  
F: +48-32 202 87 56  
info@medi-polska.pl  
www.medi-polska.pl

medi Bayreuth Unipessoal, Lda  
Rua do Centro Cultural, no. 43  
1700-106 Lisbon  
Portugal  
T +351-21 843 71 60  
F +351-21 847 08 33  
medi.portugal@medibayreuth.pt  
www.medi.pt

medi RUS LLC  
Business Center NEO GEO  
Butlerova Street 17  
117342 Moscow  
Russia  
T +7-495 374 04 56  
F +7-495 374 04 56  
info@medirus.ru  
www.medirus.ru

# Inhalt / Content

ENGLISH / ENGLISH	4
DEUTSCH / GERMAN	6
FRANCAIS / FRENCH	9
ITALIANO / ITALIAN	11
ESPAÑOL / SPANISH	14
PORTUGUES / PORTUGUESE	16
NEDERLANDS / DUTSCH	18
SVENSKA / SWEDISH	21
NORSK / NORWEGIAN	23
DANSK / DANISH	26
POLSKI / POLISH	28
ČEŠTINA / CZECH	30
MAGYAR / HUNGARIAN	33
УКРАЇНСЬКА / UKRAINIAN	35
ايربوع / ARABIC	39
РУССКИЙ / RUSSIAN	40
中文 / CHINESE	42
תירבע / HEBREW	47
EESTI / ESTONIAN	47
BAHASA MELAYU / MALAY	49

ENGLISH / ENGLISH

# circaid® juxtafit® essentials hand wrap

## Intended purpose

This garment is designed to provide compression to the hand of patients with venous and lymphatic disorders.

## Indications

- Lymphedema

## Contraindications

- Severe Peripheral Arterial Disease
- Uncontrolled Congestive Heart Failure
- Septic Phlebitis
- Phlegmasia Cerulea Dolens
- Untreated infection
- Any circumstance where increased venous and lymphatic return is undesirable
- Suspected or known untreated acute DVT

## Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes

and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

## Cautions

- Mild to Moderate Peripheral Arterial Disease
- Impaired Sensation - must be able to sense applied pressure
- Intolerance to material(s)
- Pediatric & assisted use - must be able to sense & communicate applied pressure

## This Package contains

- 1 circaid juxtafit essentials hand wrap
- 1 booklet instructions for use

## Storage instructions

Store in a dry place

Keep out of sunlight

## Lifespan

Due to material wear-and-tear the medical efficacy can only be expected for a defined lifespan. This assumes correct handling of the garment (e.g. proper care, application, and removal).

## circaid juxtafit essentials hand wrap

3 month - depends on the wear and tear

## Disposal

The product can be disposed of in the domestic waste.

## Material composition

Nylon 53%

Polyurethane 39,5

Elastane 7%

Other Fibres 0,5%

## Care instructions

The garment can be machine washed using a gentle cycle and dried on a low heat setting. If doing so, it is recommended that you place the garment in a pillow case to aid in keeping the hook and loop tap lint-free and to protect other clothing from getting damaged.

To extend the life of the garment, hand wash and drip dry.

Machine wash warm

Do not bleach

Tumble dry, low heat

Do not iron

Do not dry clean

Do not wring

---

## Wearing recommendation

Slightly loosen the bands of the garment for night-time wear. If experiencing any pain, or discomfort, reduce the compression (mmHg) applied or remove de-

vice. Inform your clinician of your action. When bands are placed correctly they should slightly overlap, leaving no space between bands.

Do not tuck bands. Hook tabs should never touch the skin!

## Donning instructions

Step 1: Slide the thumb through the hand wrap's thumb hole. The seam should run from the base of the thumb to the wrist. The curved portion of the wrap should be in the palm of the hand.

Step 2: Unroll the band closest to the fingers and secure it with the velcro tab to the hand wrap with firm and comfortable compression. The top edge of the hand wrap should line up with the base of the fingers.

Step 3: Wrap the remaining band around the wrist and secure the Velcro tab. Be sure not to over tighten the wrist band so that fluid can freely move out of the hand. The circaid juxtaft essentials hand wrap should overlap with the circaid juxtaft essentials arm.

**Note:** *Once all bands are secured inspect the garment for any gaps or creases. Adjust the bands as necessary.*

## Doffing instructions

Step 1: Detach the bands, starting with the wrist band, and fold them back onto themselves to keep the hook tabs lint-free.

Step 2: Detach the inner hook tab (palm) and slide the wrap of the hand.

#### **Note**

*If using multiple circaid juxtafit garments, only lightly apply bands overlapping the legging below to ensure gradient compression. Set the remaining bands to firm and comfortable compression or prescribed compression level.*

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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medi Trading (Shanghai) Co. Ltd.
medi Ukraine LLC