medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth
Germany
T +49 921 912-0
www.medi.de

medi Australia Pty Ltd
83 Fennell Street
North Parramatta NSW 2151
Australia
T +61-2 9890 8696
F +61-2 9890 8439
sales@mediaustralia.com.au
www.mediaustralia.com.au

medi Austria GmbH
Adamgasse 16/7
6020 Innsbruck
Austria
T +43 512 57 95 15
F +43 512 57 95 15 45
vertrieb@medi-austria.at
www.medi-austria.at

medi Belgium NV
Staatsbaan 77/0099
3945 Ham
Belgium
T: +32-11 24 25 60
F: +32-11 24 25 64
info@medibelgium.be
www.medibelgium.be

medi Canada Inc / médi Canada Inc
104-1375 Lionel-Boulet,
Varenes, Québec,
QC Canada J3X 1P7
T +1 450-583-3317 / +1 800 361 3153
F +1 888-583-6827
service@medicanada.ca
www.medicanada.ca

MAXIS a.s.,
medi group company
Slezská 2127/13
120 00 Prague 2
Czech Republic
T: +420 571 633 510
F: +420 571 616 271
info@maxis-medica.com
www.maxis-medica.com

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Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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

medi Trading (Shanghai) Co. Ltd.
Room 1703, T1, CES WEST BUND CENTER,
No. 277, Long Lan Road
Xuhui District
200232 Shanghai
People's Republic of China
T: +86-21 50582319
F: +86-21 50582319

Medi Turk Ortopedi Medikal İc ve Dis Tic. AS
Mustafa Kemal Mah. 2156 Sokak No: 16/4-5
06510 Çankaya Ankara
Turkey
T: +90 312 435 20 26
F: +90 312 434 17 67
info@medi-turk.com
www.medi-turk.com

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circaid® juxtalite®
ankle foot wrap

Intended purpose

The compression system is designed to provide compression to the foot for pati-

ents with venous disorders. The device consists of a series of juxtaposed, inelastic bands extending from a central portion of the wrap. The wrap contours the foot and is designed to cover from the base of the toes to just above the ankle, overlapping with the legging.

Indications

- Chronic venous insufficiency
- Venous stasis ulcer/mixed leg ulcer
- DVT/Thrombosis Prevention
- Acute deep vein thrombosis
- Post Sclerotherapy
- Post Thrombotic Syndrome
- Dependent oedema
- Varicose veins

Contraindications

- Severe Peripheral Arterial Disease
- Uncontrolled Congestive Heart Failure
- Septic Phlebitis
- Phlegmasia Cerulea Dolens
- Untreated infection
- Any circumstance where increased venous and lymphatic return is undesirable

Risks / Side Effects

- Mild to Moderate Peripheral Arterial Disease
- Impaired Sensation - must be able to sense applied pressure
- Intolerance to material(s)
- Pediatric & assisted use - must be able to sense & communicate applied pressure

Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Note on wearing

Slightly loosen the bands of the garment for night-time wear. If experiencing any pain, or discomfort, reduce the compression (mmHg) applied or remove device. Inform your clinician of your action.

When bands are placed correctly they should slightly overlap, leaving no space between bands.

Do not tuck bands. Hook tabs should never touch the skin!

Customising instructions

circaid juxtalite ankle foot wrap can be tailored to fit the patient's foot length. The ankle foot wrap should be fitted by a medical professional. Doing other than instructed will void the three month limited warranty.

See more details in the section ankle-band and foot-length adjustment.

Donning instructions

Step 1: Don optional circaid sock or undersleeve. The circaid juxtalite ankle foot wrap (afw) can also be worn directly against the skin.

Step 2: First, unroll the band closest to the toes and secure it to the opposing side to create a kind of slipper. Slide the foot into the slipper shape and firmly press the heel into the heel cup.

Step 3: The front edge of the foot piece should align with the base of the toes.

Step 4: Tighten the band closest to the toes to firm and comfortable compression.

While pulling the ankle band in the opposing direction, attach the second band with firm and comfortable compression.

Step 5: Unroll the ankle-band and pull it around the ankle with firm and comfortable compression. Ensure there are no gaps and fasten the Velcro tab. Adjust the ankle-band if any gaps or pressure points occur.

Step 6: Readjust the bands as necessary throughout the day to maintain a firm and comfortable compression level.

The garment should never hurt. It should feel firm but comfortable.

If pressure increases or decreases noticeably during wear, readjust the garment as needed.

If experiencing any pain, immediately remove the garment and consult your physician.

Removal: First detach the ankle band and fold it back onto itself to keep it lint free. Lastly, loosen the front foot band. The circaid juxtalite ankle foot wrap slips off like a shoe.

Ankle-band adjustment

The following steps will permanently reduce the length of the long ankle-band.

Step 1: Determine and mark the desired length of the ankle-band.

Step 2: Remove the afw and carefully cut it to the desired length.

Step 3: Use the included Velcro tabs to secure the ankle-band in place.

Foot-length adjustment

The following steps will permanently reduce the foot-length of the circaid juxtalite ankle foot wrap. The afw can be trimmed by up to 1.25 in. (3 cm) (see notches on both sides).

Step 1: With the beige side facing up, mark the desired length.

Step 2: Draw a straight line from where you would like to trim.

Step 3: Carefully cut along the drawn line, assuring that the reinforced seam section of the Velcro material is not cut.

Please note: No other modifications can be made to the garment.

Be careful that the reinforced seam within the Velcro material is not cut while shortening the garment

Care instructions

The garment can be machine washed using a gentle cycle and dried on a low heat setting. We recommend using a laundry bag. To extend the life of the garment, hand wash and drip dry.

-
- | | |
|--|----------------------|
|  | Machine wash warm |
|  | Do not bleach |
|  | Tumble dry, low heat |
|  | Do not iron |
|  | Do not dry clean |
|  | Do not wring |
-

Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

Lifespan:

3 months – due to material wear and/or slackening, medical efficacy can only be guaranteed for a defined lifespan. This assumes correct handling of the product (e.g. during care and when putting on and taking off the product).

Material composition

6% Nylon, 39% Polyurethane, 5% Elastane

Disposal

The product can be disposed of in the domestic waste.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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