**medi GmbH & Co. KG**  
**Medicusstraße 1**  
**95448 Bayreuth**  
**Germany**  
**T +49 921 912-0**  
**www.medi.de**

medi Australia Pty Ltd  
 Unit 4/13-21 Hallmark Street  
 Pendle Hill NSW 2145  
 Australia  
 T +61 2 9890 8696  
 F +61 2 9890 8439  
 sales@medi.au  
 www.mediaustralia.com.au

medi Austria GmbH  
 Adamgasse 16/7  
 6020 Innsbruck  
 Austria  
 T +43-512 57 95 15  
 F +43-512 57 95 15 45  
 vertrieb@medi-austria.at  
 www.medi-austria.at

medi Belgium NV  
 Staatsbaan 77/0099  
 3945 Ham  
 Belgium  
 T: +32-11 24 25 60  
 F: +32-11 24 25 64  
 info@medibelgium.be  
 www.medibelgium.be

medi Canada Inc / médi Canada Inc  
 104-1375 Lionel-Boulet,  
 Varennes, Québec,  
 QC Canada J3X 1P7  
 T +1 450-583-3317 / +1 800 361 3153  
 F +1 888-583-6827  
 service@medicanada.ca  
 www.medicanada.ca

MAXIS a.s.,  
 medi group company  
 Slezská 2127/13  
 120 00 Prague 2  
 Czech Republic  
 T: +420 571 633 510  
 F: +420 571 616 271  
 info@maxis-medica.com  
 www.maxis-medica.com

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medi Trading (Shanghai) Co. Ltd.  
 Room 613, HnR Group Building  
 No.1178, Tian Yao Qiao Road  
 Xuhui District, Shanghai, 200232  
 People's Republic of China  
 T +86-21 50582318  
 infocn@medi.cn  
 www.medi.cn

medi Sweden AB  
 Box 6034  
 192 06 Sollentuna  
 Sweden  
 T +46 8 96 97 98  
 F +46 8 626 68 70  
 info@medi.se  
 www.medi.se

Medi Turk Ortopedi Medikal İc ve Dis Tic. AS  
 Mustafa Kemal Mah. 2156 Sokak No: 16 /4-5  
 06510 Çankaya Ankara  
 Turkey  
 T: +90 312 435 20 26  
 F: +90 312 434 17 67  
 info@medi-turk.com  
 www.medi-turk.com

medi Ukraine LLC  
 Evhena Sverstiuka str, 11.  
 Kiev 02002  
 Ukraine  
 T: +380 44 591 11 63  
 F: +380 44 392 73 73  
 info@medi.ua  
 http://medi.ua

medi USA L.P.  
 6481 Franz Warner Parkway  
 Whitsett, N.C. 27377-3000  
 USA  
 T +1-336 4 49 44 40  
 F +1-888 5 70 45 54  
 info@mediusa.com  
 www.mediusa.com

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
### Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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medi Norway AS  
Vestheimvegen 35  
4250 Kopervik  
Norway  
T: +47 52 84 45 00  
F: +47 52 84 45 09  
post@medinorway.no  
www.medinorway.no

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## circaid® compression anklet

### Intended purpose

The garment is designed to provide therapeutic compression to the foot and ankle region for patients with venous and lymphatic disorders and is always to be worn in combination with lower leg compression garments and if necessary a foot garment.

### Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

**Important to note:** *Independently of medical diagnosis and prescribed compression treatment, circaid compression anklet is to*

*be worn in combination with foot, lower leg and if necessary upper leg compression garments. The circaid compression anklet is not defined as "medical compression stockings" according to RAL-GZ 387:2000-09 Medizinische Kompressionsstrümpfe. Detailed information on the indications can be found in the respective instruction manuals for the main products circaid juxtacures, circaid juxtalite or circaid juxtafit".*

### Note on wearing

The circaid compression anklet may be worn day and night, but removal is recommended if experiencing discomfort when lying down.

### Donning

If using the circaid compression anklet make sure that you have already applied a circaid leg garment and a circaid undersleeve. To avoid uncomfortable wearing due to wrinkles in the undersleeve ensure that the bottom part of the circaid undersleeve is folded over the bottom part of the circaid leg garment.

- For an easier application, slide your hand into the circaid compression anklet and turn it inside out just up to the heel.
- Slide the circaid compression anklet onto the foot and ensure that it is pulled high enough that it overlaps with the circaid leg garment.

### Care instructions

The garment can be machine washed using a gentle cycle and dried on a low heat setting. To extend the life of the garment, hand wash and drip dry.

Machine wash warm

Do not bleach

Tumble dry, low heat

Do not iron

Do not dry clean

Do not wring

---

### Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

---

This garment is not designed to be cut or modified in any way.

### Lifespan

3 months – due to material wear and/or slackening, medical efficacy can only be guaranteed for a defined lifespan. This assumes correct handling of the product (e.g. during care and when putting on and taking off the product).

### Material composition

Nylon, elastane

### Disposal

Dispose of properly after use.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufac-

turer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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