medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth
Germany
T +49 921 912-0
www.medi.de

medi Australia Pty Ltd
 83 Fennell Street
 North Parramatta NSW 2151
 Australia
 T +61 2 9890 8696
 F +61 2 9890 8439
 sales@mediaustralia.com.au
 www.mediaustralia.com.au

medi Austria GmbH
 Adamgasse 16/7
 6020 Innsbruck
 Austria
 T +43-512 57 95 15
 F +43-512 57 95 15 45
 vertrieb@medi-austria.at
 www.medi-austria.at

medi Belgium NV
 Staatsbaan 77/0099
 3945 Ham
 Belgium
 T: +32-11 24 25 60
 F: +32-11 24 25 64
 info@medi-belgium.be
 www.medi-belgium.be

medi Canada Inc / médi Canada Inc
 104-1375 Lionel-Boulet,
 Varennes, Québec,
 QC Canada J3X 1P7
 T +1 450-583-3317 / +1 800 361 3153
 F +1 888-583-6827
 service@medicanada.ca
 www.medicanada.ca

MAXIS a.s.,
 medi group company
 Slezská 2127/13
 120 00 Prague 2
 Czech Republic
 T: +420 571 633 510
 F: +420 571 616 271
 info@maxis-medica.com
 www.maxis-medica.com

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medi

circaid® compression anklet elastic compression system

Gebrauchsanweisung. Instructions for use. Mode d'emploi. Instrucciones de uso. Instruções para aplicação. Istruzioni per l'uso. Gebruiksaanwijzing. Brugsvejledning. Bruksanvisning. Návod k použití. Инструкция по использованию. Instrukcja zakładania. إرشادات الاستخدام. Használati útmutató. Інструкція з використання. הוראות שימוש. Kasutusjuhend. Käyttöohje.

circaid compression anklet

medi. I feel better.

E012043/01.2024

4 049772 365708

medi Polska Sp. z o.o.
 ul. Łąbiedzka 22
 44-121 Gliwice
 Poland
 T: +48-32 230 60 21
 F: +48-32 202 87 56
 info@medi-polska.pl
 www.medi-polska.pl

medi Bayreuth Unipessoal Lda.
 Rua do Centro Cultural, no. 43
 1700-106 Lisbon
 Portugal
 T +351 21 843 71 60
 F +351 21 847 08 33
 medi.portugal@medibayreuth.pt
 www.medi.pt

medi RUS LLC
 Business Center NEO GEO
 Butlerova Street 17
 117342 Moscow
 Russia
 T +7-495 374 04 56
 F +7-495 374 04 56
 info@medirus.ru
 www.medirus.ru

medi Trading (Shanghai) Co. Ltd.
 Room 1703, T1, CES WEST BUND
 CENTER,
 No. 277, Long Lan Road
 Xuhui District
 200232 Shanghai
 People's Republic of China
 T: +86-21 50582319
 F: +86-21 50582319

medi Sweden AB
 Box 6034
 192 06 Sollentuna
 Sweden
 T +46 8 96 97 98
 F +46 8 626 68 70
 info@medi.se
 www.medi.se


Medi Turk Ortopedi Medikal İc ve Dis Tic. AS
 Mustafa Kemal Mah. 2156 Sokak No: 16 /4-5
 06510 Çankaya Ankara
 Turkey
 T: +90 312 435 20 26
 F: +90 312 434 17 67
 info@medi-turk.com
 www.medi-turk.com

medi Ukraine LLC
 Evhena Sverstiuka str, 11.
 Kiev 02002
 Ukraine
 T: +380 44 591 11 63
 F: +380 44 392 73 73
 info@medi.ua
 http://medi.ua

medi USA L.P.
 6481 Franz Warner Parkway
 Whitsett, N.C. 27377-3000
 USA
 T +1-336 4 49 44 40
 F +1-888 5 70 45 54
 info@mediusa.com
 www.mediusa.com

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
Important notes

This medical device **MD** is made only for multiple uses on a single patient.  If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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medi Norway AS
Vestheimvegen 35
4250 Korpervik
Norway
T: +47 52 84 45 00
F: +47 52 84 45 09
post@medinorway.no
www.medinorway.no

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circaid® compression anklet

Intended purpose

The garment is designed to provide therapeutic compression to the foot and ankle region for patients with venous and lymphatic disorders and is always to be worn in combination with lower leg compression garments and if necessary a foot garment.

Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Important to note: *Independently of medical diagnosis and prescribed compression treatment, circaid compression anklet is to be worn in combination with foot, lower leg and if necessary upper leg compression garments. The circaid compression anklet is not defined as "medical compression stockings" according to RAL-GZ 387:2000-09 Medizinische Kompressionsstrümpfe. Detailed information on the indications can be found in the respective instruction manuals for the main products circaid juxtacures, circaid juxtalite or circaid juxtafit".*

Note on wearing

The circaid compression anklet may be worn day and night, but removal is recommended if experiencing discomfort when lying down.

Donning

If using the circaid compression anklet make sure that you have already applied a circaid leg garment and a circaid undersleeve. To avoid uncomfortable wearing due to wrinkles in the undersleeve ensure that the bottom part of the circaid undersleeve is folded over the bottom part of the circaid leg garment.

- For an easier application, slide your hand into the circaid compression anklet and turn it inside out just up to the heel.
- Slide the circaid compression anklet onto the foot and ensure that it is pulled high enough that it overlaps with the circaid leg garment.

Care instructions

The garment can be machine washed

using a gentle cycle and dried on a low heat setting. To extend the life of the garment, hand wash and drip dry.

	Machine wash warm
	Do not bleach
	Tumble dry, low heat
	Do not iron
	Do not dry clean
	Do not wring

Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

This garment is not designed to be cut or modified in any way.

Lifespan

<3 month - depends on the wear and tear
Due to material wear-and-tear and the medical efficacy can only be guaranteed for a defined lifespan. This assumes correct handling of the garment (e.g. proper care, donning, and doffing).

Material composition

Nylon, elastane

Disposal

The product can be disposed of in the domestic waste.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your

specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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