**medi GmbH & Co. KG**  
Medicusstraße 1  
95448 Bayreuth  
Germany  
T +49 921 912-0  
www.medi.de

medi Australia Pty Ltd  
83 Fennell Street  
North Parramatta NSW 2151  
Australia  
T +61-2 9890 8696  
F +61-2 9890 8439  
sales@mediaustralia.com.au  
www.mediaustralia.com.au

medi Austria GmbH  
Adamgasse 16/7  
6020 Innsbruck  
Austria  
T +43 512 57 95 15  
F +43 512 57 95 15 45  
vertrieb@medi-austria.at  
www.medi-austria.at

medi Belgium NV  
Staatsbaan 77/0099  
3945 Ham  
Belgium  
T: +32-11 24 25 60  
F: +32-11 24 25 64  
info@medibelgium.be  
www.medibelgium.be

medi Canada Inc / médi Canada Inc  
104-1375 Lionel-Boulet,  
Varenes, Québec,  
QC Canada J3X 1P7  
T +1 450-583-3317 / +1 800 361 3153  
F +1 888-583-6827  
service@medicanada.ca  
www.medicanada.ca

MAXIS a.s.,  
medi group company  
Slezská 2127/13  
120 00 Prague 2  
Czech Republic  
T: +420 571 633 510  
F: +420 571 616 271  
info@maxis-medica.com  
www.maxis-medica.com

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medi Trading (Shanghai) Co. Ltd.  
Room 1703, T1, CES WEST BUND CENTER,  
No. 277, Long Lan Road  
Xuhui District  
200232 Shanghai  
People's Republic of China  
T: +86-21 50582319  
F: +86-21 50582319

UK Responsible Person  
medi UK Ltd.  
Plough Lane  
Hereford HR4 0EL  
Great Britain  
T +44-1432 37 35 00  
F +44-1432 37 35 10  
enquiries@mediuk.co.uk  
www.mediuk.co.uk

**medi**

## circaid® juxtafit® premium inelastic compression system

Instructions for use / Gebrauchsanweisung / Instructions d'utilisation / Istruzioni per l'uso / Instrucciones de uso / Instruções de uso / Gebruiksaanwijzing / Användningsinstruktioner / Brugsanvisning / دليل الاستخدام / Instrukcja użytkowania / Инструкция по использованию / הוראות שימוש / Kasutusjuhend / Arahan penggunaan. Pukeutumisohte.




circaid juxtafit premium ankle foot wrap

circaid juxtafit premium interlocking ankle foot wrap

circaid juxtafit premium lower leg

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### Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If undue pain or an unpleasant sensation develops while you are wearing the product, please consult your doctor or orthotist immediately. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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## circaid® juxtafit® premium

### Intended purpose

The compression system is designed to provide compression to the leg for patients with venous and lymphatic disorders. The device consists of a series of juxtapositioned, inelastic bands extending from the central portion of the garment. The wrap contours to the leg and is designed to cover some or all of the leg, depending on the custom design.

### Indications

- Lymphedema
- other forms of edema
- Lipedema
- Chronic venous insufficiency
- Varicose veins
- DVT/Thrombosis Prevention
- Acute deep vein thrombosis
- Post Thrombotic Syndrome
- Venous Stasis Ulcer
- Post Sclerotherapy
- Lipodermatosclerosis
- 

### Contraindications

- Severe Peripheral Arterial Disease
- Uncontrolled Congestive Heart Failure
- Septic Phlebitis
- Phlegmasia Cerulea Dolens
- Untreated infection
- Any circumstance where increased venous and lymphatic return is undesirable

### Risks / Side Effects

- Mild to Moderate Peripheral Arterial Disease
- Impaired Sensation - must be able to sense applied pressure
- Intolerance to material(s)
- Pediatric & assisted use - must be able to sense & communicate applied pressure

### Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

### Wearing recommendation

**Note:** *Unless your treating physician orders or recommends otherwise, we recommend wearing the leg fitting together with the respective foot options enclosed with the set. In case of doubt, please consult your attending physician.*

Always ensure that the hook tabs are placed on the corresponding landing pads. Always ensure that the correct prescribed pressure range is being applied. Slightly loosen the bands of the garment for night-time wear. Each vertical BPS line on the bands may be at a different location on the

leg depending on the size and shape of the leg. They do not need to be aligned down the front center of the leg or in any other specific location on the leg. If experiencing any pain, or discomfort, reduce the compression (mmHg) applied or remove device. Inform your clinician of your action. When bands are placed correctly they should slightly overlap, leaving no space between bands. Do not tuck bands. The Hook tabs should never touch the skin!

### Donning instructions

Step 1: Slide the enclosed circaid undersleeve or undersock onto the leg up to the knee (for lower leggings) or to the groin (for upper leg and whole leg).

Step 2: Position the legging with the black side of the material facing the skin and the seam running along the back of the leg.

Step 3: Unroll the top two bands and loosely secure them to the outside of the legging to anchor it in place. This enables you to easily raise or lower the legging for correct positioning

**Note:** *If you are using your legging with a circaid foot piece, e.g. pac band or juxtafit ankle foot wraps, follow the application instructions included with your foot piece.*

*Do not secure the ankle band of the foot piece at this time as it be applied later to overlap the legging.*

*If using the circaid single band ankle foot wrap, wait until after the legging has been applied to apply the wrap.*

Step 4: Unroll the bottom two bands and secure the bottommost band to the outside of the legging to a firm and comfortable compression level.

The bottom most band can be identified by the care label.

Step 5: While holding the second band, detach the next band. Secure the second band and continue this process with the other bands, going up the leg. Ensure that the bands are placed in an alternating order.

Step 6: The legging should now lay flat and wrinkle-free against the leg.

Adjust each band to a firm and comfortable compression. See section "Using the Built-In-Pressure System" for instructions on using the BPS card to set the compression to a prescribed level.

Step 7: Once all the bands are secured, inspect the garment for any gaps or creases. Adjust the bands as necessary. You may now fold the bottom and top parts of the circaid undersleeve or sock over the top and bottom edge of the legging if desired.

Step 8 (optional): Finish applying foot compression. At this time, finish applying your foot compression as directed in the separate instructions. The ankle band will overlap with the lower portion of the juxtafit legging.

Step 9 (optional): Apply shelf straps, if included with your custom compression system. Use the shelf strap to provide additional support or coverage to areas in need, pulling up and around the area on the limb until firm and comfortable compression is achieved.

Step 10 (optional): Apply knee cover, if your compression system covers the knee. Place the knee cover over the knee cap and pull the individual tabs at each corner to firm and comfortable compression. Secure the tabs to the legging once the knee cover is securely against the knee.

**Note:** *If using a circaid cover up, apply it at this time. Slide the*

cover up over the legging and secure it using the hook tabs along the top edge. Refer to the instructions for use included with the cover up for additional information.

#### Doffing instructions

To remove the legging, first remove the cover up, if using. Then remove the knee cover, shelf straps, and/or foot compression, if using any of these. Then unfold the circaid undersleeve.

Detach all bands of the legging starting with the top band. Fold each band onto itself to keep the hook and loop tabs lint-free.

Remove the legging and then remove the circaid undersleeve or sock.

#### Using the Built-In-Pressure System (BPS)

Step 1: Locate the Built-In-Pressure card in your packaging.

Step 2: Identify your prescribed pressure scale using the color-coded system.

Step 3: Starting with the bottom band, line up the black triangle on the card with one of the BPS lines on the bottom band.

Step 4: Note where the second BPS line lines up with the card's compression ranges (A).

Step 5: If the BPS line on the garment either falls short (B) or goes beyond the correct compression range (C), readjust the band as necessary so that the second BPS line on the garment is aligned with the prescribed compression range (A). Completely loosen the band before reapplying it if too much tension has been applied.

Step 6: Repeat the steps 3 - 5 for each band going up the leg. Be sure to adjust all bands to the same tension to ensure gradient compression.

Bands may need adjusting throughout the day to maintain the prescribed compression range as edema is reduced.

**Note:** Please note it is recommended that the limb is always treated in its entirety with medical compression.

*This means that compression products for the leg begin at the foot or the ankle and extend up to below the hollow of the knee or beyond to include the thigh.*

#### Care instructions

The garment can be machine washed using a gentle cycle and dried on a low heat setting. We recommend using a laundry bag. To extend the life of the garment, hand wash and drip dry.

 Machine wash warm

 Do not bleach

 Tumble dry, low heat

 Do not iron

 Do not dry clean

 Do not wring

#### Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

#### Lifespan

12 months - due to material wear-and-tear the medical efficacy can only be expected for a defined lifespan. This assumes correct handling of the garment (e.g. proper care, application, and removal).

circaid juxtafit premium shelf strap:

<3 month - depends on the wear and tear

#### Material composition

- 49 % Polyurethane
- 46 % Nylon
- 5 % Elastane

#### Disposal

The product can be disposed of in the domestic waste

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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