medi GmbH & Co. KG
Medicusstraße 1
D-95448 Bayreuth
Germany
T +49 921 912-0
www.medi.de

medi Australia Pty Ltd
 83 Fennell Street
 North Parramatta NSW 2151
 Australia
 T +61 2 9890 8696
 F +61 2 9890 8439
 sales@mediaustralia.com.au
 www.mediaustralia.com.au

medi Austria GmbH
 Adamgasse 16/7
 6020 Innsbruck
 Austria
 T +43-512 57 95 15
 F +43-512 57 95 15 45
 vertrieb@medi-austria.at
 www.medi-austria.at

medi Belgium NV
 Posthoornstraat 13/1
 3582 Koersel
 Belgium
 T: + 32-11 24 25 60
 F: +32-11 24 25 64
 info@medibelgium.be
 www.medibelgium.be

medi Brasil
 Rua Cristovan de Vita, n.º 260, Un. 30,
 Centro Logístico Raposo Tavares,
 Vila Camargo
 Vargem Grande Paulista - SP
 CEP 06730-000
 Brasil
 T: +55-11-3500 8005
 sac@medibrasil.com
 www.medibrasil.com

medi Canada Inc / médi Canada Inc
 597, Rue Duvernay, Verchères
 QC Canada J0L 2R0
 T: +1 450-583-3317 / +1 800 361 3153
 F: +1 888-583-6827
 service@medicanada.ca
 www.medicanada.ca

MAXIS a.s.,
 medi group company
 Slezská 2127/13
 120 00 Prague 2
 Czech Republic
 T: +420 571 633 510
 F: +420 571 616 271
 info@maxis-medica.com
 www.maxis-medica.com

medi Danmark ApS
 Vejlegårdsvej 59
 2665 Vallensbaek Strand
 Denmark
 T +45 70 25 56 10
 F +45 70 25 56 20
 kundeservice@sw.dk
 www.medidanmark.dk

medi

circaid® customizable hand wrap inelastic compression system

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

circaid® customizable hand wrap

Purpose

The compression system is designed to provide compression to the hand for patients with venous and lymphatic disorders. The wrap consists of a hand piece that can be individually adjusted to the patient's size and needs as edema reduces. The hand component covers from the base of the fingers to the lower arm.

Indications

- Lymphedema
- other forms of edema

Contraindications

- Severe Peripheral Arterial Disease
- Uncontrolled Congestive Heart Failure
- Septic Phlebitis
- Phlegmasia Cerulea Dolens
- Untreated infection
- Any circumstance where increased venous and lymphatic

return is undesirable

- Suspected or known untreated acute DVT

Cautions

- Mild to Moderate Peripheral Arterial Disease
- Impaired Sensation - must be able to sense applied pressure
- Intolerance to material(s)
- Pediatric & assisted use - must be able to sense & communicate applied pressure

Storage instructions

Store in a dry place

Keep out of sunlight

Lifespan

Due to material wear-and-tear and the medical efficacy can only be guaranteed for a defined lifespan. This assumes correct handling of the garment (e.g. proper care, donning, and doffing).

circaid customizable hand wrap:
2 months

medi Bayreuth España SL
C/Canigó 2 – 6 bajos
Hospitalet de Llobregat
08901 Barcelona
Spain
T +34 93 260 04 00
F +34 93 260 23 14
medi@mediespana.com
www.mediespana.com

medi France
Z.I. Charles de Gaulle
25, rue Henri Farman
93297 Tremblay en France Cedex
France
T +33 1 48 61 76 10
F +33 1 49 63 33 05
infos@medi-france.com
www.medi-france.com

medi Hungary Kft.
Bokor u. 21.
1037 Budapest
Hungary
T +36 1371 0090
F +36 1371 0091
info@medi.hu
www.medi.hu

medi Italia s.r.l.
Via Giuseppe Ghedini, 2
40069 Zola Predosa (BO)
Italy
T: +39-051 6 13 24 84
F: +39-051 6 13 29 56
info@medi-italia.it
www.medi-italia.it

medi Japan K.K.
5-1 Nihonbashi Kabutocho Chuoku
103-0026 Tokyo
Japan
T: +81 3 6778 2590
F: +81 3 5847 7901
info@medi-japan.jp
www.medi-japan.co.jp

Tachezy Sanit Srl.
Via Druso 48/a
39100 Bolzano
Italy
T: +39 0471 288 068
F: +39 0471 282 773
info@tachezysanit.com
www.tachezysanit.com

medi Nederland BV
Heusing 5
4817 ZB Breda
The Netherlands
T +31 76 57 22 555
F +31 76 57 22 565
info@medi.nl
www.medi.nl

medi Norway AS
Vestheimvegen 35
4250 Kopervik
Norway
T: +47 52 84 45 00
F: +47 52 84 45 09
post@medinorway.no
www.medinorway.no

The hand wrap is not designed to be cut or modified in any manner other than that outlined in the directions for use. Doing so in any way other than instructed will void the warranty.

Disposal

Please dispose of the devices in household waste. No special criteria for disposal apply for medical compression garments.

Material composition

circaid compression garment
Nylon 74%
Polyurethane 20%
Elastane 6%

This garment does not contain latex.

Care instructions

The garment can be machine washed using a gentle cycle and

dried on a low heat setting. If doing so, it is recommended that you place the garment in a pillow case to aid in keeping the hook and loop tap lint-free and to protect other clothing from getting damaged. To extend the life of the garment, hand wash and drip dry.

 Machine wash warm - gentle cycle

 Do not bleach

 Tumble dry, low heat

 Do not iron

 Do not dry clean

 Do not wring

Wearing recommendations

Slightly loosen the bands of the garment for night-time wear. If experiencing any pain, or discomfort, reduce the compression (mmHg) applied or remove device. Inform your clinician of your action.

When bands are placed correctly they should slightly overlap, leaving no space between bands. Do not tuck bands. Hook tabs should never touch the skin!

Note:

The hand wrap can be worn on the left or right hand.

Apply an appropriate finger compression, if needed.

Apply a circaid undersleeve for added comfort and exudate, if desired.

Donning/Customizing Instructions

Step 1:

Measure the circumference around the mid-palm in centimetre. (1)

Step 2:

Identify the corresponding measurement mark on the hand wrap. Cut along the indicated white cut lines. (2)

Step 3:

Identify which hand you want to put the hand wrap on, then identify the corresponding side of the hand wrap. Apply the small piece of black hook tab to the palm portion of the hand wrap, placing it securely inside the indicated rectangle outline.

- If the patient is wearing the hand wrap on the left hand, apply the hook tab to the side that says “LEFT PALM.” (3a)
 - If the patient is wearing the hand wrap on the right hand, apply the hook tab to the side that says “RIGHT PALM.” (3b)
-

Note:

Apply the black hook tab so that the side with longer, stronger, more pronounced hooks and rougher feel is applied directly to the hand wrap. The smooth side with less pronounced hooks should face outwards. (3c-d) This ensures that the hook tab will stay firmly secured to the hand wrap throughout wearing.

Step 4:

Instruct the patient to slide the thumb through the thumbhole. The black hook tab should not touch the skin, and you should be able to read “RIGHT SIDE FACE OUT” or “LEFT SIDE FACE OUT”, depending on which hand the hand wrap is applied worn on. (4a-b)

Step 5:

With the free hand, the patient pull the body of the hand wrap towards the palm and secure it firmly but comfortably to the black hook tab. (5)

Step 6:

Instruct the patient to wrap the band closest to the fingers around the hand and secure it with one of the included hook tabs. (6) The band can be trimmed in length if desired.

Step 7:

Instruct the patient to wrap the band closest to the wrist around the hand and wrist and secure it with one of the included hook tabs. (7) The band can be trimmed in length if desired.

The circaid customizable hand wrap has now been applied. (8a-b)

Note:

If the wrist band does not contour comfortably to the base of the hand, additional material can be trimmed from the middle of the hand wrap (see Step 2). You will achieve an optimal fit if the body of the hand wrap covers the hook

tab with just a small amount of excess fabric (9). Never trim so much material that the body of the hand wrap does not cover the hook tab.

Doffing instructions

Instruct the patient to detach all bands starting at the wrist band. Fold each band back onto itself to keep the hook material lint free. Gently unfasten the material hooked together at the palm and remove the hand wrap completely. Leave all hook tabs in place to ease the next application.

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medi UK Ltd.
medi Ukraine LLC