**medi GmbH & Co. KG**  
**Medicusstraße 1**  
**D-95448 Bayreuth**  
**Germany**  
**T +49 921 912-0**  
**F +49 921 912-780**  
**ortho@medi.de**  
**www.medi.de**

medi Australia Pty Ltd  
 Unit 4/13-21 Hallmark Street  
 Pendle Hill NSW 2145  
 Australia  
 T +61 2 9890 8696  
 F +61 2 9890 8439  
 sales@medi.au  
 www.mediaustralia.com.au

medi Austria GmbH  
 Adamgasse 16/7  
 6020 Innsbruck  
 Austria  
 T +43 512 57 95 15  
 F +43 512 57 95 15 45  
 vertrieb@medi-austria.at  
 www.medi-austria.at

medi Belgium NV  
 Staatsbaan 77/0099  
 3945 Ham  
 Belgium  
 T: +32-11 24 25 60  
 F: +32-11 24 25 64  
 info@medibelgium.be  
 www.medibelgium.be

medi Canada Inc / médi Canada Inc  
 104-1375 Lionel-Boulet,  
 Varennes, Québec,  
 QC Canada J3X 1P7  
 T +1 450-583-3317 / +1 800 361 3153  
 F +1 888-583-6827  
 service@medicanada.ca  
 www.medicanada.ca

MAXIS a.s.,  
 medi group company  
 Slezská 2127/13  
 120 00 Prague 2  
 Czech Republic  
 T: +420 571 633 510  
 F: +420 571 616 271  
 info@maxis-medica.com  
 www.maxis-medica.com

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
medi Trading (Shanghai) Co. Ltd.  
 Room 1703, T1, CES WEST BUND CENTER,  
 No. 277, Long Lan Road  
 Xuhui District  
 200232 Shanghai  
 People's Republic of China  
 T: +86-21 50582319  
 F: +86-21 50582319

UK Responsible Person  
 medi UK Ltd.  
 Plough Lane  
 Hereford HR4 0EL  
 Great Britain  
 T +44-1432 37 35 00  
 F +44-1432 37 35 10  
 enquiries@mediuk.co.uk  
 www.mediuk.co.uk

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#### Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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Aluminium, Velours, Polyester, PU, Polyamid

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# M.4s<sup>®</sup> comfort

## Intended Purpose

M.4s comfort is a brace with a frame that is intended to reinforce and stabilise  the knee joint using limitation of extension and flexion.

## Indications

All indications that require physiological stabilisation and reinforcement of the knee joint and/or stress relief of the joint in at least two planes, such as:

- Following rupture of the cruciate ligament (conservative / post-operative)
- Following rupture of the collateral ligament (conservative / post-operative)
- For knee instability
- Following surgical intervention on the meniscus
- Following cartilage surgery

## Contraindications

None known at present.

## Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in conjunction with sweating) or due to the composition of the material. The immobilisation which is achievable by the device presents a risk for development of thrombosis in the lower veins of the legs, especially in the context of surgical interventions or trauma.

## Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

## Setting the flexion / extension limits

Extension limit wedge, gold:

5° Hypext. , 0°, 10°, 20°, 30°, 45°

Flexion limit wedge, silver:

0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°, 120°

Can be immobilised at:

0°, 10°, 20°, 30°, 45°

For immobilising the orthosis first insert the extension wedge and then the corresponding flexion wedge. The medial (inside) and the lateral (outside) wedge settings must always be the same.

- 1 0° extension and 90° flexion wedges are pre-installed in your knee brace as standard (In the UK only 10° extension wedges fitted, remainder included as spare parts).
- 2 In order to change the wedges, press the joint cover on the rear side (see profile) outwards.
- 3 You can change the wedges by removing the tweezers from the joint cover and pulling the wedges forwards (extension) or backwards (flexion) out of the joint using the tweezers. The wedges are accordingly labelled with „FLEX“ or „EXT“ and the corresponding degree. The wedges may also be distinguished by the colour. Extension wedges are golden, flexion wedges are silver. Push the new wedge from the front side (extension) or from the back (flexion) into the joint, so the markings „UP“ on the joint and the wedge are matched accordingly. Upon proper locking you will hear a click sound.
4. Move the joints of the orthosis to check the

position of the wedges. Put the tweezers back into the mounting and attach the joint covers by hooking the two noses of the joint cover into the front side of the joint and then pushing the cover onto the joint. The cover is firmly attached if you hear a click sound.

### Important notes

Only change the extension and flexion limits on instruction from your doctor or orthotist. In order to prevent hyperextension, always insert extension wedges or 0° wedges. The wedges define the exact excursion of the orthosis. In difficult soft-tissue conditions, e.g. where the soft-tissues are under compression, it may be necessary to insert the larger wedge to limit the joint movement to the desired excursion. In addition, make sure that the orthosis joints are correctly positioned: the middles of the joints should be level with the upper border of the patella. The joint must lie behind the midline of the leg.

### Shaping the orthosis

When shaping the orthosis, please ensure that the joints stay parallel. This guarantees that the orthosis will function well and run easily and that there is no excessive wear.

### Fastening the belt eyes

The belt eyes have been adequately fastened and are, to some extent, deliberately mobile. Please do not tighten the screws; they are tight enough. If you apply too much force, you can break off the screws.

### Fitting instructions

- Open all the straps. Sit on the edge of a chair, pull the brace up the affected leg and flex your knee to 45°.
- Position the joints so that the centre of the joint (extension wedges' height) is positioned at the level of the top edge of the kneecap and slightly behind the leg's side central line. Please make sure that the brace has been placed on the leg straight and is not twisted. Fasten Belt No. 1 first and then Belt No. 2.
- Keep your leg bent and check again that you have put the brace on properly.

- First fasten the lower front calf strap (strap #1), then the upper back thigh strap (#2). Keep your knee in a flexed position and push both joints towards the back of the leg.
- Now fasten the lower back thigh strap (#3).
- Now fasten the upper front calf strap (#4). Then the upper back calf strap (#5), and finally the front thigh strap (#6).
- Now move the knee lightly a few times and check again that you have put the brace on properly. Make any necessary corrections.

### Care instructions

The brace cannot be worn in saltwater. If the brace is worn in chlorinated water, it should be rinsed out thoroughly with clean water afterwards. Soap residues, lotions and ointments can cause skin irritation and material wear.

- Wash the product by hand, preferably using medi clean washing agent.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

With normal use, the hinge mechanism should not require any maintenance, but it can be lubricated with PTFE spray. Check regularly that the joint screws are secure, tighten them if necessary.

### Storage instructions

Keep the brace in a dry place and do not expose to direct sunlight.

### What can I do if the brace slips?

- Check that the brace is fitted properly. The centre of the brace joints should be level with the upper border of the kneecap.
- Make sure that the back thigh strap (#5) is above the belly of the calf muscle
- Take the covers off the knee pads and wear the pads directly on your skin. Watch for skin reactions; skin irritation may occur in people

- with very sensitive skin.
- Do not wear the brace over clothes.

**Material composition**

Aluminum, velours, polyester, PU, polyamide

**Liability**

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

**Disposal**

Dispose of properly after use.

**Your medi Team**

Wishes you a speedy recovery!

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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