**medi GmbH & Co. KG**  
**Medicusstraße 1**  
**D-95448 Bayreuth**  
**Germany**  
**T +49 921 912-0**  
**www.medi.de**

medi Australia Pty Ltd  
 83 Fennell Street  
 North Parramatta NSW 2151  
 Australia  
 T +61 2 9890 8696  
 F +61 2 9890 8439  
 sales@mediaustralia.com.au  
 www.mediaustralia.com.au

medi Austria GmbH  
 Adamgasse 16/7  
 6020 Innsbruck  
 Austria  
 T +43-512 57 95 15  
 F +43-512 57 95 15 45  
 vertrieb@medi-austria.at  
 www.medi-austria.at

medi Belgium NV  
 Posthoornstraat 13/1  
 3582 Koersel  
 Belgium  
 T: +32-11 24 25 60  
 F: +32-11 24 25 64  
 info@medibelgium.be  
 www.medibelgium.be

medi Brasil  
 Rua Cristovan de Vita, n.º 260, Un. 30,  
 Centro Logístico Raposo Tavares,  
 Vila Camargo  
 Vargem Grande Paulista - SP  
 CEP 06730-000  
 Brasil  
 T: +55-11-3500 8005  
 sac@medibrasil.com  
 www.medibrasil.com

medi Canada Inc / médi Canada Inc  
 597, Rue Duvernay, Verchères  
 QC Canada J0L 2R0  
 T: +1 450-583-3317 / +1 800 361 3153  
 F: +1 888-583-6827  
 service@medicanada.ca  
 www.medicanada.ca

MAXIS a.s.,  
 medi group company  
 Slezská 2127/13  
 120 00 Prague 2  
 Czech Republic  
 T: +420 571 633 510  
 F: +420 571 616 271  
 info@maxis-medica.com  
 www.maxis-medica.com

medi Danmark ApS  
 Vejlegardsvej 59  
 2665 Vallensbaek Strand  
 Denmark  
 T +45 70 25 56 10  
 F +45 70 25 56 20  
 kundeservice@sw.dk  
 www.medit danmark.dk

medi Bayreuth España SL  
 C/Canigó 2 – 6 bajos  
 Hospitalet de Llobregat  
 08901 Barcelona  
 Spain  
 T +34 93 260 04 00  
 F +34 93 260 23 14  
 medi@mediespana.com  
 www.mediespana.com

medi France  
 Z.I. Charles de Gaulle  
 25, rue Henri Farman  
 93297 Tremblay en France Cedex  
 France  
 T +33 1 48 61 76 10  
 F +33 1 49 63 33 05  
 infos@medi-france.com  
 www.medi-france.com

medi Hungary Kft.  
 Bokor u. 21.  
 1037 Budapest  
 Hungary  
 T +36 1371 0090  
 F +36 1371 0091  
 info@medi.hu  
 www.medi.hu

medi Italia s.r.l.  
 Via Giuseppe Ghedini, 2  
 40069 Zola Predosa (BO)  
 Italy  
 T: +39-051 6 13 24 84  
 F: +39-051 6 13 29 56  
 info@medi-italia.it  
 www.medi-italia.it

medi Japan K.K.  
 5-1 Nihonbashi Kabutocho Chuoku  
 103-0026 Tokyo  
 Japan  
 T: +81 3 6778 2590  
 F: +81 3 5847 7901  
 info@medi-japan.jp  
 www.medi-japan.co.jp

Tachezy Sanit Srl.  
 Via Druso 48/a  
 39100 Bolzano  
 Italy  
 T: +39 0471 288 068  
 F: +39 0471 282 773  
 info@tachezysanit.com  
 www.tachezysanit.com

medi Nederland BV  
 Heusing 5  
 4817 ZB Breda  
 The Netherlands  
 T +31 76 57 22 555  
 F +31 76 57 22 565  
 info@medi.nl  
 www.medi.nl

medi Norway AS  
 Vestheimvegen 35  
 4250 Kopervik  
 Norway  
 T: +47 52 84 45 00  
 F: +47 52 84 45 09  
 post@medinorway.no  
 www.medinorway.no

medi Polska Sp. z o.o.  
 ul. Łąbędzka 22  
 44-121 Gliwice  
 Poland  
 T: +48-32 230 60 21  
 F: +48-32 202 87 56  
 info@medi-polska.pl  
 www.medi-polska.pl

medi Bayreuth Unipessoal Lda.  
 Rua do Centro Cultural, no. 43  
 1700-106 Lisbon  
 Portugal  
 T +351 21 843 71 60  
 F +351 21 847 08 33  
 medi.portugal@medibayreuth.pt  
 www.medi.pt

**medi**

## circaid® reduction kit™ knee inelastic compression system

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




## circaid® reduction kit knee

### Purpose

The compression system is designed to provide compression to the leg for patients with venous and lymphatic disorders.

The device consists of a series of juxtapositioned, inelastic bands extending from the central portion of the garment. The wrap contours the leg and is designed to cover from just below the knee to the mid of upper leg.

### Indications

- Lymphedema
- other forms of edema
- Lipedema
- Chronic venous insufficiency
- Varicose veins
- DVT/Thrombosis Prevention
- Post Thrombotic Syndrome
- Venous Stasis Ulcer
- Post Sclerotherapy
- Lipodermatosclerosis

### Contraindications

- Severe Peripheral Arterial Disease
- Uncontrolled Congestive Heart Failure
- Septic Phlebitis
- Phlegmasia Cerulea Dolens
- Untreated infection
- Any circumstance where increased venous and lymphatic return is undesirable
- Suspected or known untreated acute DVT

### Cautions

- Mild to Moderate Peripheral Arterial Disease
- Impaired Sensation - must be able to sense applied pressure
- Intolerance to material(s)

- Pediatric & assisted use - must be able to sense & communicate applied pressure

### Washing instructions

- Machine wash warm - gentle cycle
- Do not bleach
- Tumble dry, low heat
- Do not iron
- Do not dry clean
- Do not wring

When washing it is recommended that the circaid reduction kit component be placed in a pillow case to keep the hook material lint-free and to prevent possible damage to other garments. Should lint collect on the tabs, it can be removed with tweezers.

### Storage instructions

- Store in a dry place
- Keep out of sunlight

### Lifespan

circaid reduction kit has a lifespan of 2 months.

Due to material wear-and-tear the medical efficacy can only be expected for a defined lifespan. This assumes correct handling of the garment (e.g. proper care, application, and removal).

### Disposal

Please dispose of the garment in household waste. No special criteria for disposal apply for medical compression garments.

medi RUS LLC  
Business Center NEO GEO  
Butlerova Street 17  
117342 Moscow  
Russia  
T +7-495 374 04 56  
F +7-495 374 04 56  
info@medirus.ru  
www.medirus.ru

medi Trading (Shanghai) Co. Ltd.  
Room 1703, T1, CES WEST BUND CENTER,  
No. 277, Long Lan Road  
Xuhui District  
200232 Shanghai  
People's Republic of China  
T: +86-21 50582319  
F: +86-21 50582319

medi Sweden AB  
Box 6034  
192 06 Sollentuna  
Sweden  
T +46 8 96 97 98  
F +46 8 626 68 70  
info@medi.se  
www.medi.se

Medi Turk Ortopedi Medikal Ic ve Dis Tic. AS  
Mustafa Kemal Mah. 2156 Sokak No: 16 /4-5  
06510 Çankaya Ankara  
Turkey  
T: +90 312 435 20 26  
F: +90 312 434 17 67  
info@medi-turk.com  
www.medi-turk.com

medi UK Ltd.  
Plough Lane  
Hereford HR4 OEL  
Great Britain  
T +44 1432 37 35 00  
F +44 1432 37 35 10  
enquiries@mediuk.co.uk  
www.mediuk.co.uk

medi Ukraine LLC  
Evhena Sverstiuka str, 11.  
Kiev 02002  
Ukraine  
T: +380 44 591 11 63  
F: +380 44 392 73 73  
info@medi.ua  
http://medi.ua

medi USA L.P.  
6481 Franz Warner Parkway  
Whitsett, N.C. 27377-3000  
USA  
T +1-336 4 49 44 40  
F +1-888 5 70 45 54  
info@mediusa.com  
www.mediusa.com

**Note:**

Undersleeves are provided with the lower leg components only or are available for purchase separately.

**Material composition**

circaid compression garment

56% Polyurethane

40% Nylon

4% Elastane

For all other components, please refer to the inserts included with each item.

This garment does not contain latex.

**Individualisation****Step 1:**

This knee component consists of three pieces: two “body” portions with printed numbers along the black interior and a “spine” portion that will be used to secure the two body portions together. The additional spine is a spare part for use when needed. Lay the body portions out flat so that the number markings on the black side of the material face up.

**Step 2:**

Measure the patient’s circumference at the largest part of the calf C, at the patella E, and 15 cm above the patella E2, in centimeters.

**Step 3:**

Locate the C measurement within the number sequences along the bottom edges of both body portions. Locate the E measurement within the number sequences along the center of both body portions. Locate the E2 measurement within the number sequences long the top edges of both body portions.

**Step 4:**

Using a white marker, or visible alternative, connect the plotted points

to create a smooth curve. Each patient’s leg measurements will produce a different spine curve. There is no one correct spine shape.

**Step 5:**

Cut along the plotted line on both portions.

**Tip:**

Cut one of the body portions first, lay the cut piece on top of the other body piece so black surfaces face each other, and use the cut edge to trace an identical curve onto the other body portion. Cut the second piece.

**Step 6:**

Flip the two pieces so that the beige side of the material is visible.

**Step 7:**

Starting at the top, push the cut edges of the body portions together without overlapping the material. While pushing together the edges, center the spine with the seam and secure it to both body portions.

**Tip:**

- Secure only a small section of the spine at a time
- While attaching the spine, curl the remainder of the spine back to keep it out of the way and ease the securing process

**Step 8:**

Once the two body portions have been secured together, trim off the excess spine.

---

**Tip:**

Once the patient reduces significantly, all necessary circumference measurements can be retaken and the knee component can be adjusted to the patient's new size by cutting and resetting the spine following steps 1- 7. The additional spine may be used for these adjustments if desired.

---

**Step 9:**

The undersleeve is worn for added comfort and when exudate is expected. Instruct the patient to slide the undersleeve onto the leg, covering the knee area. The undersleeve may be trimmed in length if desired, keeping in mind that if this is a whole leg application, then it should remain intact until the last circaid reduction kit component is customized.

**Step 10:**

To apply the circaid reduction kit knee, the patient should position it around the leg with the black side of the material facing the skin and with the center of the knee component, identified by the middle number sequences, in the knee crease. Instruct the patient to unroll the top two bands, to pull them above and below each other across the front of the leg, and to loosely secure them to the outside of the knee component to anchor the knee component in place. Be sure the patient does not tuck the bands under the knee component. The hook tabs should never make contact with the skin or undersleeve.

**Step 11:**

Once the knee component is secure, instruct the patient to unroll the bottom two bands of the knee component and to pull them above and below each other across the front of the leg. Once firm and comfortable compression is achieved, instruct the patient to secure the hook

tab of the bottom band. Before securing the second band from the bottom, instruct the patient to disengage the third band from the bottom. While holding the second and third bands taut, the patient should pull the second band from the bottom to firm and comfortable compression and secure it to the outside of the knee component. Repeat this process for the remaining bands. Instruct the patient to angle each band individually to achieve a snug and wrinkle-free fit while ensuring that there are no gaps between any of the bands of the knee component.

**Step 12:**

If parts of the knee require additional support or a gap exists, use the provided shelf strap as needed. The shelf strap can be customized into as many as three different straps. Once cut to meet the patients needs, attach the provided tabs to both ends of the strap and secure it to the outside of the knee component.

**Step 13:**

If patient is not going to wear a circaid reduction kit upper leg component, then instruct them to roll the top part of the undersleeve over the top edge of the knee component.

**Reminders for the patient:**

- The knee component should feel firm and comfortable.
- For improved comfort, the patient may slightly loosen the bands of the knee component for nighttime wear.
- When bands are placed in parallel they should slightly overlap, leaving no space between bands.
- The knee component should never hurt. If the patient experiences any pain, the knee component should be immediately removed and the patient should consult a physician.

**Removal of the circaid reduction kit knee:** Instruct the patient to detach all bands starting at the top band. Fold each band back onto itself to keep the hook material lint free. It is important that the spine not be removed or repeatedly folded and creased. Lastly, slide the undersleeve off the leg.

#### **circaid reduction kit shelf strap**

##### **Material composition**

55% Nylon  
38% Polyurethane  
4% Polyethylene  
3% Elastane

This garment does not contain latex.

##### **Donning the circaid reduction kit shelf strap:**

**Step 1:**  
Apply the appropriate circaid compression system as directed.

**Step 2:**  
If using the reduction kit shelf strap, trim the strap to the desired lengths or keep as one long band, and secure hook tabs to the strap ends.

**Step 3:**  
If there is an overhang of oedematous skin and additional support or coverage is desired, use the shelf strap to wrap up underneath the area in need, pulling up and around the limb until firm and comfortable compression is achieved and secure to the circaid compression system of choice.

**Step 4:**  
If a gap between the bands occurs and cannot be corrected by adjusting the bands then cover the gap with the shelf

strap pulling to firm and comfortable compression and securing the two hook tabs to the circaid compression system of choice.

##### **Removing the circaid reduction kit shelf straps:**

**Step 1:** Detach the two hook tabs and remove the shelf strap.

---

**Note:**  
Do not tightly roll the hook tabs. This will cause the tabs to curl and to no longer lay flat while in use.

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