medi GmbH & Co. KG  
Medicusstraße 1  
95448 Bayreuth  
Germany  
T +49 921 912-0  
www.medi.de

medi Australia Pty Ltd  
83 Fennell Street  
North Parramatta NSW 2151  
Australia  
T +61 2 9890 8696  
F +61 2 9890 8439  
sales@mediaustralia.com.au  
www.mediaustralia.com.au

medi Austria GmbH  
Adamgasse 16/7  
6020 Innsbruck  
Austria  
T +43-512 57 95 15  
F +43-512 57 95 15 45  
vertrieb@medi-austria.at  
www.medi-austria.at

medi Belgium NV  
Staatsbaan 77/0099  
3945 Ham  
Belgium  
T: +32-11 24 25 60  
F: +32-11 24 25 64  
info@medibelgium.be  
www.medibelgium.be

medi Canada Inc / médi Canada Inc  
104-1375 Lionel-Boulet,  
Varenes, Québec,  
QC Canada J3X 1P7  
T +1 450-583-3317 / +1 800 361 3153  
F +1 888-583-6827  
service@medicanada.ca  
www.medicanada.ca

MAXIS a.s.,  
medi group company  
Slezská 2127/13  
120 00 Prague 2  
Czech Republic  
T: +420 571 633 510  
F: +420 571 616 271  
info@maxis-medica.com  
www.maxis-medica.com

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medi Norway AS  
Vestheimvegen 35  
4250 Kopervik  
Norway  
T: +47 52 84 45 00  
F: +47 52 84 45 09  
post@medinorway.no  
www.medinorway.no

medi Polska Sp. z o.o.  
ul. Łabędzka 22  
44-121 Gliwice  
Poland  
T: +48-32 230 60 21  
F: +48-32 202 87 56  
info@medi-polska.pl  
www.medi-polska.pl

medi Bayreuth Unipessoal Lda.  
Rua do Centro Cultural, no. 43  
1700-106 Lisbon  
Portugal  
T +351 21 843 71 60  
F +351 21 847 08 33  
medi.portugal@medibayreuth.pt  
www.medi.pt


UK Responsible Person  
medi UK Ltd.  
Plough Lane  
Hereford HR4 OEL  
Great Britain  
T +44-1432 37 35 00  
F +44-1432 37 35 10  
enquiries@mediuk.co.uk  
www.mediuk.co.uk

medi

circaid® reduction kit  
upper leg


1


2



3

4

5

## circaid® reduction kit™ upper leg inelastic compression system

Gebrauchsanweisung. Instructions for use. Mode d'emploi. Instrucciones de uso.  
Instruções de uso. Istruzioni per l'uso. Gebruiksaanwijzing. Bruksanvisning.  
Инструкция по использованию. Instrukcja załączania. دليل الاستخدام.  
Használati útmutató. 使用说明. הוראות שימוש. Bruksanvisning. Kasutusjuhend.

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#### Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If undue pain or an unpleasant sensation develops while you are wearing the product, please consult your doctor or orthotist immediately. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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#### circaid reduction kit shelf strap

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55% Nylon  
38% Polyurethan  
4% Polyethylen  
3% Elasthan

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##### Ablegen des circaid reduction kit shelf strap:

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# circaid® reduction kit upper leg

### Intended purpose

The compression system is designed to provide compression to the leg for patients with venous and lymphatic disorders. The device consists of a series of juxtapositioned, inelastic bands extending from the central portion of the garment. The wrap contours the leg and is designed to cover from just above the ankle bone to just below the knee.

### Indications

- Lymphedema
- other forms of edema
- Lipedema
- Chronic venous insufficiency
- Varicose veins
- DVT/Thrombosis Prevention
- Acute DVT
- Post Thrombotic Syndrome
- Venous Stasis Ulcer
- Post Sclerotherapy
- Lipodermatosclerosis

### Contraindications

- Severe Peripheral Arterial Disease
- Uncontrolled Congestive Heart Failure
- Septic Phlebitis
- Phlegmasia Cerulea Dolens
- Untreated infection
- Any circumstance where increased venous and lymphatic return is undesirable

### Risks / Side effects

- Mild to Moderate Peripheral Arterial Disease
- Impaired Sensation - must be able to sense applied pressure
- Intolerance to material(s)
- Pediatric & assisted use - must be able to sense & communicate applied pressure

### Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

### Individualisation

Step 1: This upper leg component consists of two pieces: a large "body" portion that has numbers printed along the top and bottom edge and a smaller "spine" portion. Separate the two pieces from one another and lay them out flat so that the number markings along the bottom and top edge of the body portion face up. Disengage the bottom most band to view all of the number markings along the bottom edge. Do not flip the spine portion upside down.

Step 2: Measure the patient's circumference 15 cm above the patella E2, in centimeters. Measure the patient's circumference at the groin G.

If the E2 and G measurements exceed 70 cm and 80 cm respectively, a wide legging will be needed, which includes an extension piece. To use the extension, align the 70 cm mark on the base of the component body with the 70 cm mark along the body of the extension piece. Do the same thing with the 80 cm mark along the top edge of the component.

Step 3: Using the E2 circumference measurement first, place the bottom left edge of the spine portion on top of the corresponding number marking along the bottom edge of the body portion. Repeat this step with the top left edge of the spine portion for the G measurement. Once the top and bottom spine tabs have been secured, adjust the center spine tab to create a smooth curve.

Each patient's leg measurements will produce different spine angles. There is no one correct spine shape. The spine tabs do not need to reach the top or bottom edge of the body portion. In those cases, the top and bottom tabs should visually connect to the correct number markings. The garment may no longer lay flat once the spine has been set. This is to be expected when dealing with more extreme spine angles.

Step 4: Cut the excess material off of the body portion of the garment leaving about 5 cm of excess. Ensure that enough of an overlap is left as not to expose any of the spine hook and to allow for potential spine adjusting.

**Note:** When dealing with more extreme spine angles, excess material extending beyond the top of the spine portion is to be expected. Using the top of the spine portion as a guide, cut off the excess material.

Step 5: If accidental removal of the spine tabs is a concern, place the stays along the edge of the spine where the patient feels the spine needs reinforcing or trace around the spine tabs with a permanent marker for easy re-application.

**Tip:** Once the patient reduces significantly, E2 and G circumference measurements can be re-taken and the upper leg component can be adjusted to the patient's new size by resetting the spine per steps 1-7.

Step 6: The undersleeve is worn for added comfort and when exudate is expected. Instruct the patient to slide the undersleeve onto the leg, up to the groin. The undersleeve may be trimmed in length if desired, keeping in mind the patients expected use of circaid reduction kit lower leg and/or circaid reduction kit knee in conjunction with the upper leg component.

Step 7: To apply the upper leg component, the patient should position it around the leg with the black side of the material facing the skin and with the top edge just below the groin. Instruct the patient to unroll the top two bands, pull them above and below each other across the front of the leg, and loosely secure them to the outside of the upper leg component to anchor the garment in place. Be sure the patient does not tuck the bands under the upper leg component. The hook tabs should never make contact with the skin or undersleeve.

**Tip:** If the upper leg component is too long for the patient, it can be trimmed by 3 cm. Trim along the indicated cut line on the black side of the material. Trim both the body and the spine portion.

Step 8: Once the upper leg component is secure, instruct the patient to unroll the bottom two bands of the upper leg component and to pull them above and below each other across the front of the leg. Once firm and comfortable compression is achieved, instruct the patient to secure the hook tab of the bottom band. Before securing the second band from the bottom, instruct the patient to disengage the third band from the bottom. While holding the second and third bands taut, the patient should pull the second band to firm and comfortable compression and secure it to the outside of the upper leg component. Repeat this process for the remaining bands. Instruct the patient to angle each band individually to achieve a snug and wrinkle-free fit while ensuring that there are no gaps between any of the bands of the upper leg component.

Step 9: Bands should be pulled to a firm but comfortable tension to provide therapeutic compression within the tolerance level of the patient. If the patient and therapist desire, the BTS card can be used to identify and record the appropriate tension settings. It is not required to have all bands set to the same tension level. If more than one level of tension is appropriate for treatment, use the BTS card to identify and record the appropriate tension levels. To do so, align one of the vertical guidelines on the band being calibrated with the triangle on the card. Then record where the second vertical guideline falls. This will allow the patient to safely repeat the settings on their own, away from the clinic.

Record tension setting for each set of markings so the patient can easily repeat the process and ensure it is donned to the same or higher tension.

Step 10: If parts of the leg require additional support or a gap exists, use the provided shelf strap as needed. The shelf strap can be customized into as many as three different straps. Attach the provided tabs to both ends of the strap and instruct the patient to secure it to the outside of the legging.

Step 11: Instruct the patient to roll the top part of the undersleeve over the top edge of the upper leg component.

#### Reminders for the patient:

- The upper leg component should feel firm and comfortable.
- The patient may readjust the bands as necessary throughout the day, but should ensure that all vertical guidelines are set to the appropriate level of tension.
- Each set of vertical guidelines may be at a different location on the leg depending on the size and shape of the leg. They do not need to be aligned down the front center of the leg or in any other specific location on the leg.
- For improved comfort, the patient may slightly loosen the bands of the upper leg component for nighttime wear.
- When bands are placed in parallel they should slightly overlap, leaving no space between bands.
- The upper leg component should never hurt. If the patient experiences any pain, the upper leg component should be immediately removed and the patient should consult a physician.

#### Removal of the circaid reduction kit upper leg component:

Instruct the patient to detach all bands starting at the top band. Fold each band back onto itself to keep the hook material lint free. It is important that none of the spine tabs be removed. Lastly, slide the undersleeve off the leg.

#### circaid reduction kit shelf strap

##### Material composition

55% Nylon  
38% Polyurethane  
4% Polyethylene  
3% Elastane

#### Donning the circaid reduction kit shelf strap:

Step 1: Apply the appropriate circaid compression system as directed.

Step 2: If using the reduction kit shelf strap, trim the strap to the desired lengths or keep as one long band, and secure hook tabs to the strap ends.

Step 3: If there is an overhang of oedematous skin and additional support or coverage is desired, use the shelf strap to wrap up underneath the area in need, pulling up and around the limb until firm and comfortable compression is achieved and secure to the circaid compression system of choice.

Step 4: If a gap between the bands occurs and cannot be corrected by adjusting the bands then cover the gap with the shelf strap pulling to firm and comfortable compression and securing the two hook tabs to the circaid compression system of choice.

#### Removing the circaid reduction kit shelf straps:

Step 1: Detach the two hook tabs and remove the shelf strap.

**Note:** Do not tightly roll the hook tabs. This will cause the tabs to curl and to no longer lay flat while in use.

#### Care instructions

We recommend using a laundry bag.

 Machine wash warm - gentle cycle

-  Do not bleach
-  Tumble dry, low heat
-  Do not iron
-  Do not dry clean
-  Do not wring

### Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

### Lifespan

2 months - Due to material wear-and-tear and the medical efficacy can only be guaranteed for a defined lifespan. This assumes correct handling of the garment (e.g. proper care, donning, and doffing).

### Material composition

54% Polyurethane  
42% Nylon  
4% Elastane

### Disposal

The product can be disposed of in the domestic waste.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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