medi GmbH & Co. KG  
Medicusstraße 1  
95448 Bayreuth  
Germany  
T +49 921 912-0  
F +49 921 912-780  
ortho@medi.de  
www.medi.de

medi Australia Pty Ltd  
83 Fennell Street  
North Parramatta NSW 2151  
Australia  
T +61-2 9890 8696  
F +61-2 9890 8439  
sales@medi australia.com.au  
www.medi australia.com.au

medi Austria GmbH  
Adamgasse 16/7  
6020 Innsbruck  
Austria  
T +43 512 57 95 15  
F +43 512 57 95 15 45  
vertrieb@medi-austria.at  
www.medi-austria.at

medi Belgium NV  
Staatsbaan 77/0099  
3945 Ham  
Belgium  
T: +32-11 24 25 60  
F: +32-11 24 25 64  
info@medibelgium.be  
www.medibelgium.be

medi Canada Inc / médi Canada Inc  
104-1375 Lionel-Boulet,  
Varennes, Québec,  
QC Canada J3X 1P7  
T +1 450-583-3317 / +1 800 361 3153  
F +1 888-583-6827  
service@medicanada.ca  
www.medicana.ca

MAXIS a.s.,  
medi group company  
Slezská 2127/13  
120 00 Prague 2  
Czech Republic  
T: +420 571 633 510  
F: +420 571 616 271  
info@maxis-medica.com  
www.maxis-medica.com

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**Important notes**  
This medical device **[MD]** is made only for multiple uses on a single patient . If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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Polyamid, Polyester, Elastan, Viskose

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# Manumed active

## Intended purpose

Manumed active is a wrist compression support.

---

## Indications

All indications in which support of the wrist by a strengthening or pad is necessary, such as:

- Wrist injuries, e.g. sprains, distortions
- Osteoarthritis of the wrist
- Tenosynovitis
- Instability

## Contraindications

None known at present.

## Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in conjunction with sweating) or due to the composition of the material.

## Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

## Note on wearing

In principle, the support can be worn for the whole day but should be taken off

during longer periods of rest (e.g. sleep). But you should make this dependent on how you feel.

The support is suitable for use in water. After use, please rinse out thoroughly.

### Putting on / Taking off

- Before you put the bandage on, stick the attached strap to the middle of the cloth-covered support. The wide end of the strap should be attached to the cloth-covered support so that the thumb opening in the bandage and the strap are on the same side (Fig. 1).
- Now pull the bandage over your hand. When put on correctly, the “m” can be seen on the back of your hand (Fig. 2).
- Fasten the strap so that the opening in the strap goes over your outside wristbone (Fig. 3).
- Then attach it by affixing the narrow end of the strap to the strap (Fig. 4).
- Ensure that the support is positioned correctly when you put the bandage on. The support should be lined up with the underside of your forearm as straight as possible.

To take it off, undo the strap. Push the palm of your hand out of the bandage. Apply pressure to the inside pockets and guides, and pull the bandage off your hand.

To remove the splint, fold the knitted fabric back at the pink edging (Fig. 6). The splint can be pushed out of its pocket and through the slit in the fabric.

### Adjusting the splint

The splint can be individually adjusted at a temperature of around 200 °C.

---

### Care instructions

Hook and loop fasteners should be closed for washing. The splint must be removed before washing. We recommend using a laundry bag. Fabric softeners, fats, oils, lotions, ointments and soap residues can affect the material, hinder the Clima Comfort effect and cause skin irritations.

- Wash the product by hand, preferably using a medi clean detergent, or in delicate cycle at 30°C using a mild detergent without fabric conditioners.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

**Storage instructions**

Keep the product in a dry place and do not expose to direct sunlight.

---

**Composition**

Polyamide, polyester, elastane, viscose

**Liability**

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

**Disposal**

The product can be disposed of in the domestic waste.

**Your medi Team**

Wishes you a rapid recovery!

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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