medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth
Germany
T +49 921 912-0
F +49 921 912-780
ortho@medi.de
www.medi.de

medi Australia Pty Ltd
83 Fennell Street
North Parramatta NSW 2151
Australia
T +61-2 9890 8696
F +61-2 9890 8439
sales@mediaustralia.com.au
www.medi.australia.com.au

medi Austria GmbH
Adamgasse 16/7
6020 Innsbruck
Austria
T +43 512 57 95 15
F +43 512 57 95 15 45
vertrieb@medi-austria.at
www.medi-austria.at

medi Belgium NV
Staatsbaan 77/0099
3945 Ham
Belgium
T: +32-11 24 25 60
F: +32-11 24 25 64
info@medibelgium.be
www.medibelgium.be

medi Canada Inc / médi Canada Inc
104-1375 Lionel-Boulet,
Varenes, Québec,
QC Canada J3X 1P7
T +1 450-583-3317 / +1 800 361 3153
F +1 888-583-6827
service@medicanada.ca
www.medicana.ca

MAXIS a.s.,
medi group company
Slezská 2127/13
120 00 Prague 2
Czech Republic
T: +420 571 633 510
F: +420 571 616 271
info@maxis-medica.com
www.maxis-medica.com

medi Danmark ApS
Vejløgardsvej 59
2665 Vallensbaek Strand
Denmark
T +45-70 25 56 10
kundeservice@sw.dk
www.meditandmark.dk

medi Bayreuth Espana SL
C/Canigo 2-6 bajos
Hospitalet de Llobregat
08901 Barcelona
Spain
T +34-932 60 04 00
F +34-932 60 23 14
medi@mediespana.com
www.mediespana.com

medi France
Z.I. Charles de Gaulle
25, rue Henri Farman
93297 Tremblay en France Cedex
France
T +33-1 48 61 76 10
F +33-1 49 63 33 05
infos@medi-france.com
www.medi-france.com

medi Hungary Kft.
Bokor u. 21.
1037 Budapest
Hungary
T +36 1 371-0090
F +36 1 371-0091
info@medi.hu
www.medi.hu

medi Japan K.K.
5-1 Nihonbashi Kabutocho Chuoku
103-0026 Tokyo
Japan
T: +81 3 6778 2590
F: +81 3 5847 7901
info@medi-japan.jp
www.medi-japan.co.jp

UK Responsible Person
medi UK Ltd.
Plough Lane
Hereford HR4 0EL
Great Britain
T +44-1432 37 35 00
F +44-1432 37 35 10
enquiries@mediuk.co.uk
www.mediuk.co.uk

medi

medi Epibrace

Gebrauchsanweisung. Instructions for use. Mode d'emploi. Instrucciones de uso. Instruções para aplicação. Istruzioni per l'uso. Gebruiksaanwijzing. Bruksanvisning. Návod k použití. Uputa za uporabu. Инструкция по использованию. Kullanna kilavuzu. Instrukcja zakładowa. Οδηγία εφαρμογής. Használati útmutató. Інструкція з використання. Navodila za uporabo. دليل الاستخدام. 使用说明. Instrucțiuni de utilizare. שימוש הוראות שימוש. Bruksanvisning. Kasutusjuhend. Uzvilksšanas instrukcija. Užsimovimo instrukcija. 使用説明書.

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Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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medi Trading (Shanghai) Co. Ltd.
Room 1703, T1, CES WEST BUND CENTER,
No. 277, Long Lan Road
Xuhui District
200232 Shanghai
People's Republic of China
T: +86-21 50582319
F: +86-21 50582319

Medi Turk Ortopedi Medikal Ic ve Dis Tic. AS
Mustafa Kemal Mah. 2156 Sokak No: 16 /4-5
06510 Çankaya Ankara
Turkey
T: +90 312 435 20 26
F: +90 312 434 17 67
info@medi-turk.com
www.medi-turk.com

medi Ukraine LLC
Evhen Sverstiuka str, 11.
Kiev 02002
Ukraine
T: +380 44 591 11 63
F: +380 44 392 73 73
info@medi.ua
http://medi.ua

medi USA L.P.
6481 Franz Warner Parkway
Whitsett, N.C. 27377-3000
USA
T +1-336 4 49 44 40
F +1-888 5 70 45 54
info@mediusa.com
www.mediusa.com

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Polyester, Polyamid, Polyurethan, Polypropylen

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medi Epibrace

Intended purpose

medi Epibrace is an epicondylitis brace for stress relief  on the muscle origins.

Indications

All indications where it is necessary to relieve muscle origins at the ulnar and / or radial epicondyle, such as:
Epicondylitis (e.g. tennis or golfer's elbow)

Contraindications

Impingement syndrome (pressure irritation of the nerves supplying the arm)

Risks / Side Effects

If the aids are too tight, this can cause pressure sores or constriction of blood vessels or nerves. Consult your doctor before using this device if you have the following conditions:

- A condition affecting the application area
- Paraesthesia, circulatory disorders
- Lymph drainage disorders

Tightly fitting aids can cause skin irritation

(pressure, sweat, material).

Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Note on wearing

During longer periods of rest (e.g. when sleeping) you should take the medi Epibrace off.

Applying the strap

- Open the strap on the brace and lie it down, fully unrolling it, in front of you. Position your forearm on the brace with your palm facing upwards so that the light grey positioning mark on the pad is positioned on the side that

you feel pain. Then attach the pad to the positioning mark so that the tip of the pad is parallel to the elbow joint (Fig. 1). Example: If you have tennis elbow (Epicondylitis humeri radialis) in your right arm, the positioning mark should now be on the right side of your right arm.

- Thread the strap through the plastic loop and attach the fastening. The distance between the edge of the fastening and the outer elbow should be approx. 3 cm when the elbow is bent (Fig. 2).
- Refer to the table for guidance values on how tight the fastening can be done up (Fig. 3). Your forearm circumference (measured 5 cm below the elbow joint) should be used as the reference value. Tighten the strap until you attain the value shown on the chart. You can adjust the pressure according to your wishes.
- The clasp should now be firmly positioned on your arm, without slipping or constricting the arm. Form a fist with your hand. You should be able to feel your muscle under the pad. When your muscles are relaxed, the brace should sit against the pad without a noticeable feeling of pressure.

- In principle, the product should be used in accordance with medical supervision and should be fitted by a trained specialist, i.e. an orthopaedic technician, the first time that you put it on. If you wish to use other products with this product, then this should be discussed and agreed with your GP.
- If your arm or hand becomes numb, loosen the aid or remove it if necessary. If your condition persists, consult a doctor.

Care instructions

Hook and loop fasteners should be closed for washing. Soap residues, lotions and ointments can cause skin irritation and material wear.

- Wash the product by hand, preferably using medi clean washing agent.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

Storage instructions

Keep the product in a dry place and do not

expose to direct sunlight.

Material composition

Polyester, polyamide, polyurethane,
polypropylene

Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

Disposal

The product can be disposed of in the domestic waste.

Your medi team

wishes you a quick return to full fitness.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration

in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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