medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth
Germany
T +49 921 912-0
F +49 921 912-780
ortho@medi.de
www.medi.de

medi Australia Pty Ltd
 83 Fennell Street
 North Parramatta NSW 2151
 Australia
 T +61-2 9890 8696
 F +61-2 9890 8439
 sales@mediaustralia.com.au
 www.mediaustralia.com.au

medi Austria GmbH
 Adamgasse 16/7
 6020 Innsbruck
 Austria
 T +43 512 57 95 15
 F +43 512 57 95 15 45
 vertrieb@medi-austria.at
 www.medi-austria.at

medi Belgium NV
 Staatsbaan 77/0099
 3945 Ham
 Belgium
 T: + 32-11 24 25 60
 F: +32-11 24 25 64
 info@medibelgium.be
 www.medibelgium.be

medi Canada Inc / médi Canada Inc
 104-1375 Lionel-Boulet,
 Varennes, Québec,
 QC Canada J3X 1P7
 T +1 450-583-3317 / +1 800 361 3153
 F +1 888-583-6827
 service@medicanada.ca
 www.medicanada.ca

MAXIS a.s.,
 medi group company
 Slezská 2127/13
 120 00 Prague 2
 Czech Republic
 T: +420 571 633 510
 F: +420 571 616 271
 info@maxis-medica.com
 www.maxis-medica.com

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medi Trading (Shanghai) Co. Ltd.
 Room 1703, T1, CES WEST BUND CENTER,
 No. 277, Long Lan Road
 Xuhui District
 200232 Shanghai
 People's Republic of China
 T: +86-21 50582319
 F: +86-21 50582319

UK Responsible Person
 medi UK Ltd.
 Plough Lane
 Hereford HR4 0EL
 Great Britain
 T +44-1432 37 35 00
 F +44-1432 37 35 10
 enquiries@mediuk.co.uk
 www.mediuk.co.uk

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Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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Aluminium, Polyamid, Elasthan, Polyurethan

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medi Epico active

Intended purpose

medi Epico active is an elbow brace for guidance with extension- and/or flexion regulation.

Indications

All indications that require reinforcement of the physiological stabilisation of the elbow joint and/or stress relief of the joint, such as:

- Conservative implementation for early functional treatment following elbow dislocation which only impacts the ligaments
- For medial and/or lateral elbow instability
- Post-operative implementation for osteochondrosis dissecans
- Post-operative implementation following operative cartilage replacement therapy
- Post-operative implementation following operative ligament suturing or ligament reconstruction

Contra-indications

- Concomitant neurological damage, e.g. ulnar nerve
- Compartment syndrome

Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paraesthesia and circulatory disorders

- (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in conjunction with sweating) or due to the composition of the material.

Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Setting the restrictions for flexion and extension and adjustment (can only be carried out by an orthopaedic technician)

Extension restriction wedges, E (stretch restriction): 10°, 20°, 30°, 45°

Flexion restriction wedges, F (bend restriction): 45°, 60°, 75°, 90°

- Your support comes with 20° extension wedges pre-installed.
- To change the wedges, pull the joint splints at the top end (strap 2) out of the joint pocket on the side (fig. 1).
- Insert the flexion wedges (F) and then the desired extension wedges (E) into the side of the articulated joints and fasten them with the attached screws (fig. 2).
- The joint splints on the side are made

of aluminium and can be positioned crosswise if required.

- Place the joint splints, narrow side first, back in the joint pockets and ensure that the straps for the pockets enclose the upper end of the articulated joint (fig. 3).

Putting on the support

- Open and loosen all straps. Pull the support on. The arm that you are using it on should be bent slightly. Your palm should be facing upwards (fig. 4).
- The guide line should be level with the bend in your elbow (fig. 5). Check that the support is positioned correctly by bending and straightening your elbow.
- Fasten all of the ends of the straps in the correct order (see numbering) and check that the support is positioned correctly again (fig. 6).
- The strap may be cut to the length required by the individual. To this end, open the hook and loop fasteners at the ends of the straps and shorten the straps accordingly (fig. 7).
- If necessary, the enclosed cushion can be attached to the external side, between the joint pocket and arm (fig. 8).

Care instructions

Please close the hook and loop fasteners before washing and remove the splint (fig. 1). Soap residues, lotions and ointments can cause skin irritation and material wear.

- Wash the product by hand, preferably using medi clean washing agent.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

Material composition

Aluminum, polyamide, elastane, PU

Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

Disposal

The product can be disposed of in the domestic waste.

Your medi team

Wishes you a speedy recovery!

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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