**medi GmbH & Co. KG**  
**Medicusstraße 1**  
**D-95448 Bayreuth**  
**Germany**  
**T +49 921 912-0**  
**www.medi.de**

medi Australia Pty Ltd  
 83 Fennell Street  
 North Parramatta NSW 2151  
 Australia  
 T +61-2 9890 8696  
 F +61-2 9890 8439  
 sales@mediaustralia.com.au  
 www.mediaustralia.com.au

medi Austria GmbH  
 Adamgasse 16/7  
 6020 Innsbruck  
 Austria  
 T +43 512 57 95 15  
 F +43 512 57 95 15 45  
 vertrieb@medi-austria.at  
 www.medi-austria.at

medi Belgium bvba  
 Posthoornstraat 13/1  
 3582 Koersel  
 Belgium  
 T: + 32-11 24 25 60  
 F: +32-11 24 25 64  
 info@medibelgium.be  
 www.medibelgium.be

medi Brasil  
 Rua Cristovan de Vita, n.º 260, Un. 30,  
 Centro Logístico Raposo Tavares,  
 Vila Camargo  
 Vargem Grande Paulista - SP  
 CEP 06730-000  
 Brasil  
 T: +55-11-3500 8005  
 sac@medibrasil.com  
 www.medibrasil.com

medi Canada Inc / médi Canada Inc  
 597, Rue Duvernay, Verchères  
 QC Canada J0L 2R0  
 T: +1 450-583-3317 / +1 800 361 3153  
 F: +1 888-583-6827  
 service@medicanada.ca  
 www.medicanada.ca

MAXIS a.s.,  
 medi group company  
 Slezská 2127/13  
 120 00 Prague 2  
 Czech Republic  
 T: +420 571 633 510  
 F: +420 571 616 271  
 info@maxis-medica.com  
 www.maxis-medica.com

medi Danmark ApS  
 Vejlegardsvej 59  
 2665 Vallensbaek Strand  
 Denmark  
 T +45-70 25 56 10  
 F +45-70 25 56 20  
 kundeservice@sw.dk  
 www.medit danmark.dk

medi Bayreuth Espana SL  
 C/Canigo 2-6 bajos  
 Hospitalet de Llobregat  
 08901 Barcelona  
 Spain  
 T +34-932 60 04 00  
 F +34-932 60 23 14  
 medi@mediespana.com  
 www.mediespana.com

medi France  
 Z.I. Charles de Gaulle  
 25, rue Henri Farman  
 93297 Tremblay en France Cedex  
 France  
 T +33-1 48 61 76 10  
 F +33-1 49 63 33 05  
 infos@medi-france.com  
 www.medi-france.com

medi Hungary Kft.  
 Bokor u. 21.  
 1037 Budapest  
 Hungary  
 T +36 1 371-0090  
 F +36 1 371-0091  
 info@medi.hu  
 www.medi.hu

medi Nederland BV  
 Heusing 5  
 4817 ZB Breda  
 The Netherlands  
 T +31-76 57 22 555  
 F +31-76 57 22 565  
 info@medi.nl  
 www.medi.nl






**circaid® customizable  
 interlocking ankle foot wrap  
 inelastic compression system**

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## circaid® customizable interlocking afw

### Purpose

The compression system is designed to provide compression to the foot for patients with venous and lymphatic disorders. The device consists of a series of juxtapositioned, inelastic bands extending from the central portion of the garment. The wrap contours the foot and is designed to cover from the base of the toes to just above the ankle, overlapping with the legging.

### Indications

- Lymphedema
- other forms of edema
- Lipedema
- Chronic venous insufficiency
- Varicose veins
- DVT/Thrombosis Prevention
- Post Thrombotic Syndrome
- Venous Stasis Ulcer
- Post Sclerotherapy
- Lipodermatosclerosis

### Contraindications

- Severe Peripheral Arterial Disease
- Uncontrolled Congestive Heart Failure
- Septic Phlebitis
- Phlegmasia Cerulea Dolens
- Untreated infection
- Any circumstance where increased venous and lymphatic return is undesirable
- Suspected or known untreated acute DVT

### Cautions

- Mild to Moderate Peripheral Arterial Disease
- Impaired Sensation - must be able to sense applied pressure
- Intolerance to material(s)
- Pediatric & assisted use - must be able to sense & communicate applied pressure

### Storage instructions

-  Store in a dry place
-  Keep out of sunlight

### Lifespan

Due to material wear-and-tear and the medical efficacy can only be guaranteed for a defined lifespan. This assumes correct handling of the garment (e.g. proper care, donning, and doffing).

circaid customizable  
interlocking afw  
<3 months - depends on the wear and tear

### Disposal

Please dispose of the garment in household waste. No special criteria for disposal apply for medical compression garments.

### Material composition

circaid compression garment  
Polyurethane 49%  
Nylon 44%  
Elastane 4%

For all other components, please refer to the inserts included with each item.

medi Polska Sp. z o.o.  
ul. Łabędzka 22  
44-121 Gliwice  
Poland  
T: +48-32 230 60 21  
F: +48-32 202 87 56  
info@medi-polska.pl  
www.medi-polska.pl

medi Bayreuth Unipessoal, Lda  
Rua do Centro Cultural, no. 43  
1700-106 Lisbon  
Portugal  
T +351-21 843 71 60  
F +351-21 847 08 33  
medi.portugal@medibayreuth.pt  
www.medi.pt

medi RUS LLC  
Business Center NEO GEO  
Butlerova Street 17  
117342 Moscow  
Russia  
T +7-495 374 04 56  
F +7-495 374 04 56  
info@medirus.ru  
www.medirus.ru

medi Trading (Shanghai) Co. Ltd.  
Room 105, Building No. 4,  
Lane 500 ShengXia Road  
Pudong New District  
201210 Shanghai  
People's Republic of China  
T: +86-21 61761988  
F: +86-21 61769607  
infocn@medi.cn

Medi Turk Ortopedi Medikal Ic ve Dis Tic.  
AS  
Mustafa Kemal Mah. 2156 Sokak No: 16  
/4-5  
06510 Çankaya Ankara  
Turkey  
T: +90 312 435 20 26  
F: +90 312 434 17 67  
info@medi-turk.com  
www.medi-turk.com

medi UK Ltd.  
Plough Lane  
Hereford HR4 OEL  
Great Britain  
T +44-1432 37 35 00  
F +44-1432 37 35 10  
enquiries@mediuk.co.uk  
www.mediuk.co.uk

medi Ukraine LLC  
Evhena Sverstiuka str, 11.  
Kiev 02002  
Ukraine  
T: +380 44 591 11 63  
F: +380 44 392 73 73  
info@medi.ua  
http://medi.ua

medi USA L.P.  
6481 Franz Warner Parkway  
Whitsett, N.C. 27377-3000  
USA  
T +1-336 4 49 44 40  
F +1-888 5 70 45 54  
info@mediusa.com  
www.mediusa.com

This garment does not contain latex.

### Care instructions

The garment can be machine washed using a gentle cycle and dried on a low heat setting. If doing so, it is recommended that you place the garment in a pillow case to aid in keeping the hook and loop tab lint-free and to protect other clothing from getting damaged. To extend the life of the garment, hand wash and drip dry.

Machine wash warm

Do not bleach

Tumble dry, low heat

Do not iron

Do not dry clean

Do not wring

### Wearing recommendation

Slightly loosen the bands of the garment for night-time wear.

If experiencing any pain, or discomfort, reduce the compression (mmHg) applied or remove device. Inform your clinician of your action.

When bands are placed correctly they should slightly overlap, leaving no space between bands.

Do not tuck bands. Hook tabs should never touch the skin!

### Individualisation

The ankle foot wrap is not designed to be cut or modified except as outlined in the directions for use.

Doing so other than instructed will void the warranty.

Be careful that the reinforced seam within the hook material is not cut while shortening the afw.

The circaid customizable interlocking ankle foot wrap can be tailored to fit the patient's foot size and contour.

The ankle foot wrap should be fitted by a medical professional.

### Donning instructions

Step 1: Apply a circaid sock first. Ensure that there are no wrinkles in the sock. (1)

Step 2: Identify the two interlocking heel bands. (2a) Fasten the interlocking heel bands by sliding the shorter band through the hole in the other and loosely securing them to the outside of the AFW. (2b)

Step 3: Slide the foot in the AFW so that the heel is securely placed in the heel pocket. (3) Ensure that the foot is properly positioned. The AFW may need to be pulled forward to ensure a proper fit.

Step 4: Unroll the foot band closest to the toes to the opposite side of the foot wrap and trim to desired length. Then secure it with firm and comfortable compression using one of the included hook tabs. (4a) Repeat with the second foot band. (4b) Adjust the interlocking heel bands so that they fit snugly around the back of the heel.

Step 5: Take the long ankle strap and pull it across the front of the foot. Wrap it around the ankle and trim to the desired length before securing it using one of the included hook tabs. (5a-b) Ensure that there are no gaps.

---

Note:

The straps may wrap from different sides or fall on different positions on the foot, depending on which foot the AFW is being worn and the size of the foot.

---

The circaid customizable ankle foot wrap has now been applied. (6)

If using a circaid legging in addition to the AFW, apply the legging first and then apply the AFW over it.

If pressure increases or decreases noticeably during wear, readjust the garment as needed.

### **Doffing instructions**

First detach the long ankle band and fold it back to the body of the AFW to keep it lint free. Lastly, loosen the foot bands. The ankle foot wrap slips off like a shoe.

### **Adjusting the foot length:**

The following steps will reduce the foot-length of the circaid ankle foot wrap (7).

Step 1:

With the beige side facing up, mark the desired length.

Step 2:

Draw a straight line connecting the notches.

Step 3:

Carefully cut along the drawn line, assuring that the reinforced seam section of the Velcro material is not cut.

---

Note:

Do not trim the component while it is on the wearer's leg.

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