medi GmbH & Co. KG
 Medicusstraße 1
 95448 Bayreuth
 Germany
 T +49 921 912-0
 www.medi.de

medi Australia Pty Ltd
 83 Fennell Street
 North Parramatta NSW 2151
 Australia
 T +61-2 9890 8696
 F +61-2 9890 8439
 sales@mediaustralia.com.au
 www.mediaustralia.com.au

medi Austria GmbH
 Adamgasse 16/7
 6020 Innsbruck
 Austria
 T +43 512 57 95 15
 F +43 512 57 95 15 45
 vertrieb@medi-austria.at
 www.medi-austria.at

medi Belgium NV
 Staatsbaan 77/0099
 3945 Ham
 Belgium
 T: +32-11 24 25 60
 F: +32-11 24 25 64
 info@medibelgium.be
 www.medibelgium.be

medi Canada Inc / médi Canada Inc
 104-1375 Lionel-Boulet,
 Varennes, Québec,
 QC Canada J3X 1P7
 T +1 450-583-3317 / +1 800 361 3153
 F +1 888-583-6827
 service@medicanada.ca
 www.medicanada.ca

MAXIS a.s.,
 medi group company
 Slezská 2127/13
 120 00 Prague 2
 Czech Republic
 T: +420 571 633 510
 F: +420 571 616 271
 info@maxis-medica.com
 www.maxis-medica.com

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
medi Trading (Shanghai) Co. Ltd.
 Room 1703, T1, CES WEST BUND CENTER,
 No. 277, Long Lan Road
 Xuhui District
 200232 Shanghai
 People's Republic of China
 T: +86-21 50582319
 F: +86-21 50582319

UK Responsible Person
 medi UK Ltd.
 Plough Lane
 Hereford HR4 0EL
 Great Britain
 T +44-1432 37 35 00
 F +44-1432 37 35 10
 enquiries@mediuk.co.uk
 www.mediuk.co.uk

**circaid® juxtafit® essentials lower leg
 inelastic compression system**

Gebrauchsanweisung. Instructions for use. Mode d'emploi. Instrucciones de uso. Instruções para aplicação. Istruzioni per l'uso. Gebruiksaanwijzing. Brugsvejledning. Bruksanvisning. Návod k použití. Инструкция по использованию. Instrukcja zakładania. Használati útmutató. Інструкції з використання. دليل الاستخدام. Brugsanvisning. Kasutusjuhend. הוראות שימוש. Käyttöohje. 使用説明書.

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


Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If undue pain or an unpleasant sensation develops while you are wearing the product, please consult your doctor or orthotist immediately. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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circaid® juxtafit® essentials lower leg

Intended purpose

The compression system is designed to provide compression to the leg for patients with venous and lymphatic disorders. The device consists of a series of juxtapositioned, inelastic bands extending from the central portion of the garment. The wrap contours the leg and is designed to cover from just above the ankle bone to just below the knee.

Indications

- Lymphedema
- Lipedema
- Chronic venous insufficiency
 - Varicose Veins
 - Lipodermatosclerosis
 - Venous Stasis Ulcer
- Deep vein thrombosis/Thrombosis Prevention
- Acute deep vein thrombosis
- Post Thrombotic Syndrome
- Post Sclerotherapy
- Dependent oedema

Contraindications

- Severe Peripheral Arterial Disease
- Uncontrolled Congestive Heart Failure
- Septic Phlebitis
- Phlegmasia Cerulea Dolens
- Untreated infection
- Any circumstance where increased venous and lymphatic return is undesirable

Risks/side effects

- Mild to Moderate Peripheral Arterial Disease
- Impaired Sensation - must be able to sense applied pressure
- Intolerance to material(s)
- Pediatric & assisted use - must be able to sense & communicate applied pressure

Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this

from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Wearing recommendation

Note: Unless your treating physician orders or recommends otherwise, we recommend wearing the the legwear together with the respective footwear included in the set. If in doubt, please consult your treating physician.

Always ensure that the correct prescribed pressure range is being applied. Slightly loosen the bands of the garment for night-time wear. Each vertical BPS line on the bands may be at a different location on the leg depending on the size and shape of the leg. They do not need to be aligned down the front center of the leg or in any other specific location on the leg. If experiencing any pain, or discomfort, reduce the compression (mmHg) applied or remove device. Inform your clinician of your action. When bands are placed correctly they should slightly overlap, leaving no space between bands. Do not tuck bands. The Velcro tabs should never touch the skin!

Donning instructions

- Slide the circaid undersleeve onto the lower leg. The cuffed edge of the undersleeve should cover the heel. Ensure that there are no wrinkles in the undersleeve.
- Align the garment just above the ankle bone and below the knee crease. The black side of the material should face the skin. Unroll the top two bands and loosely secure them to the outside of the garment to anchor the legging in place. This enables you to easily raise or lower the garment for correct positioning.

Note: *The bottom most band can be identified by the care label.*

- Unroll the bottom two bands and secure the bottom most band to the outside of the legging to a firm and comfortable compression level.
- While holding the second band, detach the next band. Secure the second band and continue this process with the other bands, going up the leg.
- The garment should lay flat and wrinkle free against the leg. Adjust each band to a firm and comfortable compression level. Ensure that the bands are placed in an alternating order.
- Once all bands are secured, inspect the garment for any gaps or creases. Adjust the bands as necessary.

Note: *Once the correct pressure is set using the BPS card, you may fold the top and bottom part of the circaid undersleeve over the top and bottom edge of the circaid juxtaft garment.*

Donning the circaid compression ankle

First fold the bottom part of the circaid undersleeve over the bottom edge of the circaid garment. Slide the circaid compression ankle onto the foot and ensure that it is pulled high enough so that it overlaps with the circaid juxtaft garment. Ensure that there are no wrinkles in the ankle.

Doffing instructions

To remove the legging, first remove the cover up, if using. Then remove the knee cover, shelf straps, and/or foot

compression, if using any of these. Then unfold the circaid undersleeve. Detach all bands of the legging starting with the top band. Fold each band onto itself to keep the hook and loop tabs lint-free. Remove the legging and then remove the circaid undersleeve or sock.

Using the Built-In-Pressure card

The set compression value can vary due to movement, specific material characteristics and individual body anatomy.

- Locate the Built-In-Pressure card in your packaging.
- Identify your prescribed pressure scale using the color-coded system.
- Starting with the bottom band, line up the black triangle on the card with one of the BPS lines on the bottom band.
- Note where the second BPS line lines up with the card's compression ranges (A).
- If the BPS line on the garment either falls short (B) or goes beyond the correct compression range (C), readjust the band as necessary so that the second BPS line on the garment is aligned with the prescribed compression range (A). Completely loosen the band before reapplying it if too much tension has been applied.
- Repeat the steps 3 - 5 for each band going up the leg.
- Be sure to adjust all bands to the same tension to ensure gradient compression.

Bands may need adjusting throughout the day to maintain the prescribed compression range as edema is reduced.

Note:

It is recommended that the limb is always treated in its entirety with medical compression. This means that compression products for the leg begin at the foot or the ankle and extend up to below the hollow of the knee or beyond to include the thigh.

Care instructions

This compressive product can be washed on a delicate cycle and dried in a tumble dryer set to a low heat. We recommend using a mesh laundry bag. To extend the life of the product, hand wash and drip dry.

-
-  Machine wash warm
 -  Do not bleach
 -  Tumble dry, low heat
 -  Do not iron
 -  Do not dry clean
 -  Do not wring
-

Storage instructions

Store the product in a dry place and do not expose it to direct sunlight.

Lifespan

6 months – due to material wear and/or slackening, medical efficacy can only be guaranteed for a defined lifespan. This assumes correct handling of the product (e.g. during care and when putting on and taking off the product).

Material composition

53% Polyurethane, 42% Nylon, 5% Elastane

Disposal

You can dispose of this product with your household waste

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

FRANÇAIS / FRENCH

circaid® juxtafit® essentials lower leg

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