medi GmbH & Co. KG
 Medicusstraße 1
 95448 Bayreuth
 Germany
 T +49 921 912-0
 www.medi.de

medi Australia Pty Ltd
 Unit 4/13-21 Hallmark Street
 Pendle Hill NSW 2145
 Australia
 T +61 2 9890 8696
 F +61 2 9890 8439
 sales@medi.au
 www.mediaustralia.com.au

medi Austria GmbH
 Adamgasse 16/7
 6020 Innsbruck
 Austria
 T +43-512 57 95 15
 F +43-512 57 95 15 45
 vertrieb@medi-austria.at
 www.medi-austria.at

medi Belgium NV
 Staatsbaan 77/0099
 3945 Ham
 Belgium
 T: + 32-11 24 25 60
 F: +32-11 24 25 64
 info@medibelgium.be
 www.medibelgium.be

medi Canada Inc / médi Canada Inc
 104-1375 Lionel-Boulet,
 Varennes, Québec,
 QC Canada J3X 1P7
 T +1 450-583-3317 / +1 800 361 3153
 F +1 888-583-6827
 service@medicanada.ca
 www.medicanada.ca

MAXIS a.s.,
 medi group company
 Slezská 2127/13
 120 00 Prague 2
 Czech Republic
 T: +420 571 633 510
 F: +420 571 616 271
 info@maxis-medica.com
 www.maxis-medica.com

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medi

circaid® power added compression™ (pac band™) inelastic compression system

Instructions for use. Gebrauchsanweisung. Instructions
 d'utilisation. Istruzioni per l'uso. Instrucciones de uso.
 Instruções de uso. Användningsinstruktioner. Brugsanvisning.
 Instrukcja użytkowania. Návod k použití. Használati útmutató.
 دليل الاستخدام. Інструкції з використання. Інструкція по
 использованию. הוראות שימוש. Kasutusjuhend. Käyttöohje.

**power added compression™
 (pac band™)**

4 051526 949508

E013733/01.2025

medi. I feel better.

medi Bayreuth España SL
 C/Canigó 2 – 6 bajos
 Hospitalet de Llobregat
 08901 Barcelona
 Spain
 T +34 93 260 04 00
 F +34 93 260 23 14
 medi@mediespana.com
 www.mediespana.com

medi France
 Z.I. Charles de Gaulle
 25, rue Henri Farman
 93297 Tremblay en France Cedex
 France
 T +33 1 48 61 76 10
 F +33 1 49 63 33 05
 infos@medi-france.com
 www.medi-france.com

medi Hungary Kft.
 Bokor u. 21.
 1037 Budapest
 Hungary
 T +36 1371 0090
 F +36 1371 0091
 info@medi.hu
 www.medi.hu




medi Italia s.r.l.
 Via Giuseppe Ghedini, 2
 40069 Zola Predosa (BO)
 Italy
 T: +39-051 6 13 24 84
 F: +39-051 6 13 29 56
 info@medi-italia.it
 www.medi-italia.it

medi Japan K.K.
 5-1 Nihonbashi Kabutocho Chuoku
 103-0026 Tokyo
 Japan
 T: +81 3 6778 2590
 F: +81 3 5847 7901
 info@medi-japan.jp
 www.medi-japan.co.jp

Tachezy Sanit Srl.
 Via Druso 48/a
 39100 Bolzano
 Italy
 T: +39 0471 288 068
 F: +39 0471 282 773
 info@tachezysanit.com
 www.tachezysanit.com

medi Nederland BV
 Heusing 5
 4817 ZB Breda
 The Netherlands
 T +31 76 57 22 555
 F +31 76 57 22 565
 info@medi.nl
 www.medi.nl

medi Norway AS
 Vestheimvegen 35
 4250 Kopervik
 Norway
 T: +47 52 84 45 00
 F: +47 52 84 45 09
 post@medinorway.no
 www.medinorway.no

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Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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medi Trading (Shanghai) Co. Ltd.
Room 1703, T1, CES WEST BUND
CENTER,
No. 277, Long Lan Road
Xuhui District
200232 Shanghai
People's Republic of China
T: +86-21 50582319
F: +86-21 50582319

medi Sweden AB
Box 6034
192 06 Sollentuna
Sweden
T +46 8 96 97 98
F +46 8 626 68 70
info@medi.se
www.medi.se

Medi Turk Ortopedi Medikal Iç ve
Dis Tic. AS
Mustafa Kemal Mah. 2156 Sokak
No: 16/4-5
06510 Çankaya Ankara
Turkey
T: +90 312 435 20 26
F: +90 312 434 17 67
info@medi-turk.com
www.medi-turk.com

medi Ukraine LLC
Evhena Sverstiuka str, 11.
Kiev 02002
Ukraine
T: +380 44 591 11 63
F: +380 44 392 73 73
info@medi.ua
http://medi.ua

medi USA L.P.
6481 Franz Warner Parkway
Whitsett, N.C. 27377-3000
USA
T +1-336 4 49 44 40
F +1-888 5 70 45 54
info@mediusa.com
www.mediusa.com

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circaid[®] pac band[™]

Intended purpose

The garment is designed to provide compression  to the foot and ankle region for patients with venous and lymphatic disorders.

Indications

- Lymphedema
- Chronic venous insufficiency
 - Varicose veins
 - Lipodermatosclerosis
 - Venous leg ulcer/Mixed leg ulcer
- Deep vein thrombosis/Thrombosis Prevention
- Acute deep vein thrombosis
- Post Thrombotic Syndrome
- Post Sclerotherapy
- Dependent oedema

Contraindications

- Severe Peripheral Arterial Disease
- Uncontrolled Congestive Heart Failure
- Septic Phlebitis
- Phlegmasia Cerulea Dolens
- Untreated infection
- Any circumstance where increased venous and lymphatic return is undesirable

Risks / Side effects

- Mild to Moderate Peripheral Arterial Disease
- Impaired Sensation - must be able to sense applied pressure
- Intolerance to material(s)
- Pediatric & assisted use - must be able to sense & communicate applied pressure

Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Notes on wearing

Note: *Unless your treating physician orders or recommends otherwise, we recommend wearing the the legwear together with the*

respective footwear included in the set. If in doubt, please consult your treating physician.

The garment should never hurt. It should feel firm but comfortable. Slightly loosen the bands of the garment for night-time wear. If experiencing any pain, remove the garment and consult your physician. When bands are placed correctly they should slightly overlap, leaving no space between bands. Do not tuck bands. The Velcro tabs should never touch the skin!

Fitting instructions

The circaid power added compression band (pac band) is designed to provide additional compression to the foot, targeting the arch and ankle area regions.

For complete compression of the foot, it is recommended to apply the circaid power added compression band over the circaid compression anklet.

- Lay the pac band, out flat. If using on the left foot, the care instructions and illustration should be visible (2a). If using on the right foot, no text should be visible (2b).
- Center the arch of the foot as shown in the illustration printed on the band. (3a) shows the left foot. (3b) shows the right foot. Wrap the short arch band over the top of the foot (4).
- Secure the long arch band over the top of the foot with the large hook tab (5a-b).
- Pull the ankle band around the back of the ankle and secure it with the small hook tab to a firm but comfortable compression (6a-c).

Note: If excess material or overlap is a concern, either the long arch band (7) or the ankle band (8) may be trimmed as needed. Be careful not to over-trim the bands.

Care instructions

The garment can be machine washed using a gentle cycle and dried on a low heat setting. We recommend using a laundry bag.. To extend the life of the garment, hand wash and drip dry.

-
-  Machine wash warm
 -  Do not bleach
 -  Drip dry or line dry
 -  Do not iron
 -  Do not dry clean
 -  Do not wring
-

Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

Lifespan

3 months – due to material wear and/or slackening, medical efficacy can only be guaranteed for a defined lifespan. This assumes correct handling of the product (e.g. during care and when putting on and taking off the product).

Material composition

Nylon, Polyuretane, Elastane

Disposal

Dispose of properly after use.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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