medi GmbH & Co. KG  
 Medicusstraße 1  
 95448 Bayreuth  
 Germany  
 T +49 921 912-0  
 www.medi.de

medi Australia Pty Ltd  
 83 Fennell Street  
 North Parramatta NSW 2151  
 Australia  
 T +61 2 9890 8696  
 F +61 2 9890 8439  
 sales@mediaustralia.com.au  
 www.mediaustralia.com.au

medi Austria GmbH  
 Adamgasse 16/7  
 6020 Innsbruck  
 Austria  
 T +43-512 57 95 15  
 F +43-512 57 95 15 45  
 vertrieb@medi-austria.at  
 www.medi-austria.at

medi Belgium NV  
 Staatsbaan 77/0099  
 3945 Ham  
 Belgium  
 T: +32-11 24 25 60  
 F: +32-11 24 25 64  
 info@medibelgium.be  
 www.medibelgium.be

4 026398 068534

medi Canada Inc / médi Canada Inc  
 104-1375 Lionel-Boulet,  
 Varennes, Québec,  
 QC Canada J3X 1P7  
 T +1 450-583-3317 / +1 800 361 3153  
 F +1 888-583-6827  
 service@medicanada.ca  
 www.medicanada.ca

MAXIS a.s.,  
 medi group company  
 Slezská 2127/13  
 120 00 Prague 2  
 Czech Republic  
 T: +420 571 633 510  
 F: +420 571 616 271  
 info@maxis-medica.com  
 www.maxis-medica.com

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### Important notes

This medical device **[MD]** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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## **circaid® glove/toe cap**

### **Intended purpose**

This garment is designed to provide compression to toes or fingers for patients with venous or lymphatic disorders.

### **Indications**

- Lymphedema
- Other forms of edema
- Acute deep vein thrombosis

### **Contraindications**

- Severe peripheral arterial disease
- Uncontrolled congestive heart failure
- Septic phlebitis
- Phlegmasia cerula dolens
- Untreated infection
- Any circumstance where increased venous and lymphatic return is undesirable

### **Risks / Side effects**

- Mild to moderate peripheral arterial disease
- Impaired sensation - must be able to sense applied pressure
- Intolerance to material(s)
- Pediatric & assisted use - must be able to sense & communicate applied pressure

### **Intended users and patient target groups**

Healthcare professionals and patients, including persons who play a supporting

role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

### **Wearing recommendation**

The garment can be worn during the day or as directed by your clinician. If any discomfort is experienced, remove the the garment.

### **Fitting instructions**

Step 1: Slide the circaid glove/toe cap onto the fingers/toes, ensure that each finger/toe is placed separately within the garment (1a 1b).

Step 2: Fingers (2a) or toes (2b) of the garment can be trimmend down individually until the dotted line. Furter options of individualisation are not intended.

### **Individualisation**

To assess the correct area to trim, place the garment on the hand/toe and mark the desired length. Remove the garment and trim with scissors.

**Note:** *The garment is designed for left or right limbs*

### Care instructions

The garment can be machine washed using a gentle cycle and dried on a low heat setting. To extend the life of the garment, hand wash and drip dry.

-  Machine wash warm
-  Do not bleach
-  Tumble dry, low heat
-  Do not iron
-  Do not dry clean
-  Do not wring

### Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

### Lifespan

3 months - due to material wear-and-tear the medical efficacy can only be expected for a defined lifespan. This assumes correct handling of the garment (e.g. proper care, application, and removal).

### Material composition

45% Polyamide  
55% Elastane

### Disposal

The product can be disposed of in the domestic waste.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report

to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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medi Ukraine LLC