medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth
Germany
T +49 921 912-0
www.medi.de

medi Australia Pty Ltd
83 Fennell Street
North Parramatta NSW 2151
Australia
T +61 2 9890 8696
F +61 2 9890 8439
sales@mediaustralia.com.au
www.mediaustralia.com.au

medi Austria GmbH
Adamgasse 16/7
6020 Innsbruck
Austria
T +43-512 57 95 15
F +43-512 57 95 15 45
vertrieb@medi-austria.at
www.medi-austria.at

medi Belgium NV
Staatsbaan 77/0099
3945 Ham
Belgium
T: +32-11 24 25 60
F: +32-11 24 25 64
info@medibelgium.be
www.medibelgium.be

4 026398 068534

medi Canada Inc / médi Canada Inc
104-1375 Lionel-Boulet,
Varenes, Québec,
QC Canada J3X 1P7
T +1 450-583-3317 / +1 800 361 3153
F +1 888-583-6827
service@medicanada.ca
www.medicanada.ca

MAXIS a.s.,
medi group company
Slezská 2127/13
120 00 Prague 2
Czech Republic
T: +420 571 633 510
F: +420 571 616 271
info@maxis-medica.com
www.maxis-medica.com

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Important notes

This medical device **[MD]** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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circaid® glove/toe cap

Intended purpose

This garment is designed to provide compression to toes or fingers for patients with venous or lymphatic disorders.

Indications

- Lymphedema
- Other forms of edema
- Acute deep vein thrombosis

Contraindications

- Severe peripheral arterial disease
- Uncontrolled congestive heart failure
- Septic phlebitis
- Phlegmasia cerula dolens
- Untreated infection
- Any circumstance where increased venous and lymphatic return is undesirable

Risks / Side effects

- Mild to moderate peripheral arterial disease
- Impaired sensation - must be able to sense applied pressure
- Intolerance to material(s)
- Pediatric & assisted use - must be able to sense & communicate applied pressure

Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting

role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Wearing recommendation

The garment can be worn during the day or as directed by your clinician. If any discomfort is experienced, remove the the garment.

Fitting instructions

Step 1: Slide the circaid glove/toe cap onto the fingers/toes, ensure that each finger/toe is placed separately within the garment (1a 1b).

Step 2: Fingers (2a) or toes (2b) of the garment can be trimmend down individually until the dotted line. Furter options of individualisation are not intended.

Individualisation

To assess the correct area to trim, place the garment on the hand/toe and mark the desired length. Remove the garment and trim with scissors.

Note: *The garment is designed for left or right limbs*

Care instructions

The garment can be machine washed using a gentle cycle and dried on a low heat setting. To extend the life of the garment, hand wash and drip dry.

-  Machine wash warm
-  Do not bleach
-  Tumble dry, low heat
-  Do not iron
-  Do not dry clean
-  Do not wring

Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

Lifespan

3 months - due to material wear-and-tear the medical efficacy can only be expected for a defined lifespan. This assumes correct handling of the garment (e.g. proper care, application, and removal).

Material composition

45% Polyamide
55% Elastane

Disposal

The product can be disposed of in the domestic waste.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report

to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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medi Ukraine LLC