medi GmbH & Co. KG
 Medicusstraße 1
 D-95448 Bayreuth
 Germany
 T +49 921 912-0
 F +49 921 912-780
 ortho@medi.de
 www.medi.de

medi Australia Pty Ltd
 83 Fennell Street
 North Parramatta NSW 2151
 Australia
 T +61-2 9890 8696
 F +61-2 9890 8439
 sales@medi australia.com.au
 www.medi australia.com.au

medi Austria GmbH
 Adamgasse 16/7
 6020 Innsbruck
 Austria
 T +43 512 57 95 15
 F +43 512 57 95 15 45
 vertrieb@medi-austria.at
 www.medi-austria.at

medi Belgium NV
 Staatsbaan 77/0099
 3945 Ham
 Belgium
 T: + 32-11 24 25 60
 F: +32-11 24 25 64
 info@medibelgium.be
 www.medi belgium.be

medi Canada Inc / médi Canada Inc
 597, Rue Duvernay, Verchères
 QC Canada J0L 2R0
 T: +1 450-583-3317 / +1 800 361 3153
 F: +1 888-583-6827
 service@medicanada.ca
 www.medicana da.ca

MAXIS a.s.,
 medi group company
 Slezská 2127/13
 120 00 Prague 2
 Czech Republic
 T: +420 571 633 510
 F: +420 571 616 271
 info@maxis-medica.com
 www.maxis-medica.com

medi Danmark ApS
 Vejlegardsvej 59
 2665 Vallensbaek Strand
 Denmark
 T +45-70 25 56 10
 F +45-70 25 56 20
 kundeservice@sw.dk
 www.medi danmark.dk

medi Bayreuth Espana SL
 C/Canigo 2-6 bajos
 Hospitalet de Llobregat
 08901 Barcelona
 Spain
 T +34-932 60 04 00
 F +34-932 60 23 14
 medi@mediespana.com
 www.medi espana.com

medi France
 Z.I. Charles de Gaulle
 25, rue Henri Farman
 93297 Tremblay en France Cedex
 France
 T +33-1 48 61 76 10
 F +33-1 49 63 33 05
 infos@medi-france.com
 www.medi-france.com

medi Hungary Kft.
 Bokor u. 21.
 1037 Budapest
 Hungary
 T +36 1 371-0090
 F +36 1 371-0091
 info@medi.hu
 www.medi.hu

medi Japan K.K.
 5-1 Nihonbashi Kabutocho Chuoku
 103-0026 Tokyo
 Japan
 T: +81 3 6778 2590
 F: +81 3 5847 7901
 info@medi-japan.jp
 www.medi-japan.co.jp

medi Nederland BV
 Heusing 5
 4817 ZB Breda
 The Netherlands
 T +31-76 57 22 555
 F +31-76 57 22 565
 info@medi.nl
 www.medi.nl

medi Polska Sp. z o.o.
 ul. Łabędzka 22
 44-121 Gliwice
 Poland
 T: +48-32 230 60 21
 F: +48-32 202 87 56
 info@medi-polska.pl
 www.medi-polska.pl

medi Bayreuth Unipessoal, Lda
 Rua do Centro Cultural, no. 43
 1700-106 Lisbon
 Portugal
 T +351-21 843 71 60
 F +351-21 847 08 33
 medi.portugal@medibayreuth.pt
 www.medi.pt

medi RUS LLC
 Business Center NEO GEO
 Butlerova Street 17
 117342 Moscow
 Russia
 T +7-495 374 04 56
 F +7-495 374 04 56
 info@medirus.ru
 www.medi rus.ru

medi Sweden AB
 Box 6034
 192 06 Sollentuna
 Sweden
 T +46 8 96 97 98
 F +46 8 626 68 70
 info@medi.se
 www.medi.se

medi Rehab one

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Important notes

This medical device **MD** is made for single patient use only. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If undue pain or an unpleasant sensation develops while you are wearing the product, please consult your doctor or orthotist immediately. Do not wear the product over open wounds, and use it only as instructed by your doctor or orthotist.

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medi Rehab one

Intended purpose

medi Rehab one is a round-knit medical compression stocking for compressive care of the lower limbs, mainly for prevention and treatment of post-operative and post-traumatic oedema and general swelling.

Indications

Oedema

- Post-traumatic oedema
- Post-operative oedema

Contraindications

- Advanced peripheral arterial occlusive disease (if one of these parameters is present: ABPI < 0.5, ankle arterial pressure < 60 mmHg, toe pressure < 30 mmHg or TcPO₂ < 20 mmHg on dorsum of foot).
- Decompensated heart failure (NYHA III + IV)
- Septic phlebitis
- Phlegmasia cerulea dolens

Please note:

In the following instances, the therapeutic decision should be made by a doctor, who weighs up the benefits and risks and selects the most suitable compressive equipment:

- Severe weeping dermatosis
- Patient unable to tolerate compressive stocking material
- Severe paraesthesia in the limb
- Advanced peripheral neuropathy (e.g. diabetes mellitus)
- Rheumatoid arthritis

Risks and side effects

Medical compression stockings may cause skin necrosis and pressure damage to the peripheral nerves, especially if it is handled incorrectly. If your skin is sensitive, medical compression stockings may cause itching, peeling or signs of inflammation. Therefore, adequate skincare applied under the compression stocking is a must.

Intended patient groups

Healthcare professionals should provide care to the adults and children, applying the available information on the measurements/sizes and necessary functions/indications and in line with the information provided by the manufacturer, and acting under their own responsibility.

Note on wearing

Unless otherwise prescribed, you should wear medi Rehab one from morning to evening. If this duration is significantly exceeded in exceptional cases, the top band should be repositioned several times within this period.

Daily wearing and washing of your compression stockings can cause the compression and elasticity of your stockings, which are required to make them medically effective, to deteriorate over time. The stockings should therefore only be used for a maximum of six months.

Special instructions

- If pain or increased irritation of the skin occurs while you are wearing the stocking, remove it immediately and consult your doctor or the medical specialists.

- Appropriate skincare must be maintained during compression therapy. This is essential. Creams or ointments that are particularly greasy and also soap residue can cause skin irritation and material wear and can thus impair the efficacy of the compression stocking. Therefore, medi also offers a skincare range (medi day, medi night, medi soft foam lotion) that is specially formulated for compression stockings. Please visit a medical specialist retailer to discuss which product is right for you.
- Sharp fingernails or toenails, faults in your shoe's lining, jewellery or rough or calloused heels may cause damage to the knit fabric.
- The sewn-in fabric label is an important component in ensuring that the identification and traceability of the product is guaranteed. Please do not remove it from the stocking by detaching it or cutting it out.

Fitting instructions

- Put the compression stocking on immediately after you get up. If you have a shower in the morning and then apply skin care products, then you should wait a while until your skin is completely dry. This will make it far easier to put on the compression stocking.
- Reach into the stocking. Grasp the heel and pull the stocking inside out. (Figs. 1 + 2)
- Pull the stocking over the foot, up to the heel. (Figs. 3 + 4)
- Now grip the top layer of the stocking as shown in Fig. 5 and pull it onto your lower leg. (Fig. 6)
- Repeat the process step by step until

the stocking is fully on. The green border should come to just under the knee. (Fig. 7-10)

- Check the entire stocking to ensure that it is seated correctly.

Please note the following:

The stocking must lie correctly and should not be twisted (the ribbed pattern provides orientation for this) or stretched too long.

Care instructions

- Wash your stocking every day after you have worn it.
- Ideally, turn the stocking inside out and wash it separately or with items of clothing that are the same colour.
- It is recommended that you use a laundry bag.
- Soap residue, creams or ointments can lead to irritation of the skin and deterioration of the material.
- Please only use medi spot ex to remove stains.
- Wash the product by hand, preferably with medi clean soap, or in a delicates wash cycle at 40°C with mild detergent and without fabric conditioner.
- Do not bleach.
- Dry at a low temperature.
- Do not iron.
- Do not dry-clean.

- After washing, you can roll the compression stocking up in a towel and squeeze it to dry it. Do not leave the damp stocking in the wet towel, instead let it dry in the air.
- Do not place it in sunlight or on a radiator/heating appliance.

Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

Material composition

64% polyamide
27% elastane
5% polyester
4% wool (merino)

Latex free

Disposal

The product can be disposed of in the domestic waste.

Your medi team

Wishes you a speedy recovery!

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR). The traceability of this product is assured via a UDI code **UDI**.

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