1


5

9

13

2


6


10

14


3


7

11

15

4

8

12

16

Deutsch

## medi Soft OA

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---

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Polyamid, Polyurethan, Polyester, Lycra, Aluminium

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# medi Soft OA

## Intended purpose

medi Soft OA is a knee brace for relief and guidance of the knee joint.

## Indications

All indications that require load relief of a compartment via the 3-point pressure relief principle, additional physiological support and, where required, limitation of the range of motion of the knee joint, such as:

- Unilateral osteoarthritis
- Cartilage surgery or chondroplasty
- Meniscus intervention

## Contra-indications

Combined lateral and medial osteoarthritis

## Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in conjunction with sweating) or due to the composition of the material.

## Intended patient groups

Healthcare professionals should provide care to the adults and children, applying the available information on the measurements/sizes and necessary functions/indications and in line with the information provided by the manufacturer, and acting under their own responsibility.

## Setting the restrictions for flexion and extension (can only be carried out by an orthopaedic technician)

Extension restriction wedges, silver (stretch restriction): -5°, 0°, 10°, 20°, 30°, 45°

Flexion restriction wedges, black (bend restriction) 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°, 120°

Can be immobilised at: 0°, 10°, 20°, 30°, 45°

## Adjusting the level of flexion and extension restriction

1. Your knee brace comes with 0° extension wedges pre-installed.
2. To change the wedges, remove the brace joint from the brace.
3. Loosen the screw on the inside section and remove the screw and wedge.
4. Place the appropriate wedge for the front and/or back in the side opening and screw the wedge in tight.
5. Move the brace joint a few times to ensure that it can move freely and the wedges are secure.
6. Mount the brace joint on the brace.

Only change the extension and flexion restrictions under the orders of the doctor treating you. To avoid hyperextension, extension wedges or 0° wedges must always be used. The wedges set the exact radius of movement for the brace. If soft tissue conditions pose issues, then it may be necessary to use a larger wedge, due to

soft tissue compression, in order to limit movement of the knee to the desired radius.

### Putting on the brace

#### Before use

- Ensure that both load relief mechanisms are in the starting position (close to the knee joint) and that the switch is set to the white dot. (Fig. 1)
- Open all quick-release clasps on the front of the brace, by turning them slightly. (Fig. 2)
- Open the hook and loop fastening on the brace. (Fig. 3)

#### Putting the brace on

- Bend your leg at approx. 70° and pull the brace onto the injured leg. (Fig. 4).
- Pull the brace up your leg until the fabric comes to below the kneecap and ensure that the kneecap is positioned in the centre of the opening. The centre of the brace joint between the two lateral joint splints is now sitting around the same level as the upper edge of the kneecap.
- Please make sure that the brace has been placed on the leg straight and is not twisted. (Fig. 5)
- Now close both of the hook and loop fastenings on the upper part of the brace. (Fig. 6)
- Fasten the quick-release clasp for strap 1 first and then fasten the quick-release clasp for strap 2 (Fig. 7).
- Now fasten straps 3 and 4. (Fig. 8)
- Push and/or pull the load relief mechanism up or down, moving along the increasing number scale. (Attention: To reduce pressure, the lock needs to be released and then refastened after setting).  
An audible click is heard for every change in level (Fig. 9 and 10) to show that the brace has engaged.  
The greatest level of load relief is

attained at level 9.

Alternatively, you can open the slider, then push the mechanism to the desired position on the number scale and then lock the slider again.

The orthopaedic technician has set the brace in advance and has established a benchmark value with you. Depending on your level of pain, you can adjust the mechanism to fit your needs.

- Take a few steps and ensure that the brace is sitting correctly on the leg. If necessary, correct the position of the brace again.

#### Taking the brace off

- Deactivate the load relief mechanism by releasing the switch and pushing the mechanism towards the knee joint. (Figs. 11 and 12)
- Then set the switch back to the white dot.
- Undo all of the quick-release clasps. (Fig. 13)
- Open the hook and loop fastening on the brace and pull it off your leg, moving downwards. (Fig. 14)
- Close the velcro fastening and all quick-release clasps to ensure that no damage is done to the brace. (Fig. 15)

#### Note for technicians

- Set the load relief system (strap system) so that the patient attains a good level of pain relief. To do this, with the brace in its starting position (control is close to joint), pull the straps tightly. Then set up the three-point system. Level 7/8 is considered the guideline starting level. The patient can increase or reduce the level of pressure. (Fig. 16)
- Allow the patient to wear the brace for a few minutes to ensure that the straps are providing an appropriate and comfortable level of tension.

---

### Care instructions

Hook and loop fasteners should be closed for washing. Soap residues, lotions and ointments can cause skin irritation and material wear.

- Wash the product by hand, preferably using medi clean washing agent.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

---

### Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

---

### Material composition

Polyamide, PU, polyester, lycra, aluminum

### Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

### Disposal

The product can be disposed of in the domestic waste.

### Your medi team

Wishes you a speedy recovery!

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for

serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR). The traceability of this product is assured via a UDI code [UDI](#).

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