medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth
Germany
T +49 921 912-0
F +49 921 912-780
ortho@medi.de
www.medi.de

medi Australia Pty Ltd
 Unit 4/13-21 Hallmark Street
 Pendle Hill NSW 2145
 Australia
 T +61 2 9890 8696
 F +61 2 9890 8439
 sales@medi.au
 www.mediaustralia.com.au

medi Austria GmbH
 Adamgasse 16/7
 6020 Innsbruck
 Austria
 T +43 512 57 95 15
 F +43 512 57 95 15 45
 vertrieb@medi-austria.at
 www.medi-austria.at

medi Belgium NV
 Staatsbaan 77/0099
 3945 Ham
 Belgium
 T: +32-11 24 25 60
 F: +32-11 24 25 64
 info@medibelgium.be
 www.medibelgium.be

medi Canada Inc / médi Canada Inc
 104-1375 Lionel-Boulet,
 Varennes, Québec,
 QC Canada J3X 1P7
 T +1 450-583-3317 / +1 800 361 3153
 F +1 888-583-6827
 service@medicanada.ca
 www.medicana.ca

MAXIS a.s.,
 medi group company
 Slezská 2127/13
 120 00 Prague 2
 Czech Republic
 T: +420 571 633 510
 F: +420 571 616 271
 info@maxis-medica.com
 www.maxis-medica.com

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medi Trading (Shanghai) Co. Ltd.
 Room 613, HnR Group Building
 No.1178, Tian Yao Qiao Road
 Xuhui District, Shanghai, 200232
 People's Republic of China
 T +86-21 50582318
 infocn@medi.cn
 www.medi.cn

UK Responsible Person
 medi UK Ltd.
 Plough Lane
 Hereford HR4 OEL
 Great Britain
 T +44-1432 37 35 00
 F +44-1432 37 35 10
 enquiries@mediuk.co.uk
 www.mediuk.co.uk

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Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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Polyamid, Polyurethan, Polyester, Lycra, Aluminium

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medi Soft OA light

Intended purpose

medi Soft OA light is a knee brace for relief of the knee joint .

Indications

All indications requiring unicompartimental relief by 3 point principle and additional physiological guidance such as:

- Unilateral osteoarthritis
- Meniscus intervention

Contra-indications

Combined lateral and medial osteoarthritis

Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in conjunction with sweating) or due to the composition of the material. The immobilisation which is achievable by the device presents a risk for development of thrombosis in the lower

veins of the legs, especially in the context of surgical interventions or trauma.

Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Setting the restrictions for flexion and extension (can only be carried out by an orthopaedic technician)

Extension restriction wedges, silver (stretch restriction): -5°, 0°, 10°, 20°, 30°, 45°

Flexion restriction wedges, black (bend restriction) 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°, 120°

Can be immobilised at: 0°, 10°, 20°, 30°, 45°

Adjusting the level of flexion and extension restriction

1. Your knee brace comes with 0° extension wedges pre-installed.
2. To change the wedges, remove the brace joint from the brace.
3. Loosen the screw on the inside section and remove the screw and wedge.
4. Place the appropriate wedge for the front and/or back in the side opening and screw the wedge in tight.
5. Move the brace joint a few times to ensure that it can move freely and the wedges are secure.

6. Mount the brace joint on the brace.

Only change the extension and flexion restrictions under the orders of the doctor treating you. To avoid hyperextension, extension wedges or 0° wedges must always be used. The wedges set the exact radius of movement for the brace. If soft tissue conditions pose issues, then it may be necessary to use a larger wedge, due to soft tissue compression, in order to limit movement of the knee to the desired radius.

Fitting instructions

Open the straps on the front and the hook and loop fastenings on the brace. (Fig. 1)

Putting the brace on

- Bend your leg at approx. 70° and pull the brace onto the injured leg. (Fig. 2)
- Pull the brace up your leg until the fabric comes to below the kneecap and ensure that the kneecap is positioned in the centre of the opening. The centre of the brace joint between the two lateral joint splints is now sitting around the same level as the upper edge of the kneecap. Please make sure that the brace has been placed on the leg straight and is not twisted. (Fig. 3)
- Now close the fasteners on the upper part of the brace. (Fig. 4)
- Once again, ensure that the brace is sitting correctly.
- First, pull the lower strap (1) through the slot on the joint and set it to the setting that was defined by the technician. Fasten the hook and loop. Then set the upper strap (2) and fasten it.
- The straps at the back have already been set by the technician. (Figs. 5 and 6)
- Take a few steps and ensure that the

brace is sitting correctly on the leg. If necessary, correct the position of the brace again. (Fig. 7)

Taking the brace off

- Open both of the straps and stick the end of the straps back onto their respective straps, to prevent them from snagging clothing. (Fig. 8)
- Open the fasteners on the brace and pull it off your leg, moving downwards. Close the hook and loop fastening to prevent damage from occurring. (Fig 9)

Note for technicians

- Set the load relief system (strap system) so that the patient attains a good level of pain relief. The base setting should be set so that the patient can still increase or decrease pressure. (Fig. 10)
- Allow the patient to wear the brace for a few minutes to ensure that the straps are providing an appropriate and comfortable level of tension.

Care instructions

Hook and loop fasteners should be closed for washing. Soap residues, lotions and ointments can cause skin irritation and material wear.

- Wash the product by hand, preferably using medi clean washing agent.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

Material composition

Polyamide, poliurethane, polyester, lycra, aluminum

Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

Disposal

Dispose of properly after use.

**Your medi team**

Wishes you a speedy recovery!

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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