**Hersteller / Manufacturer**  
**medi GmbH & Co. KG**  
**Medicusstraße 1**  
**D-95448 Bayreuth**  
**Germany**  
**T +49 921 912-0**  
**F +49 921 912-780**  
**ortho@medi.de**  
**www.medi.de**

medi Australia Pty Ltd  
Unit 4/13-21 Hallmark Street  
Pendle Hill NSW 2145  
Australia  
T +61 2 9890 8696  
F +61 2 9890 8439  
sales@medi.au  
www.medi.australia.com.au

medi Austria GmbH  
Adamgasse 16/7  
6020 Innsbruck  
Austria  
T +43 512 57 95 15  
F +43 512 57 95 15 45  
vertrieb@medi-austria.at  
www.medi-austria.at

medi Belgium NV  
Staatsbaan 77/0099  
3945 Ham  
Belgium  
T: + 32-11 24 25 60  
F: +32-11 24 25 64  
info@medibelgium.be  
www.medibelgium.be

medi Canada Inc / médi Canada Inc  
104-1375 Lionel-Boulet,  
Varenes, Québec,  
QC Canada J3X 1P7  
T +1 450-583-3317 / +1 800 361 3153  
F +1 888-583-6827  
service@medicanada.ca  
www.medicanada.ca

MAXIS a.s.,  
medi group company  
Slezská 2127/13  
120 00 Prague 2  
Czech Republic  
T: +420 571 633 510  
F: +420 571 616 271  
info@maxis-medica.com  
www.maxis-medica.com

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medi Trading (Shanghai) Co. Ltd.  
Room 1703, T1, CES WEST BUND CENTER,  
No. 277, Long Lan Road  
Xuhui District  
200232 Shanghai  
People's Republic of China  
T: +86-21 50582319  
F: +86-21 50582319

UK Responsible Person  
medi UK Ltd.  
Plough Lane  
Hereford HR4 OEL  
Great Britain  
T +44-1432 37 35 00  
F +44-1432 37 35 10  
enquiries@mediuk.co.uk  
www.mediuk.co.uk

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#### Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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# medi ROM Kidz

## Intended Purpose

The medi ROM Kidz is a universal knee splint for mobilisation .

## Indications

All indications which require early functional mobilisation and limitation of the range of motion of the knee joint, such as:

- Following ligament injury (post-traumatic/post-operative)
- Following meniscal injury (post-traumatic/post-operative)
- Following tendon injury (post-traumatic/post-operative)

## Contraindications

None known at present.

## Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paresthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in conjunction with sweating) or due to the composition of the material. The immobilisation achieved by the device presents a risk for development of

thrombosis in the lower veins of the legs, especially in the context of surgical interventions or trauma. The product is not suitable for children under 6 years of age. Danger of suffocation from small parts that can be swallowed.

## Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

## Adjusting the flexion and extension stops.

- To set the radius of movement, first loosen the screw (Fig. 1). Make sure not to let the spacer shim become lost.
- Then turn the wedge-shaped plastic opening until it aligns with one of the pins (Fig. 2).
- Now you can take the pin out, place it in a different position and thus set the desired degree on the brace. Repeat this for the other articulation (Fig. 3).
- Replace the spacer disc and turn the round opening back to its original position (Fig. 4).
- Finally, tighten the safety screw (Fig. 5)

Limitation of extension:

0°, 15°, 30°, 45°, 60°

Limitation of flexion:

0°, 15°, 30°, 45°, 60°, 75°, 90°, 105°

Immobilisation: 0°, 15°, 30°, 45°, 60°

### Fitting instructions

Open all straps and loosen the joint splints from the pad.

Wrap the pad for the lower and upper leg around the leg tightly and attach the pad to the front of the leg (Fig. 6).

First, fasten the lateral (outer) splint to the pad (Fig. 7). Ensure that the ROM articulation is aligned on the patella and the mid-line of the knee. Now fasten the medial (inner) splint, using the same method. Ensure that both of the articulations are level with each other when doing this. Now close the top lower leg strap (Fig. 8) first and then close the bottom upper leg strap (Fig. 9). Then close the rest of the straps, working from the bottom of the brace up.

To take off the brace, you simply need to open the straps and the pad and pull the brace backwards.

### Important notes

Only change the extension and flexion limits on instruction from your doctor or orthotist.

### Care instructions

Hook and loop fasteners should be closed for washing. Soap residues, lotions and ointments can cause skin irritation and material wear.

- Wash the paddings by hand, preferably using medi clean washing agent.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

### Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

### Material composition

Aluminum, PU foam, cotton

### Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

### Disposal

Dispose of properly after use.

### Your medi Team

Wishes you a quick return to full fitness.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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