**medi GmbH & Co. KG**  
**Medicusstraße 1**  
**95448 Bayreuth**  
**Germany**  
**T +49 921 912-0**  
**www.medi.de**

medi Australia Pty Ltd  
 83 Fennell Street  
 North Parramatta NSW 2151  
 Australia  
 T +61-2 9890 8696  
 F +61-2 9890 8439  
 sales@mediaustralia.com.au  
 www.mediaustralia.com.au

medi Austria GmbH  
 Adamgasse 16/7  
 6020 Innsbruck  
 Austria  
 T +43 512 57 95 15  
 F +43 512 57 95 15 45  
 vertrieb@medi-austria.at  
 www.medi-austria.at

medi Belgium NV  
 Staatsbaan 77/0099  
 3945 Ham  
 Belgium  
 T: +32-11 24 25 60  
 F: +32-11 24 25 64  
 info@medibelgium.be  
 www.medibelgium.be

medi Canada Inc / médi Canada Inc  
 104-1375 Lionel-Boulet,  
 Varennes, Québec,  
 QC Canada J3X 1P7  
 T +1 450-583-3317 / +1 800 361 3153  
 F +1 888-583-6827  
 service@medicanada.ca  
 www.medicanada.ca

MAXIS a.s.,  
 medi group company  
 Slezská 2127/13  
 120 00 Prague 2  
 Czech Republic  
 T: +420 571 633 510  
 F: +420 571 616 271  
 info@maxis-medica.com  
 www.maxis-medica.com

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medi Trading (Shanghai) Co. Ltd.  
 Room 1703, T1, CES WEST BUND CENTER,  
 No. 277, Long Lan Road  
 Xuhui District  
 200232 Shanghai  
 People's Republic of China  
 T: +86-21 50582319  
 F: +86-21 50582319

UK Responsible Person  
 medi UK Ltd.  
 Plough Lane  
 Hereford HR4 0EL  
 Great Britain  
 T +44-1432 37 35 00  
 F +44-1432 37 35 10  
 enquiries@mediuk.co.uk  
 www.mediuk.co.uk

## circaid® juxtafit® lower leg inelastic compression system

Gebrauchsanweisung. Instructions for use. Mode d'emploi. Instrucciones de uso. Instruções para aplicação. Istruzioni per l'uso. Gebruiksaanwijzing. Bruksanvisning. Kasutusjuhend. Käyttöohje. 使用説明書.

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### Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If undue pain or an unpleasant sensation develops while you are wearing the product, please consult your doctor or orthotist immediately. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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circaid juxtafit essentials 6 Monate

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53% Polyurethan, 43% Nylon, 4% Elastan

circaid juxtafit premium Unterschenkel

48% Polyurethan, 46% Nylon, 3% Elastan, 3% Polyethylen

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## circaid® juxtafit® lower leg

#### Intended Purpose

The compression system is designed to provide compression to the leg for patients with venous and lymphatic disorders. The device consists of a series of juxtaapositioned, inelastic bands extending from the central portion of the garment. The wrap contours the leg and is designed to cover from just above the ankle bone to just below the knee.

#### Indications

- Lymphedema
- Lipedema
- Chronic venous insufficiency
  - Varicose veins
  - Lipodermatosclerosis
- Venous leg ulcer
- Deep vein thrombosis/Thrombosis Prevention
- Acute deep vein thrombosis
- Post Thrombotic Syndrome
- Post Sclerotherapy
- Dependent oedema

#### Contraindications

- Severe Peripheral Arterial Disease
- Uncontrolled Congestive Heart Failure
- Septic Phlebitis
- Phlegmasia Cerulea Dolens
- Untreated infection
- Any circumstance where increased venous and lymphatic return is undesirable

#### Risks / Side effects

- Mild to Moderate Peripheral Arterial Disease
- Impaired Sensation - must be able to sense applied pressure
- Intolerance to material(s)
- Pediatric & assisted use - must be able to sense & communicate applied pressure

#### Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for

this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

**Note:** *Unless your treating physician orders or recommends otherwise, we recommend wearing the leg fitting together with the respective foot options enclosed with the set. In case of doubt, please consult your attending physician.*

#### **Wearing recommendation**

Always ensure that the hook tabs are placed on the corresponding landing pads.

Always ensure that the correct prescribed pressure range is being applied. Slightly loosen the bands of the garment for night-time wear. Each vertical BPS line on the bands may be at a different location on the leg depending on the size and shape of the leg. They do not need to be aligned down the front center of the leg or in any other specific location on the leg. If experiencing any pain, or discomfort, reduce the compression (mmHg) applied or remove device. Inform your clinician of your action. When bands are placed correctly they should slightly overlap, leaving no space between bands. Do not tuck bands. The Hook tabs should never touch the skin!

#### **Donning instructions**

Step 1: Slide the compressive undersock onto the leg.

Ensure that there are no wrinkles in the undersock.

Step 2: Align the garment just above the ankle bone and below the knee crease. The black side of the material should face the skin.

Unroll the top two bands and loosely secure them to the outside of the garment to anchor it in place. This enables you to easily raise or lower the garment for correct positioning.

**Note:** *The bottom most band can be identified by the care label.*

Step 3: Unroll the bottom two bands and secure the bottom most band to the outside of the garment to a firm and comfortable compression level.

Step 4: While holding the second band, detach the next band. Secure the second band and continue this process with the other bands, going up the leg.

Step 5: The garment should lay flat and wrinkle-free against the leg. Adjust each band to a firm and comfortable compression level.

Ensure that the bands are placed in an alternating order.

Step 6: Once all bands are secured, inspect the garment for any gaps or creases. Adjust the bands as necessary.

Step 7: Fold the top end of the compressive undersock over the lower leg compression garment.

#### **Doffing instructions**

To remove the legging, first remove the cover up, if using. Then remove the knee cover, shelf straps, and/or foot compression, if using any of these. Then unfold the circaid undersleeve. Detach all bands of the legging starting with the top band. Fold each band onto itself to keep the hook and loop tabs lint-free. Remove the compression garment and

then remove the compressive undersock.

#### **Using the Built-In-Pressure System (BPS)**

This garment is designed to provide gradient compression based on these average ankle compression levels.

Step 1: Locate the Built-In-Pressure card in your packaging.  
Step 2: Identify your prescribed pressure scale using the color-coded system.

Step 3: Starting with the bottom band, line up the black triangle on the card with one of the BPS lines on the bottom band.

Step 4: Note where the second BPS line lines up with the card's compression ranges (A).

Step 5: If the BPS line on the garment either falls short (B) or goes beyond the correct compression range (C), readjust the band as necessary so that the second BPS line on the garment is aligned with the prescribed compression range (A). Completely loosen the band before reapplying it if too much tension has been applied.

Step 6: Repeat the steps 3 - 5 for each band going up the leg. Be sure to adjust all bands to the same tension to ensure gradient compression.

Bands may need adjusting throughout the day to maintain the prescribed compression range as edema is reduced.

**Note:** *Please note it is recommended that the limb is always treated in its entirety with medical compression. This means that compression products for the leg begin at the foot or the ankle and extend up to below the hollow of the knee or beyond to include the thigh.*

#### **Care instructions**

The garment can be machine washed using a gentle cycle and dried on a low heat setting. We recommend using a laundry bag. To extend the life of the garment, hand wash and drip dry.

- 
-  Machine wash warm
  -  Do not bleach
  -  Tumble dry, low heat
  -  Do not iron
  -  Do not dry clean
  -  Do not wring
- 

#### **Storage instructions**

Keep the product in a dry place and do not expose to direct sunlight.

#### **Lifespan**

circaid juxtafit essentials 6 months.

circaid juxtafit premium 12 months.

Due to material wear-and-tear the medical efficacy can only be expected for a defined lifespan. This assumes correct handling of the garment (e.g. proper care, application, and removal).

#### **Material composition**

circaid juxtafit essentials lower leg

53% Polyurethane, 43% Nylon, 4% Elastane,

circaid juxtafit premium lower leg

48% Polyurethane, 46% Nylon, 3% Elastane, 3% Polyéthylène

### Disposal

The product can be disposed of in the domestic waste.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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