**medi GmbH & Co. KG**  
 Medicusstraße 1  
 95448 Bayreuth  
 Germany  
 T +49 921 912-0  
 www.medi.de

medi Australia Pty Ltd  
 Unit 4/13-21 Hallmark Street  
 Pendle Hill NSW 2145  
 Australia  
 T +61 2 9890 8696  
 F +61 2 9890 8439  
 sales@medi.au  
 www.mediaustralia.com.au

medi Austria GmbH  
 Adamgasse 16/7  
 6020 Innsbruck  
 Austria  
 T +43 512 57 95 15  
 F +43 512 57 95 15 45  
 vertrieb@medi-austria.at  
 www.medi-austria.at

medi Belgium NV  
 Staatsbaan 77/0099  
 3945 Ham  
 Belgium  
 T: +32-11 24 25 60  
 F: +32-11 24 25 64  
 info@medibelgium.be  
 www.medibelgium.be

medi Canada Inc / médi Canada Inc  
 104-1375 Lionel-Boulet,  
 Varennes, Québec,  
 QC Canada J3X 1P7  
 T +1 450-583-3317 / +1 800 361 3153  
 F +1 888-583-6827  
 service@medicanada.ca  
 www.medicanada.ca

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## Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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# circaid® compressive undersock

## Intended purpose

The compression system is designed to provide graduated compression to the foot and ankle region for patients with venous and lymphatic disorders.

## Indications

- Lymphedema
- Other forms of edema
- Lipedema
- Chronic venous insufficiency
- Varicose veins
- DVT/Thrombosis Prevention
- Acute deep vein thrombosis
- Post Thrombotic Syndrome
- Venous Stasis Ulcer
- Post Sclerotherapy
- Lipodermatosclerosis

## Contraindications

- Advanced peripheral arterial disease (if one of these parameters is present: ABPI < 0.5; ankle arterial pressure < 60mmHg; toe pressure < 30 mmHg or TcPO<sub>2</sub> < 20 mmHg on dorsum of foot). If non-elastic materials are used, the patient can still try out compression clothing if they have ankle arterial pressure between 50 and 60 mmHg, under close clinical supervision
- Decompensated heart failure
- Septic phlebitis
- Phlegmasia cerulea dolens
- Untreated infection of the leg and/or foot

- Any circumstances under which an increase in venous or lymphatic reflux is undesirable
- Suspected or known untreated acute deep vein thrombosis

## Risks/side effects

- Mild to Moderate Peripheral Arterial Disease
- Impaired Sensation - must be able to sense applied pressure
- Intolerance to material(s)
- Pediatric & assisted use - must be able to sense & communicate applied pressure

## Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

## Wearing recommendation

*Note: Unless your treating physician orders or recommends otherwise, we recommend wearing the the legwear together with the respective footwear included in the set. If in doubt, please consult your treating physician.*

The circaid compression garment should never hurt. It should feel firm but comfortable. The undersock may be worn day and

night, but removal is recommended if experiencing discomfort when lying down.

### Donning instructions

- Slide the compressive undersock onto the leg.
- Ensure that there are no wrinkles.
- Don the lower leg compression garment.
- Fold the top end of the compressive undersock over the lower leg compression garment

### Care instructions

The garment can be machine washed using a gentle cycle and dried on a low heat setting. To extend the life of the garment, hand wash and drip dry.

---

 Machine wash warm

 Do not bleach

 Tumble dry, low heat

 Do not iron

 Do not dry clean

---

### Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

### Lifespan

6 months – due to material wear and/or slackening, medical efficacy can only be guaranteed for a defined lifespan. This assumes correct handling of the product (e.g. during care and when putting on and taking off the product).

### Material composition

Polyester, nylon, elastane

### Disposal

Dispose of properly after use.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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Tachezy Sanit Srl.