medi GmbH & Co. KG
Medicusstraße 1
D-95448 Bayreuth
Germany
T +49 921 912-0
www.medi.de

medi Australia Pty Ltd
83 Fennell Street
North Parramatta NSW 2151
Australia
T +61 2 9890 8696
F +61 2 9890 8439
sales@mediaustralia.com.au
www.mediaustralia.com.au

medi Austria GmbH
Adamgasse 16/7
6020 Innsbruck
Austria
T +43-512 57 95 15
F +43-512 57 95 15 45
vertrieb@medi-austria.at
www.medi-austria.at

medi Belgium NV
Posthoornstraat 13/1
3582 Koersel
Belgium
T: +32-11 24 25 60
F: +32-11 24 25 64
info@medibelgium.be
www.medibelgium.be

medi Brasil
Rua Cristovan de Vita, n.º 260, Un. 30,
Centro Logístico Raposo Tavares,
Vila Camargo
Vargem Grande Paulista - SP
CEP 06730-000
Brasil
T: +55-11-3500 8005
sac@medibrasil.com
www.medibrasil.com

medi Canada Inc / médi Canada Inc
597, Rue Duvernay, Verchères
QC Canada J0L 2R0
T: +1 450-583-3317 / +1 800 361 3153
F: +1 888-583-6827
service@medicanada.ca
www.medicanada.ca

MAXIS a.s.,
medi group company
Slezská 2127/13
120 00 Prague 2
Czech Republic
T: +420 571 633 510
F: +420 571 616 271
info@maxis-medica.com
www.maxis-medica.com

medi Danmark ApS
Vejlegardsvej 59
2665 Vallensbaek Strand
Denmark
T +45 70 25 56 10
F +45 70 25 56 20
kundeservice@sw.dk
www.medit danmark.dk

medi Bayreuth España SL
C/Canigó 2 – 6 bajos
Hospitalet de Llobregat
08901 Barcelona
Spain
T +34 93 260 04 00
F +34 93 260 23 14
medi@mediespana.com
www.mediespana.com

medi France
Z.I. Charles de Gaulle
25, rue Henri Farman
93297 Tremblay en France Cedex
France
T +33 1 48 61 76 10
F +33 1 49 63 33 05
infos@medi-france.com
www.medi-france.com

medi Hungary Kft.
Bokor u. 21.
1037 Budapest
Hungary
T +36 1371 0090
F +36 1371 0091
info@medi.hu
www.medi.hu

medi

circaid® juxtalite® lower leg inelastic compression system

Instructions for use. Gebrauchsanweisung. Instruções para aplicação. Instrucciones de uso. Gebruiksaanwijzing. Istruzioni per l'uso. Bruksanvisning. Brugsvejledning. مداخلت سالات اداش رر. Інструкція з використання. Instrukcja użytkowania. תוארה שומיש.

circaid juxtalite lower leg

circaid juxtalite ankle foot wrap

circaid juxtalite with ankle foot wrap


EO15977/03.2020


4 061736 636020

medi. I feel better.

5a

5b

5c

5d


5e

5f

5g

5h

5i


5j


6a

6b

6c

English

circaid® juxtalite® lower leg

Purpose

The compression system is designed to provide compression to the lower leg for patients with venous disorders. The device consists of a series of juxtapositioned, inelastic bands extending from the central portion of the garment. The wrap contours the leg and is designed to cover from just above the ankle bone to just below the knee.

Indications

- Chronic venous insufficiency
- Venous Stasis Ulcer
- DVT/Thrombosis Prevention
- Post Sclerotherapy
- Post thrombotic syndrome
- Dependent edema
- Varicose veins

Contraindications

- Severe Peripheral Arterial Disease
- Uncontrolled Congestive Heart Failure
- Septic Phlebitis
- Phlegmasia Cerulea Dolens
- Untreated infection
- Any circumstance where increased venous and lymphatic return is undesirable
- Suspected or known untreated acute DVT

Cautions

- Mild to Moderate Peripheral Arterial Disease
- Impaired Sensation - must be able to sense applied pressure
- Intolerance to material(s)
- Pediatric & assisted use - must

be able to sense & communicate applied pressure

This Package contains

- 1 circaid juxtalite lower leg
- 2 circaid compressive undersocks
- 1 Built-In-Pressure System™ (BPS™) card
- 1 booklet instructions for use

Storage instructions

Store in a dry place

Keep out of sunlight

Lifespan

circaid juxtalite has a lifespan of 6 months.

Due to material wear-and-tear the medical efficacy can only be expected for a defined lifespan. This assumes correct handling of the garment (e.g. proper care, application, and removal).

Disposal

Please dispose of the garment in household waste. No special criteria for disposal apply for medical compression garments.

medi Italia s.r.l.
Via Giuseppe Ghedini, 2
40069 Zola Predosa (BO)
Italy
T: +39-051 6 13 24 84
F: +39-051 6 13 29 56
info@medi-italia.it
www.medi-italia.it

medi Japan K.K.
5-1 Nihonbashi Kabutocho Chuoku
103-0026 Tokyo
Japan
T: +81 3 6778 2590
F: +81 3 5847 7901
info@medi-japan.jp
www.medi-japan.co.jp

Tachezy Sanit Srl.
Via Druso 48/a
39100 Bolzano
Italy
T: +39 0471 288 068
F: +39 0471 282 773
info@tachezysanit.com
www.tachezysanit.com

medi Nederland BV
Heusing 5
4817 ZB Breda
The Netherlands
T +31 76 57 22 555
F +31 76 57 22 565
info@medi.nl
www.medi.nl

medi Norway AS
Vestheimvegen 35
4250 Kopervik
Norway
T: +47 52 84 45 00
F: +47 52 84 45 09
post@medinorway.no
www.medinorway.no

medi Polska Sp. z o.o.
ul. Łabędzka 22
44-121 Gliwice
Poland
T: +48-32 230 60 21
F: +48-32 202 87 56
info@medi-polska.pl
www.medi-polska.pl

medi Bayreuth Unipessoal Lda.
Rua do Centro Cultural, no. 43
1700-106 Lisbon
Portugal
T +351 21 843 71 60
F +351 21 847 08 33
medi.portugal@medibayreuth.pt
www.medi.pt

medi RUS LLC
Business Center NEO GEO
Butlerova Street 17
117342 Moscow
Russia
T +7-495 374 04 56
F +7-495 374 04 56
info@medirus.ru
www.medirus.ru

medi Trading (Shanghai) Co. Ltd.
Room 1703, T1, CES WEST BUND CENTER,
No. 277, Long Lan Road
Xuhui District
200232 Shanghai
People's Republic of China
T: +86-21 50582319
F: +86-21 50582319

medi Sweden AB
Box 6034
192 06 Sollentuna
Sweden
T +46 8 9697 98
F +46 8 626 68 70
info@medi.se
www.medi.se

Medi Turk Ortopedi Medikal İc ve Dis Tic. AS
Mustafa Kemal Mah. 2156 Sokak No: 16 /4-5
06510 Çankaya Ankara
Turkey
T: +90 312 435 20 26
F: +90 312 434 17 67
info@medi-turk.com
www.medi-turk.com

medi UK Ltd.
Plough Lane
Hereford HR4 OEL
Great Britain
T +44 1432 37 35 00
F +44 1432 37 35 10
enquiries@mediuk.co.uk
www.mediuk.co.uk

medi Ukraine LLC
Evhena Sverstiuka str, 11.
Kiev 02002
Ukraine
T: +380 44 591 11 63
F: +380 44 392 73 73
info@medi.ua
http://medi.ua

medi USA L.P.
6481 Franz Warner Parkway
Whitsett, N.C. 27377-3000
USA
T +1-336 4 49 44 40
F +1-888 5 70 45 54
info@mediusa.com
www.mediusa.com

Material composition

circaid compression garment

- Nylon 56%
- Polyurethane 40%
- Elastane 4%

For all other components, please refer to the inserts included with each item.

This garment does not contain latex.

Care instructions

The garment can be machine washed using a gentle cycle and dried on a low heat setting. If doing so, it is recommended that you place the garment in a pillow case to aid in keeping the Velcro tab material lint-free and to protect other clothing from getting damaged.

To extend the life of the garment, hand wash and drip dry.

Machine wash warm

Do not bleach

Tumble dry, low heat

Do not iron

Do not dry clean

Do not wring

Wearing recommendation

Always ensure that the correct prescribed pressure range is being applied.

Slightly loosen the bands of the garment for night-time wear.

Each vertical BPS line on the bands may be at a different location on the leg depending on the size and shape of the leg. They do not need to be aligned down the front center of the leg

or in any other specific location on the leg.

If experiencing any pain, or discomfort, reduce the compression (mmHg) applied or remove device. Inform your clinician of your action.

When bands are placed correctly they should slightly overlap, leaving no space between bands.

Do not tuck bands. The Velcro tabs should never touch the skin!

Donning instructions

Step 1: Slide the compressive undersock onto the leg. Ensure that there are no wrinkles in the undersleeve.

Step 2: Align the garment just above the ankle bone and below the knee crease. The black side of the material should face the skin. Unroll the top two bands and loosely secure them to the outside of the garment to anchor the legging in place. This enables you to easily raise or lower the garment for correct positioning.

Note:

The bottom most band can be identified by the care label.

Step 3: Unroll the bottom two bands and secure the bottom most band to the outside of the legging to a firm and comfortable compression level.

Step 4: While holding the second band, detach the next band.

Secure the second band and continue this process with the other bands, going up the leg.

Step 5: The garment should lay flat and wrinkle free against the leg. Adjust each band to a firm and comfortable compression

level. Ensure that the bands are placed in an alternating order.

Step 6: Once all bands are secured, inspect the garment for any gaps or creases. Adjust the bands as necessary.

Step 7: Fold the top end of the compressive undersock over the lower leg compression garment

Doffing instructions

Detach all bands of the legging starting with the top band. Fold each band onto itself to keep the hook and loop tabs lint-free.

Remove the compression garment and then remove the compressive undersock.

Using the Built-In-Pressure card:

Step 1: Locate the Built-In-Pressure card in your packaging.

Step 2: Identify your prescribed pressure scale using the color-coded system.

Step 3: Starting with the bottom band, line up the vertical line on the card with one of the BPS lines on the bottom band.

Step 4: Note where the second BPS line lines up with the card's compression ranges (6a).

Step 5: If the BPS line on the garment either falls short (6b) or goes beyond the correct compression range (6c), readjust the band as necessary so that the second BPS line on the garment is aligned with the prescribed compression range (6a).

Completely loosen the band before reapplying it if too much tension has been applied.

Step 6: Repeat the steps 3 - 5 for each band going up the leg. Be sure to adjust all bands to the same compression level. Bands may need adjusting throughout the day to maintain the prescribed compression range.

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