medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth
Germany
T +49 921 912-0
www.medi.de

medi Australia Pty Ltd
Unit 4/13-21 Hallmark Street
Pendle Hill NSW 2145
Australia
T +61 2 9890 8696
F +61 2 9890 8439
sales@medi.au
www.mediaustralia.com.au

medi Austria GmbH
Adamgasse 16/7
6020 Innsbruck
Austria
T +43-512 57 95 15
F +43-512 57 95 15 45
vertrieb@medi-austria.at
www.medi-austria.at

medi Belgium NV
Staatsbaan 77/0099
3945 Ham
Belgium
T: +32-11 24 25 60
F: +32-11 24 25 64
info@medibelgium.be
www.medibelgium.be

medi Canada Inc / médi Canada Inc
104-1375 Lionel-Boulet,
Varennes, Québec,
QC Canada J3X 1P7
T +1 450-583-3317 / +1 800 361 3153
F +1 888-583-6827
service@medicanada.ca
www.medicanada.ca

MAXIS a.s.,
medi group company
Slezská 2127/13
120 00 Prague 2
Czech Republic
T: +420 571 633 510
F: +420 571 616 271
info@maxis-medica.com
www.maxis-medica.com

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medi

circaid® juxtalite® lower leg inelastic compression system

Gebrauchsanweisung. Instructions for use. Mode d'emploi.
Instrucciones de uso. Instruções para aplicação. Istruzioni per
l'uso. Gebruiksaanwijzing. Bruksanvisning. Uputa za uporabu.
Instrukcja zakładania. Használati útmutató. Bruksanvisning.
Kasutusjuhend. Käyttöohje.

EO15977 / 01.2025

medi. I feel better.

circaid juxtalite lower leg

circaid juxtalite ankle foot wrap

circaid juxtalite with ankle foot
wrap

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Important notes
This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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



medi Norway AS
Vestheimvegen 35
4250 Kopervik
Norway
T: +47 52 84 45 00
F: +47 52 84 45 09
post@medinorway.no
www.medinorway.no

medi Polska Sp. z o.o.
ul. Łąbiedzka 22
44-121 Gliwice
Poland
T: +48-32 230 60 21
F: +48-32 202 87 56
info@medi-polska.pl
www.medi-polska.pl

medi Bayreuth Unipessoal Lda.
Rua do Centro Cultural, no. 43
1700-106 Lisbon
Portugal
T +351 21 843 71 60
F +351 21 847 08 33
medi.portugal@medibayreuth.pt
www.medi.pt

medi RUS LLC
Business Center NEO GEO
Butlerova Street 17
117342 Moscow
Russia
T +7-495 374 04 56
F +7-495 374 04 56
info@medirus.ru
www.medirus.ru

medi Trading (Shanghai) Co. Ltd.
Room 1703, T1, CES WEST BUND CENTER,
No. 277, Long Lan Road
Xuhui District
200232 Shanghai
People's Republic of China
T: +86-21 50582319
F: +86-21 50582319

medi Sweden AB
Box 6034
192 06 Sollentuna
Sweden
T +46 8 9697 98
F +46 8 626 68 70
info@medi.se
www.medi.se

Medi Turk Ortopedi Medikal İc ve Dis Tic. AS
Mustafa Kemal Mah. 2156 Sokak No: 16 /4-5
06510 Çankaya Ankara
Turkey
T: +90 312 435 20 26
F: +90 312 434 17 67
info@medi-turk.com
www.medi-turk.com

medi Ukraine LLC
Evhena Sverstiuka str, 11.
Kiev 02002
Ukraine
T: +380 44 591 11 63
F: +380 44 392 73 73
info@medi.ua
http://medi.ua

medi USA L.P.
6481 Franz Warner Parkway
Whitsett, N.C. 27377-3000
USA
T +1-336 4 49 44 40
F +1-888 5 70 45 54
info@mediusa.com
www.mediusa.com

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lower leg

Intended purpose

The compression system is designed to provide compression  to the lower leg for patients with venous disorders.

Indications

- Chronic venous insufficiency
 - Varicose Veins
 - Lipodermatosclerosis
 - Venous Stasis Ulcer
- Deep vein thrombosis/Thrombosis Prevention
- Acute DVT
- Post Thrombotic Syndrome
- Post Sclerotherapy
- Dependent oedema

Contraindications

- Advanced peripheral arterial disease (if one of these parameters is present: ABPI < 0.5; ankle arterial pressure < 60mmHg; toe pressure < 30 mmHg or TcPO₂ < 20 mmHg on dorsum of foot). If non-elastic materials are used, the patient can still try out compression clothing if they have ankle arterial pressure between 50 and 60 mmHg, under close clinical supervision
- Decompensated heart failure
- Septic phlebitis
- Phlegmasia cerulea dolens
- Untreated infection of the leg and/or foot
- Any circumstances under which an increase in venous or lymphatic reflux is undesirable
- Suspected or known untreated acute deep vein thrombosis

Risks / Side effects

- Mild to Moderate Peripheral Arterial Disease
- Impaired Sensation - must be able to sense applied pressure
- Intolerance to material(s)
- Pediatric & assisted use - must be able to sense & communicate applied pressure

Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to

adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Notes on wearing

Always ensure that the correct prescribed pressure range is being applied. Slightly loosen the bands of the garment for night-time wear. Each vertical BPS line on the bands may be at a different location on the leg depending on the size and shape of the leg. They do not need to be aligned down the front center of the leg or in any other specific location on the leg. If experiencing any pain, or discomfort, reduce the compression (mmHg) applied or remove device. Inform your clinician of your action. When bands are placed correctly they should slightly overlap, leaving no space between bands. Do not tuck bands. The Velcro tabs should never touch the skin!

Note: Unless your treating physician orders or recommends otherwise, we recommend wearing the leg fitting together with the respective foot options enclosed with the set. In case of doubt, please consult your attending physician.

Donning instructions

- Slide the compressive undersock onto the leg. Ensure that there are no wrinkles in the undersleeve.
- Align the garment just above the ankle bone and below the knee crease. The black side of the material should face the skin. Unroll the top two bands and loosely secure them to the outside of the garment to anchor the legging in place. This enables you to easily raise or lower the garment for correct positioning.

Note:

The bottom most band can be identified by the care label.

- Unroll the bottom two bands and secure the bottom most band to the outside of the legging to a firm and comfortable compression level.
- While holding the second band, detach the next band. Secure the second band and continue this process with the other bands, going up the leg.
- The garment should lay flat and

wrinkle free against the leg. Adjust each band to a firm and comfortable compression level. Ensure that the bands are placed in an alternating order.

- Once all bands are secured, inspect the garment for any gaps or creases. Adjust the bands as necessary.
- Fold the top end of the compressive undersock over the lower leg compression garment

Using the Built-In-Pressure card

This garment is designed to provide gradient compression based on these average ankle compression levels.

- Locate the Built-In-Pressure card in your packaging.
- Identify your prescribed pressure scale using the color-coded system.
- Starting with the bottom band, line up the vertical line on the card with one of the BPS lines on the bottom band.
- Note where the second BPS line lines up with the card's compression ranges (6a).
- If the BPS line on the garment either falls short (6b) or goes beyond the correct compression range (6c), readjust the band as necessary so that the second BPS line on the garment is aligned with the prescribed compression range (6a). Completely loosen the band before reapplying it if too much tension has been applied.
- Repeat the steps 3 - 5 for each band going up the leg. Be sure to adjust all bands to the same compression level. Bands may need adjusting throughout the day to maintain the prescribed compression range.

Care instructions

The garment can be machine washed using a gentle cycle and dried on a low heat setting. We recommend using a laundry bag. To extend the life of the garment, hand wash and drip dry.

 Machine wash warm

 Do not bleach

 Tumble dry, low heat

 Do not iron

 Do not dry clean

 Do not wring

Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

Lifespan

6 months - due to material wear-and-tear the medical efficacy can only be expected for a defined lifespan. This assumes correct handling of the garment (e.g. proper care, application, and removal).

Material composition

Nylon, Polyurethane, Elastane

Disposal

Dispose of properly after use.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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