medi GmbH & Co.KG  
Medicusstraße 1  
95448 Bayreuth  
Germany  
T +49 921 912-0  
F +49 921 912-57  
medi@medi.de  
www.medi.de

medi Australia Pty Ltd  
Unit 4/13-21 Hallmark Street  
Pendle Hill NSW 2145  
Australia  
T +61 2 9890 8696  
F +61 2 9890 8439  
sales@medi.au  
www.mediaustralia.com.au

medi Austria GmbH  
Adamgasse 16/7  
6020 Innsbruck  
Austria  
T +43-512 57 95 15  
F +43-512 57 95 15 45  
vertrieb@medi-austria.at  
www.medi-austria.at

medi Belgium NV  
Staatsbaan 77/0099  
3945 Ham  
Belgium  
T: + 32-11 24 25 60  
F: +32-11 24 25 64  
info@medibelgium.be  
www.medibelgium.be

medi Canada Inc / médi Canada Inc  
104-1375 Lionel-Boulet,  
Varennes, Québec,  
QC Canada J3X 1P7  
T +1 450-583-3317 / +1 800 361 3153  
F +1 888-583-6827  
service@medicanada.ca  
www.medicanada.ca

UK Responsible Person  
medi UK Ltd.  
Plough Lane  
Hereford HR4 0EL  
Great Britain  
T +44-1432 37 35 00  
F +44-1432 37 35 10  
enquiries@mediuk.co.uk  
www.mediuk.co.uk

MAXIS a.s.,  
medi group company  
Slezská 2127/13  
120 00 Prague 2  
Czech Republic  
T: +420 571 633 510  
F: +420 571 616 271  
info@maxis-medica.com  
www.maxis-medica.com

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medi Norway AS  
Vestheimvegen 35  
4250 Kopervik  
Norway  
T: +47 52 84 45 00  
F: +47 52 84 45 09  
post@medinorway.no  
www.medinorway.no

medi Polska Sp. z o.o.  
ul. Łabędzka 22  
44-121 Gliwice  
Poland  
T: +48-32 230 60 21  
F: +48-32 202 87 56  
info@medi-polska.pl  
www.medi-polska.pl

medi Bayreuth Unipessoal Lda.  
Rua do Centro Cultural, no. 43  
1700-106 Lisbon  
Portugal  
T +351 21 843 71 60  
F +351 21 847 08 33  
medi.portugal@medibayreuth.pt  
www.medi.pt

medi RUS LLC  
Business Center NEO GEO  
Butlerova Street 17  
117342 Moscow  
Russia  
T +7-495 374 04 56  
F +7-495 374 04 56  
info@medirus.ru  
www.medirus.ru

medi Trading (Shanghai) Co. Ltd.  
Room 613, HnR Group Building  
No.1178, Tian Yao Qiao Road  
Xuhui District, Shanghai, 200232  
People's Republic of China  
T +86-21 50582318  
info@medi.cn  
www.medi.cn

medi Sweden AB  
Box 6034  
192 06 Sollentuna  
Sweden  
T +46 8 96 97 98  
info@medi.se  
www.medi.se

Medi Turk Ortopedi Medikal Ic ve Dis Tic. AS  
Mustafa Kemal Mah. 2156 Sokak No: 16 /4-5  
06510 Çankaya Ankara  
Turkey  
T: +90 312 435 20 26  
F: +90 312 434 17 67  
info@medi-turk.com  
www.medi-turk.com

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Tips on putting  
on and wearing  
compression  
stockings

Your medical compression stocking is a superior quality therapeutic garment. Here are a few application and care tips to help you:

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# sheer & soft®

### Intended purpose

Round knitted medical stockings are used for compression  of the lower extremities, mainly for the treatment of disorders of the venous systems.

### Features

Medical compression stocking works using a pressure gradient (decreasing from the end of your limb towards your heart) to exert compression on your limbs.

### Indications

Please always follow your doctor's orders and/or the orders of your medical specialist.

## Chronic venous diseases

- To improve venous symptoms
- To improve quality of life in chronic venous diseases
- Prevention and therapy of venous oedema
- Prevention and therapy of venous skin changes
- Eczema and pigmentation
- Lipodermatosclerosis and atrophic blanche
- Therapy of the venous ulcer
- Therapy for mixed (caused by arterial and venous issues) ulceration (taking the contraindications into account)
- Prevention of recurrent venous ulcer
- Pain reduction for venous ulcer
- Varicose veins
- Initial phase following varicose vein therapy
- Functional venous insufficiency (with obesity or jobs that involve a lot of sitting or standing)
- Venous malformations
- Severe chronic venous insufficiency

## Thromboembolic venous diseases

- Superficial venous thrombosis
- Deep vein thrombosis in the leg
- Conditions that follow on from thrombosis
- Post-thrombotic syndrome
- Thrombotic prophylaxis in mobile patients

## Oedema

- Lymphoedema
- Oedema during pregnancy
- Post-traumatic oedema
- Post-operative oedema
- Post-operative reperfusion oedema
- Cyclic idiopathic oedema
- Lipoedema
- Conditions involving stasis due to immobility (arthrogenic congestive syn-

drome, paresis and partial paresis of the extremity)

- Oedema caused by job (jobs which involve a lot of sitting or standing)
- Oedema caused by taking medication which cannot be altered (no viable alternative)

### Other indications

- Obesity with functional venous insufficiency
- Inflammatory dermatitis of the leg
- Nausea, dizziness during pregnancy
- Complaints involving stasis during pregnancy
- Conditions resulting from burns
- Scar treatment

As a rule, medi recommends round-knit compression stockings for venous conditions and flat-knit compression stockings for disorders of the lymphatic system. However, given certain conditions, a doctor may also decide that a flat-knit compressive aid is an appropriate therapy for a venous condition, for example (e.g. if the limb is very large at one point and narrow at another or if tissue folds are very deep). Patient-dependent factors in particular, such as the patient's weight, the type and severity of the oedema and the nature of their connective tissue all play a role in this decision.

### Contraindications

- Advanced peripheral arterial occlusive disease (if one of these parameters is present: ABPI < 0.5, ankle arterial pressure < 60 mmHg, toe pressure < 30mmHg or TcPO<sub>2</sub> < 20 mmHg on dorsum of foot). When non-elastic materials are implemented, the patient can still try out compression stocking if they have ankle arterial

pressure between 50 and 60 mmHg and they are under close clinical supervision.

- Decompensated heart failure (NYHA III + IV)
- Septic thrombophlebitis
- Phlegmasia cerulea dolens

In the following instances, the therapeutic decision should be made by weighing up the benefits and risks and by selecting the most suitable compressive equipment:

- Severe weeping dermatosis
- Intolerance to compressive material
- Severe paraesthesia in the limbs
- Advanced peripheral neuropathy (e.g. in diabetes mellitus)
- Primary chronic polyarthritis

### Risks/side effects

Especially if handled incorrectly, medical compression stockings may incur

- skin necrosis and
- pressure damage to peripheral nerves.
- If your skin is sensitive, the compressive equipment may cause itchiness, peeling or signs of inflammation. Therefore, applying adequate skin-care under the compression product is vital. medi offers a skincare range (medi day, medi night) that is specially formulated for compression stockings. Please visit a medical specialist retailer to discuss which product is right for you.

If the following symptoms are experienced then the compression clothing should be removed immediately and the clinical findings should be checked: If the toes are blue or white in colour, if you experience paraesthesia or numbness, increasing pain, shortness of breath and break out in sweats or acute limitation in terms of your movement.

## **Intended users and patient target groups**

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

## **Length of time that compression clothing should be worn and used for**

Unless otherwise specified by the doctor, wear your mediven compression stockings from morning to night every day. Compression stockings with a silicone top band can cause skin irritation in people with sensitive skin. Even if the maximum wearing time is exceeded by a significant amount in exceptional cases, the top band should be moved frequently if you are wearing a thigh-length stocking, or you should consider wearing tights as an alternative. Daily wearing and washing of your compression stockings can cause the compression and elasticity of your stockings, which are required to make them medically effective, to deteriorate over time. The stockings should therefore only be used for a maximum of six months. After this period, if a repeat prescription for medical compression stockings is issued, the patient's dimensions must be checked again by the medical specialist retailer. The date next to the hourglass symbol (label) contains the maximum use period of six months.

## **Important note**

The sewn-in fabric label must not be removed, as this renders the claim to a guarantee or exchange null and void.

## **Instructions for putting the stocking on**

Put the compression stocking on directly after you get up. Fabric or rubber gloves protect the knit fabric against damage and facilitate putting on and distributing the fabric of the compressive aid considerably. Ask about donning aids such as the medi Butler and medi 2in1.

- Reach into the compression stocking and grasp the heel. Hold the heel and turn the compression stocking inside out.
- The foot part (which has remained inside the stocking) now forms an opening. Stretch the opening out slightly using both hands.
- Slide your foot into the opening and pull the compression stocking over your foot carefully, up to the heel.
- Then distribute the material evenly over the leg without causing wrinkles. Do this by pulling it up the leg bit by bit, ensuring that the stocking is not overstretched. Check that it is correctly fitted on the toe and heel. If required, correct this once you have put it on by smoothing the stocking down slightly.
- Finally, check the position of your compression stocking. A knee stocking should be positioned around 1 cm under the back of the knee, a thigh stocking should be positioned around 2 cm under the fold of the buttock. Tights should be comfortable around the waist and the gusset should fit securely around the crotch.

## **Tips**

Grip the knit fabric carefully when put-

ting the stocking on – do not grip the seam of the top band.

### Instructions for care

Wash your stocking every day after you have worn it. We recommend that you use a net washing bag. Fabric conditioners, fats, oils, lotions, balms and soap residue can damage the material and lead to skin irritation. Wash it separately or with clothing that is similar in colour.

 Wash the product by hand, preferably with medi clean soap or in a delicates wash cycle at 40°C with mild detergent and without fabric conditioner.

 Do not bleach.

 Air dry or dry in a tumble dryer using the delicates cycle.

 Do not iron.

 Do not dry-clean.

### Liability

The manufacturer's liability will become void if the product is not used as intended (including own repairs or those of external service providers).

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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### Important notes

This medical device **MD** is made only for multiple uses on a single patient . If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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medi Ukraine LLC