medi GmbH & Co. KG
 Medicusstraße 1
 95448 Bayreuth
 Germany
 F +49 921 912-0
 F +49 921 912-57
 info@medi.de
 www.medi.de

medi Australia Pty Ltd
 83 Fennell Street
 North Parramatta NSW 2151
 Australia
 T +61 2 9890 8696
 F +61 2 9890 8439
 sales@mediaustralia.com.au
 www.mediaustralia.com.au

medi Austria GmbH
 Adamgasse 16/7
 6020 Innsbruck
 Austria
 T +43-512 57 95 15
 F +43-512 57 95 15 45
 vertrieb@medi-austria.at
 www.medi-austria.at

medi Belgium NV
 Staatsbaan 77/0099
 3945 Ham
 Belgium
 T: +32-11 24 25 60
 F: +32-11 24 25 64
 info@medibelgium.be
 www.medibelgium.be

medi Canada Inc / médi Canada Inc
 104-1375 Lionel-Boulet,
 Varennes, Québec,
 QC Canada J3X 1P7
 T +1 450-583-3317 / +1 800 361 3153
 F +1 888-583-6827
 service@medicanada.ca
 www.medicanada.ca

MAXIS a.s.,
 medi group company
 Slezská 2127/13
 120 00 Prague 2
 Czech Republic
 T: +420 571 633 510
 F: +420 571 616 271
 info@maxis-medica.com
 www.maxis-medica.com

medi Danmark ApS
 Vejleåvej 66
 2635 Ishøj
 Denmark
 T +45-70 25 56 10
 F +45-70 25 56 20
 kundeservice@sw.dk
 www.meditandmark.dk

UK Responsible Person
 medi UK Ltd.
 Plough Lane
 Hereford HR4 0EL
 Great Britain
 T +44-1432 37 35 00
 F +44-1432 37 35 10
 enquiries@mediuk.co.uk
 www.mediuk.co.uk

medi Bayreuth España SL
 C/Canigó 2 – 6 bajos
 Hospitalet de Llobregat
 08901 Barcelona
 Spain
 T +34 93 260 04 00
 F +34 93 260 23 14
 medi@mediespana.com
 www.mediespana.com

medi France
 Z.I. Charles de Gaulle
 25, rue Henri Farman
 93297 Tremblay en France Cedex
 France
 T +33 1 48 61 76 10
 F +33 1 49 63 33 05
 infos@medi-france.com
 www.medi-france.com

medi Hungary Kft.
 Bokor u. 21.
 1037 Budapest
 Hungary
 T +36 1371 0090
 F +36 1371 0091
 info@medi.hu
 www.medi.hu

medi Italia s.r.l.
 Via Giuseppe Ghedini, 2
 40069 Zola Predosa (BO)
 Italy
 T: +39-051 6 13 24 84
 F: +39-051 6 13 29 56
 info@medi-italia.it
 www.medi-italia.it

medi Japan K.K.
 5-1 Nihonbashi Kabutocho Chuoku
 103-0026 Tokyo
 Japan
 T: +81 3 6778 2590
 F: +81 3 5847 7901
 info@medi-japan.jp
 www.medi-japan.co.jp

Tachezy Sanit Srl.
 Via Druso 48/a
 39100 Bolzano
 Italy
 T: +39 0471 288 068
 F: +39 0471 282 773
 info@tachezysanit.com
 www.tachezysanit.com

medi

duomed®

Gebrauchsanweisung. Instructions for use.
 Mode d'emploi. Instrucciones de uso. Instruções
 para aplicação. Istruzioni per l'uso.
 Gebruiksaanwijzing. Brugsvejledning.
 Brugsanvisning. Návod k použití. Uputa za
 uporabu. Инструкција по испoлзуванoу.
 Kullanim kılavuzu. Instrukcja zakładowania. Οδηγία
 εφαρμογής. Használati útmutató. Інструкція з
 використання. Navodila za uporabo. Návod na
 použitie. دليل الاستخدام. 使用说明. Instrucțiuni
 de utilizare. הוראות שימוש. Bruksanvisning.
 Инструкции за употреба. Jalga panemise
 juhised. Käyttööhje. Uzvilksanas
 instrukcija. Užsimovimo instrukcija.
 使用説明書. Hướng dẫn cách sử dụng.
 यूजर मैनुअल.

Trocken lagern / Store in a dry place / Stocker au sec / Almacenar seco /
 Conserve em ambiente seco / Conservare in luogo fresco e asciutto / Droog
 bewaren / Opbevaars tørt / Lagras torr / Skladujte v suchu / Skladištiti suho /
 Хранить в сухом месте / Kuru bir yerde saklayın / Przechowywać w suchym
 miejscu / Να αποθηκεύεται σε ξηρό μέρος / Száraz helyen tárolandó /
 Зберігати сухими / Hraniti na suhem / Skladovať na suchom mieste /
 יישן במקום יבש / يخبزن في مكان جاف / 存放在干燥环境中 / Depozitați uscat /
 / Oppbevaars tørt / Съхранявайте в сухо състояние / Hoida kuivas kohas /
 Säilytettävä kuivassa / Uzglabāt sausā vietā / Laikyti sausoje vietoje / 乾燥存
 放 / Bảo quản khô ráo / सुखी जगह पर स्टोर करें

Vor Sonnenlicht schützen / Keep out of sunlight / Protéger du soleil / Proteger
 de la luz solar / Proteger da luz solar / Riparare dai raggi solari / Beschermen
 tegen zonlicht / Beskyttes mod sollys / Skyddas mot direkt solljus / Chraňte
 před slunečním zářením / Čuvati od sunčeve svjetlosti / Защищать от
 воздействия солнечного света / Güneş ışınlarından koruyun / Chronić przed
 słońcem / Να προστατεύεται από το ηλιακό φως / Napfénytől védendő /
 Защищати від сонячних променів / Zaščititi pred sončno svetlobo / Chrániť
 pred slnečným svetlom / لا تعرضه لأشعة الشمس / 避免阳光直射 / Protejați de
 lumina soarelui / יש להגן מפני אור שמש / Beskyttes mot solstråler / Pazere ot
 slънчева светлина / Kaitsta päikesevalguse eest / Säilytettävä suojassa
 auringonvalolta / Sargāt no saules gaismas / Saugoti nuo saulės šviesos / 避免
 直接日晒 / Tránh ánh nắng mặt trời / सूर्य की किरणों से सुरक्षित रखें

Bitte über den Hausmüll entsorgen / Please dispose of this product with your
 household waste / Élimination avec les déchets ménagers / Rogamos que se
 elimine en la basura doméstica / Eliminar no lixo doméstico / Smaltire il
 prodotto con i rifiuti domestici / Gelieve mee te geven met het huisvuil /
 Bortskaffes sammen med husholdningsaffaldet / Släng i hushållssoporna /
 Zlikvidujte prosím jako domovní odpad / Molimo da se odloži u kućansko
 smeće / Пожаљуйста, выбрасывайте с бытовыми отходами / Lűtfen evsel
 atiklar úzerinden bertaraf edin / Wyrzucić z odpadami gospodarstwa
 domowego / Να απορρίπτεται στα οικιακά απορρίμματα / A háztartási
 hulladékkal együtt ártalmatlanítsa / Будь ласка, утилізуйте разом із
 побутовими відходами / Odložite med gospodinjske odpadke / Likvidujte,
 prosím, s domovým odpadom / يرجى التخلص من المنتج مع النفايات المنزلية / 请作为生
 活垃圾废弃处理 / xVã rugãm sã eliminați cu deșeurile menajere /
 הביתית לאשפה המוצר את השליך / Kan kastes i husholdningsavfallet /
 Моля, изхвърлете заедно с битовите отпадъци / Kõrvaldage koos
 olmejäätmetega / Hävitä kotitalousjätteen mukana / Likvidējiet, izmetot pie
 sadzīves atkritumiem / Šalinkite su buitinėmis atliekomis / 請與家庭垃圾一
 起處理 / Xin vút bỏ cùng với rác thải sinh hoạt trong nhà / कृपया घरेलू
 कचरे में निपटान करें

- duomed
- duomed smooth


E016183/01.2024


medi. I feel better.

1

4

2

5

3

DEUTSCH / GERMAN

duomed®

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medi Norway AS
Vestheimvegen 35
4250 Kopervik
Norway
T: +47 52 84 45 00
F: +47 52 84 45 09
post@medinorway.no
www.medinorway.no

medi Polska Sp. z o.o.
ul. Łabędzka 22
44-121 Gliwice
Poland
T: +48-32 230 60 21
F: +48-32 202 87 56
info@medi-polska.pl
www.medi-polska.pl

medi Bayreuth Unipessoal Lda.
Rua do Centro Cultural, no. 43
1700-106 Lisbon
Portugal
T +351 21 843 71 60
F +351 21 847 08 33
medi.portugal@medibayreuth.pt
www.medi.pt

medi RUS LLC
Business Center NEO GEO
Butlerova Street 17
117342 Moscow
Russia
T +7-495 374 04 56
F +7-495 374 04 56
info@medirus.ru
www.medirus.ru

medi Trading (Shanghai) Co. Ltd.
Room 1703, T1, CES WEST BUND CENTER,
No. 277, Long Lan Road
Xuhui District
200232 Shanghai
People's Republic of China
T: +86-21 50582319
F: +86-21 50582319

medi Sweden AB
Box 6034
192 06 Sollentuna
Sweden
T +46 8 96 97 98
F +46 8 626 68 70
info@medi.se
www.medi.se

Medi Turk Ortopedi Medikal Iç ve Dis Tic. AS
Mustafa Kemal Mah. 2156 Sokak No: 16/4-5
06510 Çankaya Ankara
Turkey
T: +90 312 435 20 26
F: +90 312 434 17 67
info@medi-turk.com
www.medi-turk.com

medi Ukraine LLC
Evhena Sverstiuka str, 11.
Kiev 02002
Ukraine
T: +380 44 591 11 63
F: +380 44 392 73 73
info@medi.ua
http://medi.ua

medi USA L.P.
6481 Franz Warner Parkway
Whitsett, N.C. 27377-3000
USA
T +1-336 4 49 44 40
F +1-888 5 70 45 54
info@mediusa.com
www.mediusa.com

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www.medi.biz/anziehen

Tips on putting on and wearing compression stockings

Your medical compression stocking is a superior quality therapeutic garment. Here are a few application and care tips to help you:

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duomed®

Intended purpose

Round knitted medical compression stockings are used for compression of the lower extremities, mainly for the treatment of disorders of the venous systems.

Features

Medical compression stocking works using a pressure gradient (decreasing from the end of your limb towards your heart) to exert compression on your limbs.

Indications

Chronic venous diseases

- To improve venous symptoms
- To improve quality of life in chronic venous diseases
- Prevention and therapy of venous oedema
- Prevention and therapy of venous skin changes
- Eczema and pigmentation
- Lipodermatosclerosis and atrophie blanche
- Therapy of the venous ulcer
- Therapy for mixed (caused by arterial and venous issues) ulceration (taking the contraindications into account)
- Prevention of recurrent venous ulcer
- Pain reduction for venous ulcer
- Varicose veins
- Initial phase following varicose vein therapy
- Functional venous insufficiency (with obesity or jobs that involve a lot of sitting or standing)
- Venous malformations
- Severe chronic venous insufficiency

Thromboembolic venous diseases

- Superficial venous thrombosis
- Deep vein thrombosis in the leg
- Conditions that follow on from thrombosis
- Post-thrombotic syndrome
- Thrombotic prophylaxis in mobile patients

Oedema

- Lymphoedema
- Oedema during pregnancy
- Post-traumatic oedema
- Post-operative oedema
- Post-operative reperfusion oedema
- Cyclic idiopathic oedema
- Lipoedema
- Conditions involving stasis due to immobility (arthrogenic congestive syndrome, paresis and partial paresis of the extremity)
- Oedema caused by job (jobs which involve a lot of sitting or standing)
- Oedema caused by taking medication which cannot be altered (no viable alternative)

Other indications

- Obesity with functional venous insufficiency
- Inflammatory dermatitis of the leg
- Nausea, dizziness during pregnancy
- Complaints involving stasis during pregnancy
- Conditions resulting from burns
- Scar treatment

As a rule, medi recommends round-knit compression stockings for venous conditions and flat-knit compression stockings for disorders of the lymphatic system. However, given certain conditions, a doctor may also decide that a flat-knit compressive aid is an appropriate therapy for a venous condition, for example (e.g. if the limb is very large at one point and narrow at another or if tissue folds are very deep). Patient-dependent factors

in particular, such as the patient's weight, the type and severity of the oedema and the nature of their connective tissue all play a role in this decision.

Contraindications

- Advanced peripheral arterial occlusive disease (if one of these parameters is present: ABPI < 0.5, ankle arterial pressure < 60 mmHg, toe pressure < 30mmHg or TcPO₂ < 20 mmHg on dorsum of foot). When non-elastic materials are implemented, the patient can still try out compression stocking if they have ankle arterial pressure between 50 and 60 mmHg and they are under close clinical supervision.
 - Decompensated heart failure (NYHA III + IV)
 - Septic thrombophlebitis
 - Phlegmasia cerulea dolens
- In the following instances, the therapeutic decision should be made by weighing up the benefits and risks and by selecting the most suitable compressive equipment:
- Severe weeping dermatosis
 - Intolerance to compressive material
 - Severe paraesthesia in the limbs
 - Advanced peripheral neuropathy (e.g. in diabetes mellitus)
 - Primary chronic polyarthritis

Risks/side effects

Especially if handled incorrectly, medical compression stockings may incur

- skin necrosis and
- pressure damage to peripheral nerves.
- If your skin is sensitive, the compressive equipment may cause itchiness, peeling or signs of inflammation. Therefore, applying adequate skincare under the compression product is vital. medi offers a skin-care range (medi day, medi night) that is specially formulated for compression stockings. Please visit a medical specialist retailer to discuss which product is right for you.

If the following symptoms are experienced then the compression clothing should be removed immediately and the clinical findings should be checked:

Blue or white discolouration of toes, paraesthesia and numbness, increasing pain, shortness of breath and sudden sweating, acute limitation of mobility, or swelling of the foot.

Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals. Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Length of time that compression clothing should be worn and used for

Unless otherwise specified by the doctor, wear your medi ven compression stockings from morning to night every day. Compression stockings with a silicone top band can cause skin irritation in people with sen-

sitive skin. Even if the maximum wearing time is exceeded by a significant amount in exceptional cases, the top band should be moved frequently if you are wearing a thigh-length stocking, or you should consider wearing tights as an alternative. Daily wearing and washing of your compression stockings can cause the compression and elasticity of your stockings, which are required to make them medically effective, to deteriorate over time. The stockings should therefore only be used for a maximum of six months. After this period, if a repeat prescription for medical compression stockings is issued, the patient's dimensions must be checked again by the medical specialist retailer. The date next to the hourglass symbol (label) contains the maximum use period of six months.

Important note

The sewn-in fabric label must not be removed, as this renders the claim to a guarantee or exchange null and void.

Instructions for putting the stocking on

Put the compression stocking on directly after you get up. Fabric or rubber gloves protect the knit fabric against damage and facilitate putting on and distributing the fabric of the compressive aid considerably. Ask about donning aids such as the medi Butler and medi 2in1.

For compression stockings with open ends, first place your foot in the attached donning aid.

- Reach into the compression stocking and grasp the heel. Hold the heel and turn the compression stocking inside out.
- The foot part (which has remained inside the stocking) now forms an opening. Stretch the opening out slightly using both hands.
- Slide your foot into the opening (if the stocking is open-ended at the toe, do this using the donning aid) and pull the compression stocking over your foot carefully, up to the heel.
- Then distribute the material evenly over the leg without causing wrinkles. Do this by pulling it up the leg bit by bit, ensuring that the stocking is not over-stretched. Check that it is correctly fitted on the toe and heel. If required, correct this once you have put it on by smoothing the stocking down slightly.

For compression stockings with open ends, remove the donning aid after this step by folding the foot part of the stocking up towards the heel and pulling the donning aid out. After this, pull the foot component over your foot again and spread the knitted fabric evenly up to the start of the toes.

- Finally, check the position of your compression stocking. A knee stocking should be positioned around 1 cm under the back of the knee, a thigh stocking should be positioned around 2 cm under the fold of the buttock. Tights should be comfortable around the waist and the gusset should fit securely around the crotch.

For stockings without a foot part (IBD): Tightly hold the opening with your hands at the textile tag of the top band and step through both openings in the weave, foot first, from above. Position the lower, narrower band of the compressive weave approx. 1 cm above the ankle. Pull up the top band until approx. 1 cm below the back

of the knee. Make sure the stocking fits without creases and smooth out any creases that might be present in the fabric.

Tips

Grip the knit fabric carefully when putting the stocking on – do not grip the seam of the top band.

Instructions for care

Wash your stocking every day after you have worn it. We recommend that you use a net washing bag. Fabric conditioners, fats, oils, lotions, balms and soap residue can damage the material and lead to skin irritation. Wash it separately or with clothing that is similar in colour.

 Wash the product by hand, preferably with mild clean soap or in a delicates wash cycle at 40°C with mild detergent and without fabric conditioner.

 Do not bleach.

 Air dry or dry in a tumble dryer using the delicates cycle.

 Do not iron.

 Do not dry-clean.

Liability

The manufacturer's liability will become void if the product is not used as intended (including own repairs or those of external service providers).

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If undue pain or an unpleasant sensation develops while you are wearing the product, please consult your doctor or orthotist immediately. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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