medi GmbH & Co.KG  
 Medicusstraße 1  
 95448 Bayreuth  
 Germany  
 T +49 921 912-0  
 F +49 921 912-57  
 medi@medi.de  
 www.medi.de

medi Australia Pty Ltd  
 Unit 4/13-21 Hallmark Street  
 Pendle Hill NSW 2145  
 Australia  
 T +61 2 9890 8696  
 F +61 2 9890 8439  
 sales@medi.au  
 www.mediaustralia.com.au

medi Austria GmbH  
 Adamgasse 16/7  
 6020 Innsbruck  
 Austria  
 T +43-512 57 95 15  
 F +43-512 57 95 15 45  
 verttrieb@medi-austria.at  
 www.medi-austria.at

medi Belgium NV  
 Staatsbaan 77/0099  
 3945 Ham  
 Belgium  
 T: +32-11 24 25 60  
 F: +32-11 24 25 64  
 info@medibelgium.be  
 www.medibelgium.be

medi Canada Inc / médi Canada Inc  
 104-1375 Lionel-Boulet,  
 Varennes, Québec,  
 QC Canada J3X 1P7  
 T +1 450-583-3317 / +1 800 361 3153  
 F +1 888-583-6827  
 service@medicanada.ca  
 www.medicanada.ca

UK Responsible Person  
 medi UK Ltd.  
 Plough Lane  
 Hereford HR4 OEL  
 Great Britain  
 T +44-1432 37 35 00  
 F +44-1432 37 35 10  
 enquiries@mediuk.co.uk  
 www.mediuk.co.uk

MAXIS a.s.,  
 medi group company  
 Slezská 2127/13  
 120 00 Prague 2  
 Czech Republic  
 T: +420 571 633 510  
 F: +420 571 616 271  
 info@maxis-medica.com  
 www.maxis-medica.com

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## medi travel® women / medi travel® men Instructions for Use

### 1. Intended use

Round-knit compression  stocking for the lower limbs which provides physical prophylaxis of thrombosis to ensure that your veins remain healthy while travelling.

### 2. Features of the product

medi travel® women / medi travel® men stocking works using a graduated pressure pattern (decreasing from the end of your limb towards your heart) to exert compression on your limbs.

### 3. Indications

Please always follow your doctor's orders and/or the orders of your medical specialist. Your medi travel women/men can be used for physical thrombosis prophylaxis of the lower limbs during travel, as long as the user has no other conditions affecting the venous or lymphatic systems or other conditions that require compressive therapy to be applied under the responsibility and supervision of medical specialists.

### 4. Contraindications

Medical compression stocking may not be worn if the patient has the following conditions:

- Advanced peripheral arterial occlusive disease (if one of these parameters is present: ABPI < 0.5, ankle arterial pressure < 60 mmHg, toe pressure < 30mmHg or TcPO2 < 20 mmHg on dorsum of foot). When non-elastic materials are implemented, the patient can still try out compression clothing if they have ankle

arterial pressure between 50 and 60 mmHg and they are under close clinical supervision.

- Decompensated heart failure (NYHA III + IV)
- Septic thrombophlebitis
- Phlegmasia cerulea dolens

In the following instances, the therapeutic decision should be made by weighing up the benefits and risks and by selecting the most suitable compressive equipment:

- Severe weeping dermatosis
- Intolerance to compressive material
- Severe paraesthesia in the limbs
- Advanced peripheral neuropathy (e.g. in diabetes mellitus)
- Primary chronic polyarthritis

If you are not sure whether one or more of these statements applies to you, speak to your doctor and/or the medical personnel treating you.

The following risks and the following side effects may occur:

Medical compression stocking may incur the following, especially if handled incorrectly:

- Skin necrosis and
- Pressure damage to peripheral nerves.

If your skin is sensitive, the compressive equipment may cause itchiness, peeling or signs of inflammation. Adequate skin-care under the compression stocking is thus vital. Therefore, please take particular care to observe our important notes and the fitting instructions for the compression stocking (sections 7 and 8).

If the following symptoms are experienced then the compression stocking should be removed immediately and the clinical findings should be checked:

If the toes are blue or white in colour, if

you experience paraesthesia or numbness, increasing pain, shortness of breath and break out in sweats or acute limitation in terms of your movement.

## **5. Intended users and patient target groups**

The intended users include healthcare professionals and people who are both planning to travel and have healthy veins.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

## **6. Length of time that the product should be worn and used for**

For best possible results, the travel stocking should be worn throughout the entire duration of the travel period (i.e. during the trip there and the return trip). A precisely graduated compression pattern (high pressure at bottom, decreasing towards the top) is required to ensure that these stockings work optimally. Wearing and washing your compression stockings frequently can cause the compression and elasticity of your stockings, which are required to make them medically effective, to deteriorate over time.

## **7. Instructions for use**

- The compression stocking should be issued by trained specialists with a briefing on how to handle it and how to care for it (see also sections 8 and 10 on this).
- If pain or increased irritation of the skin occurs while you are wearing the stocking, remove it immediately and consult your doctor or the medical specialists.
- Appropriate skincare must be applied when wearing the travel stockings. Creams

or salves that are particularly greasy and also soap residue can cause skin irritation and material wear and can thus impair the efficacy of the stocking. Therefore, medi also offers a skincare range (medi day, medi night, medi soft foam lotion) that is specially formulated for stockings. Please visit a medical specialist retailer to discuss which product is right for you.

- Sharp fingernails or toenails, faults in your shoe's lining, jewellery or rough or calloused heels may cause damage to the knit fabric.
- Do not attempt to repair your compression stocking yourself or to have it repaired by an external service provider, as this will render the guarantee void, without exception. Furthermore, this can impair the quality, safety and efficacy of the product. You can obtain information on correct repair by medi – if possible – from your medical specialist retailer or our services.
- In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).
- The sewn-in fabric label is an important component in ensuring that the identification and traceability of the product is guaranteed. Therefore, ensure that it is not removed from the stocking under any circumstances. Do not tear it out or cut it out, as this renders any claim to the guarantee, repair or exchange null and void.

## 8. Fitting instructions

Explanatory videos with tips and tricks for putting on the clothing correctly can be found here:

<https://vimeo.com/377253363> (Putting on the AD with closed toe: instructive video)

Ideally, you should put the stockings on before you begin your trip. If you have issues putting them on, ask for medi fabric or rubber gloves in order to make it easier to distribute the fabric over the leg, or ask for a handy donning aid such as medi Butler or medi 2in1.

When you are putting the stocking on, ensure that your fingernails have been filed into a smooth, rounded shape. It is also important to ensure that your heels and toenails are well cared for to avoid causing damage to the knit fabric. Only put your jewellery or watch on once you have put the compression stocking on. To make it easy to put your stockings on, sit down and place your foot on an elevated support, such as a small stool.

**A.** Reach into the stocking and grasp the heel. Hold the heel and turn the stocking inside out.

The foot part (which has remained inside the stocking) now forms an opening. Stretch the opening out slightly using both hands.

Slide your foot into the opening (if the stocking has an opening in the toe, do this using the donning aid) and pull the stocking over your foot carefully, up to the heel.

**B.** Then distribute the material evenly over the leg without causing wrinkles. Do this by pulling it up the leg bit by bit, ensuring that the stocking is not overstretched. In the interim, check that it is correctly fitted on the toe and heel. If

required, correct this once you have put it on by pulling the stocking down slightly (tip: this is easier if you use the gloves).

C. Finally, check that your stocking is sitting correctly. medi travel® women / medi travel® men should sit around 2 cm under the crease of the knee.

**Tip for taking off the clothing:** Fold the stocking down in stages until you reach the heel. When you reach the heel, grip the knit fabric and pull it over the foot.

## 9. Storage instructions

The products must be stored in a dry place and out of direct sunlight.

- The hourglass symbol shown on the box label shows the latest date that the compression clothing can be used up until. For medi travel® women / medi travel® men, this is 60 months in total from the production date.

## 10. Care instructions

- Wash your stocking every day after you have worn it. We recommend that you use medi clean or a specially formulated hand wash detergent to ensure gentle care. Alternatively, you can use washing detergent for delicates to wash the clothing by hand or in the machine (delicates cycle at 40°C) without optical bleach or conditioner.

**Tip:** A net washing bag offers extra protection for the knit fabric.

- We recommend that you turn your stocking inside out when washing it. Wash it separately or with clothing that is similar in colour.
- If the compression stocking is stained, we recommend that you use medi spot ex. Please do not dry clean medi travel® women / medi travel® men or clean it with household agents.

- After washing, you can roll the compression stocking up in a towel and press it to dry it. Do not leave the damp stocking in the wet towel, instead let it dry in the air. Do not place it in sunlight or on a radiator/heating appliance.
- Alternatively, you could also dry your medi travel® women / medi travel® men in your tumble dryer using the delicates cycle.

- 
- Wash up to 40°C (delicates cycle)
  - Do not bleach
  - Tumble dry at low temperatures (delicates cycle)
  - Do not iron
  - Do not dry-clean
- 

## 11. Material composition

- All medi compression stocking contains polyamide and elastane. You can find more detailed information on the fabric label sewn into your compression stocking.

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medi Ukraine LLC