**medi GmbH & Co. KG**  
**Medicusstraße 1**  
**95448 Bayreuth**  
**Germany**  
**T +49 921 912-0**  
**F +49 921 912-780**  
**ortho@medi.de**  
**www.medi.de**

**medi Australia Pty Ltd**  
**83 Fennell Street**  
**North Parramatta NSW 2151**  
**Australia**  
**T +61-2 9890 8696**  
**F +61-2 9890 8439**  
**sales@mediaustralia.com.au**  
**www.mediaustralia.com.au**

**medi Austria GmbH**  
**Adamgasse 16/7**  
**6020 Innsbruck**  
**Austria**  
**T +43 512 57 95 15**  
**F +43 512 57 95 15 45**  
**vertrieb@medi-austria.at**  
**www.medi-austria.at**

**medi Belgium NV**  
**Staatsbaan 77/0099**  
**3945 Ham**  
**Belgium**  
**T: + 32-11 24 25 60**  
**F: +32-11 24 25 64**  
**info@medibelgium.be**  
**www.medibelgium.be**

**medi Canada Inc / médi Canada Inc**  
**104-1375 Lionel-Boulet,**  
**Varenes, Québec,**  
**QC Canada J3X 1P7**  
**T +1 450-583-3317 / +1 800 361 3153**  
**F +1 888-583-6827**  
**service@medicanada.ca**  
**www.medicana.ca**

**MAXIS a.s.,**  
**medi group company**  
**Slezská 2127/13**  
**120 00 Prague 2**  
**Czech Republic**  
**T: +420 571 633 510**  
**F: +420 571 616 271**  
**info@maxis-medica.com**  
**www.maxis-medica.com**

**medi Danmark ApS**  
**Vejlegardsvej 59**  
**2665 Vallensbaek Strand**  
**Denmark**  
**T +45-70 25 56 10**  
**kundeservice@sw.dk**  
**www.meditandmark.dk**

**medi Bayreuth Espana SL**  
**C/Canigo 2-6 bajos**  
**Hospitalet de Llobregat**  
**08901 Barcelona**  
**Spain**  
**T +34-932 60 04 00**  
**F +34-932 60 23 14**  
**medi@mediespana.com**  
**www.mediespana.com**

**medi France**  
**Z.I. Charles de Gaulle**  
**25, rue Henri Farman**  
**93297 Tremblay en France Cedex**  
**France**  
**T +33-1 48 61 76 10**  
**F +33-1 49 63 33 05**  
**infos@medi-france.com**  
**www.medi-france.com**

**medi Hungary Kft.**  
**Bokor u. 21.**  
**1037 Budapest**  
**Hungary**  
**T +36 1 371-0090**  
**F +36 1 371-0091**  
**info@medi.hu**  
**www.medi.hu**

**medi Japan K.K.**  
**5-1 Nihonbashi Kabutocho Chuoku**  
**103-0026 Tokyo**  
**Japan**  
**T: +81 3 6778 2590**  
**F: +81 3 5847 7901**  
**info@medi-japan.jp**  
**www.medi-japan.co.jp**

**medi Nederland BV**  
**Heusing 5**  
**4817 ZB Breda**  
**The Netherlands**  
**T +31-76 57 22 555**  
**F +31-76 57 22 565**  
**info@medi.nl**  
**www.medi.nl**

**medi Polska Sp. z o.o.**  
**ul.Łabędzka 22**  
**44-121 Gliwice**  
**Poland**  
**T: +48-32 230 60 21**  
**F: +48-32 202 87 56**  
**info@medi-polska.pl**  
**www.medi-polska.pl**

**medi Bayreuth Unipessoal, Lda**  
**Rua do Centro Cultural, no. 43**  
**1700-106 Lisbon**  
**Portugal**  
**T +351-21 843 71 60**  
**F +351-21 847 08 33**  
**medi.portugal@medibayreuth.pt**  
**www.medi.pt**

**medi RUS LLC**  
**Business Center NEO GEO**  
**Butlerova Street 17**  
**117342 Moscow**  
**Russia**  
**T +7-495 374 04 56**  
**F +7-495 374 04 56**  
**info@medirus.ru**  
**www.medirus.ru**

**medi Sweden AB**  
**Box 6034**  
**192 06 Sollentuna**  
**Sweden**  
**T +46 8 96 97 98**  
**F +46 8 626 68 70**  
**info@medi.se**  
**www.medi.se**

**medi Trading (Shanghai) Co. Ltd.**  
**Room 1703, T1, CES WEST BUND CENTER,**  
**No. 277, Long Lan Road**  
**Xuhui District**  
**200232 Shanghai**  
**People's Republic of China**  
**T: +86-21 50582319**  
**F: +86-21 50582319**

**UK Responsible Person**  
**medi UK Ltd.**  
**Plough Lane**  
**Hereford HR4 OEL**  
**Great Britain**  
**T +44-1432 37 35 00**  
**F +44-1432 37 35 10**  
**enquiries@mediuk.co.uk**  
**www.mediuk.co.uk**

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### Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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# Collamed®

## Intended purpose

Collamed is a knee support brace which functions based on the 4-point principle using limitation of extension / flexion .

## Indications

All indications that require reinforcement of the physiological stabilisation of the knee joint and/or stress relief of the joint, such as:

- Following injury to the collateral ligament (post-traumatic/post-operative)
- Following injury to the cruciate ligament (post-traumatic/post-operative)
- For knee instability
- Following meniscal injury
- For knee osteoarthritis

## Contraindications

None known at present.

## Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in

conjunction with sweating) or due to the composition of the material. The immobilisation which is achievable by the device presents a risk for development of thrombosis in the lower veins of the legs, especially in the context of surgical interventions or trauma.

## Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

## Setting the flexion / extension limits Only to be carried out by an orthopedic technician.

Extension limit wedge, EXT (limits extension): -5°, 0°, 10°, 20°, 30°

Flexion limit wedge, FLEX (limits flexion): 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°

Can be immobilised at:  
0°, 10°, 20°, 30°

The medial (inside) and the lateral (outside) wedge settings must always be the same.

1. Only change the wedges if the orthoses are in the lateral supports.
2. 0° extension wedges and 90° flexion wedges are pre-installed in your knee brace.
3. To change the wedges, open the top and bottom latches. Now, simply remove the wedges.

4. In the wedge set, the replacement wedges are marked R – FLEX, L – FLEX or R – EXT, L – EXT and with the respective angles.
5. Insert the new flexion wedge first and the new extension wedge next into the joint. Then close the two latches to secure the wedges. Repeat the process with the other joint.  
**WARNING:** The extension and flexion wedges always have to be inserted.
6. Move the brace joints to check if the wedges are securely seated.

### Important notes

Only change the extension and flexion limits on instruction from your doctor or orthotist. In order to prevent hyperextension, always insert extension stops or 0° stops. The stops define the exact excursion of the orthosis. In difficult soft-tissue conditions, e.g. where the soft-tissues are under compression, it may be necessary to insert a larger stop to limit the joint movement to the desired excursion. In addition, make sure that the orthosis hinges are correctly positioned: the middles of the hinges should be in one line with the upper border of the patella. The hinge must lie behind the midline of the leg.

### Fitting instructions

- Undo all straps and the two hook and loop strips. Pull the straps back toward the back so that they will not interfere with you putting on the brace (Image 1).
- Now put your leg into a 45-degree angle flexion position. Wrap the brace tightly around your leg from the back and close it with the two hook and loop strips (Image 2).
- Make sure that the kneecap is positioned centrally in the space for the patella. Position the brace so that the middles of the hinges are level with the upper border of the kneecap and

- just behind the lateral centre line of the leg (Image 3).
- First fasten the lower calf strap in the front (Image 4) and then the lower calf strap in the back (Image 5).
- Now fasten the upper thigh strap and remove the small synthetic plates (Images 6+7).
- Now fasten all remaining straps, following the number order.
- If possible, walk a few steps to check that the brace is in the correct position. Please adjust the position as required (Image 8).

### Care instructions

Please close the hook and loop fasteners before washing and remove the orthoses (Images 9,10,11). Soap residues, lotions and ointments can cause skin irritation and material wear.

- Wash the product by hand, preferably using medi clean washing agent.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

With normal use, the hinge mechanism should not require any maintenance, but it can be lubricated with PTFE spray.

### Storage instructions

Keep the brace in a dry place and do not expose to direct sunlight.

### Material composition

Avional, polyamide, elasthane, polyester

### Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and

instructions in this manual.

**Disposal**

The product can be disposed of in the domestic waste.

**Your medi team**

wishes you a quick return to full fitness

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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