medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth / Germany
T +49 921 912-0
F +49 921 912-780
ortho@medi.de
www.medi.de

medi Australia Pty Ltd
83 Fennell Street
North Parramatta NSW 2151
Australia
T +61 2 9890 8696
F +61 2 9890 8439
sales@mediaustralia.com.au
www.mediaustralia.com.au

medi Austria GmbH
Adamgasse 16/7
6020 Innsbruck
Austria
T +43-512 57 95 15
F +43-512 57 95 15 45
vertrieb@medi-austria.at
www.medi-austria.at

medi Belgium NV
Staatsbaan 77/0099
3945 Ham
Belgium
T: + 32-11 24 25 60
F: +32-11 24 25 64
info@medibelgium.be
www.medibelgium.be

medi Canada Inc / médi Canada Inc
104-1375 Lionel-Boulet,
Varennes, Québec,
QC Canada J3X 1P7
T +1 450-583-3317 / +1 800 361 3153
F +1 888-583-6827
service@medicanada.ca
www.medicanada.ca

MAXIS a.s.,
medi group company
Slezská 2127/13
120 00 Prague 2
Czech Republic
T: +420 571 633 510
F: +420 571 616 271
info@maxis-medica.com
www.maxis-medica.com

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medi Norway AS
Vestheimvegen 35
4250 Kopervik
Norway
T: +47 52 84 45 00
F: +47 52 84 45 09
post@medinorway.no
www.medinorway.no

medi Polska Sp. z o.o.
ul.Łabędzka 22
44-121 Gliwice
Poland
T: +48-32 230 60 21
F: +48-32 202 87 56
info@medi-polska.pl
www.medi-polska.pl

medi Bayreuth Unipessoal Lda.
Rua do Centro Cultural, no. 43
1700-106 Lisbon
Portugal
T +351 21 843 71 60
F +351 21 847 08 33
medi.portugal@medibayreuth.pt
www.medi.pt

medi RUS LLC
Business Center NEO GEO
Butlerova Street 17
117342 Moscow
Russia
T +7-495 374 04 56
F +7-495 374 04 56
info@medirus.ru
www.medirus.ru

UK Responsible Person
medi UK Ltd.
Plough Lane
Hereford HR4 0EL
Great Britain
T +44-1432 37 35 00
F +44-1432 37 35 10
enquiries@mediuk.co.uk
www.mediuk.co.uk

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Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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Collamed® long

Intended purpose

Collamed is a knee support brace which functions based on the 4-point principle using limitation of extension / flexion .

Indications

All indications that require reinforcement of the physiological stabilisation of the knee joint and/or stress relief of the joint, such as:

- Following injury to the collateral ligament (post-traumatic/postoperative)
- Following injury to the cruciate ligament (post-traumatic/postoperative)
- For knee instability
- Following meniscal injury
- For knee osteoarthritis

Contraindications

None known at present.

Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in

conjunction with sweating) or due to the composition of the material. The immobilisation which is achievable by the device presents a risk for development of thrombosis in the lower veins of the legs, especially in the context of surgical interventions or trauma.

Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Setting the flexion / extension limits Only to be carried out by an orthopedic technician.

Extension limit wedge, EXT (limits extension): -5°, 0°, 10°, 20°, 30°

Flexion limit wedge, FLEX (limits flexion): 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°

Can be immobilised at:
0°, 10°, 20°, 30°

The medial (inside) and the lateral (outside) wedge settings must always be the same.

1. Your knee brace is equipped with 0° extension wedges and 90° flexion wedges.
2. To change the wedges, open the top and bottom sliding lock. Now, the wedges can simply be taken out.
3. The replacement wedges are included in the set of wedges labelled with R

- FLEX, L - FLEX or R - EXT, L - EXT and the respective degree settings.
- 4. Firstly, insert the new flexion wedge and then the extension wedge into the joint. Afterwards, close the two sliders to secure the wedges. Repeat the same process at the other joint.
WARNING: Extension and flexion wedges must always be inserted together.
- 5. Move the brace joints to check that the wedges are positioned securely.

Important notes

Only change the extension and flexion limits on instruction from your doctor or orthotist.

In order to prevent hyperextension, always insert extension stops or 0° stops. The stops define the exact excursion of the orthosis. In difficult soft-tissue conditions, e.g. where the soft-tissues are under compression, it may be necessary to insert a larger stop to limit the joint movement to the desired excursion. In addition, make sure that the orthosis hinges are correctly positioned: the middles of the hinges should be in one line with the upper border of the patella. The hinge must lie behind the midline of the leg.

Fitting instructions

- Open all straps and the two front hook and loop fasteners. Release the straps from one of the two hook and loop points (Image 1).
- Put your leg into a 45-degree angle flexion position. Wrap the brace closely around the leg from behind and close the two hook and loop fasteners (Image 2).
- Make sure that the kneecap is positioned centrally in the space for the patella. Position the brace in such a way that the centres of the hinges are level with the upper edge of the kneecap and just behind the middle of

- the side of the leg. (Image 3).
- Now, fasten the bottom calf strap (Image 4) and then the top thigh strap (Image 5).
- Now fasten all remaining straps, following the number order.
- If possible, walk a few steps in order to check whether the brace is in the correct position (Image 6). Please adjust the position as required (Image 7).

Care instructions

Please close the hook and loop fasteners before washing and remove the orthoses. Soap residues, lotions and ointments can cause skin irritation and material wear.

- Wash the product by hand, preferably using medi clean washing agent.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

Storage instructions

Keep the brace in a dry place and do not expose to direct sunlight.

Material composition

Avional, polyamide, elasthane

Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

Disposal

The product can be disposed of in the domestic waste.

Your medi team

wishes you a quick return to full fitness

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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