**medi GmbH & Co. KG**  
**Medicusstraße 1**  
**95448 Bayreuth**  
**Germany**  
**T +49 921 912-0**  
**F +49 921 912-780**  
**ortho@medi.de**  
**www.medi.de**

medi Australia Pty Ltd  
 Unit 4/13-21 Hallmark Street  
 Pendle Hill NSW 2145  
 Australia  
 T +61 2 9890 8696  
 F +61 2 9890 8439  
 sales@medi.au  
 www.medi.australia.com.au

medi Austria GmbH  
 Adamgasse 16/7  
 6020 Innsbruck  
 Austria  
 T +43 512 57 95 15  
 F +43 512 57 95 15 45  
 vertrieb@medi-austria.at  
 www.medi-austria.at

medi Belgium NV  
 Staatsbaan 77/0099  
 3945 Ham  
 Belgium  
 T: + 32-11 24 25 60  
 F: +32-11 24 25 64  
 info@medibelgium.be  
 www.medibelgium.be

medi Canada Inc / médi Canada Inc  
 104-1375 Lionel-Boulet,  
 Varennes, Québec,  
 QC Canada J3X 1P7  
 T +1 450-583-3317 / +1 800 361 3153  
 F +1 888-583-6827  
 service@medicanada.ca  
 www.medicanada.ca

MAXIS a.s.,  
 medi group company  
 Slezská 2127/13  
 120 00 Prague 2  
 Czech Republic  
 T: +420 571 633 510  
 F: +420 571 616 271  
 info@maxis-medica.com  
 www.maxis-medica.com

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medi Trading (Shanghai) Co. Ltd.  
 Room 1703, T1, CES WEST BUND CENTER,  
 No. 277, Long Lan Road  
 Xuhui District  
 200232 Shanghai  
 People's Republic of China  
 T: +86-21 50582319  
 F: +86-21 50582319

Medi Turk Ortopedi Medikal Iç ve Dis Tic. AS  
 Mustafa Kemal Mah. 2156 Sokak No: 16 /4-5  
 06510 Çankaya Ankara  
 Turkey  
 T: +90 312 435 20 26  
 F: +90 312 434 17 67  
 info@medi-turk.com  
 www.medi-turk.com

UK Responsible Person  
 medi UK Ltd.  
 Plough Lane  
 Hereford HR4 0EL  
 Great Britain  
 T +44-1432 37 35 00  
 F +44-1432 37 35 10  
 enquiries@mediuk.co.uk  
 www.mediuk.co.uk

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### Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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Aluminium, Polyamid, Elasthan, Polyurethan

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# medi Epico ROMs PRU

## Intended purpose

medi Epico ROMs PRU is an elbow brace for mobilisation  with adjustable immobilisation  of the proximal radioulnar articulation.

## Indications

All indications that require early functional mobilisation of the elbow joint with limitation of the range of motion and adjustable immobilisation of the proximal radioulnar articulation, such as:

- Surgical fixation of medial and/or lateral collateral ligaments
- Surgically treated bicep or tricep tears
- Conservative therapy of dislocated elbows
- Surgery on dislocated elbow fractures
- Medial or lateral collateral ligament reconstructive surgery
- Treatment of radial head fracture (nonsurgical, postoperative)
- Treatment following insertion of a prothesis
- Distal humeral fractures (nonsurgical, postoperative)

## Contraindications

- Concomitant neurological damage, e.g. ulnar nerve
- Compartment syndrome

## Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with

infectious indications (excessive heat, swelling or reddening)

- Paresthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in conjunction with sweating) or due to the composition of the material.

## Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

## QuickLock mechanism

The QuickLock mechanism allows a quick fixation or release of the range of motion (ROM). In order to fixate the joint, move the bar at the hinge to the closed symbol (). To release the pre-set ROM, move the bar to the open symbol () again.

## Setting the level of flexion and extension limitation

- The medi Epico ROMs proximal radioulnar joint brace comes with a pre-installed extension and flexion angle of 0° extension and 90° flexion. In order to adjust the flexion or

extension, first release the QuickLock mechanism (Fig. 1). Then, turn the indicator until the symbol  aligns with the pin for limiting flexion.

- Note: The points of the indicator and the flexion pin must be exactly opposite each other (Fig. 2).
- Now you can pull out the light grey pin and adjust the flexion (Fig. 3). To do this, move the flexion pin onto the desired number of degrees (Fig. 4).
- Then push the pin towards the centre of the joint again to secure the setting (Fig. 5). In order to adjust the extension, turn the indicator to the light grey pin for limiting the extension and proceed as you would for adjusting the amount of flexion.
- To ensure that the angle of the articulated joint of the brace does not change unless desired (e.g., by being bumped), make sure that the indicator is at the closed symbol  again (Fig. 6).
- Limitation of flexion:  $-10^{\circ}$ - $120^{\circ}$
- Limitation of extension:  $-10^{\circ}$ - $90^{\circ}$

### Adjusting the hand support

The hand support of the brace is fixed in a neutral position. In order to adjust the proximal radioulnar joint, loosen the two hexagon socket screws (Fig. 7). Now the joint is freely mobile. Depending on the specification, you can make the desired adjustments to pronation/supination by screwing the screws into the provided positions 0-3 ( $1=30^{\circ}$ ;  $2=60^{\circ}$ ;  $3=90^{\circ}$ ).

### Instructions for putting the brace on

- Open all of the hook-and-loop straps and the shoulder strap and put on the brace so that the articulated joint of the brace is positioned exactly over the lateral epicondyle (Fig. 8).
- The length of the hand support can be adjusted using the two screws on the lower arm splint (Fig. 9).

- First close the two hook-and-loop fasteners on the straps on the lower arm and tighten the straps (Fig. 10). Then proceed in the same manner for the straps on the upper arm. If necessary, the padding along the raised seams can be cut to size using scissors (Fig. 11). The straps can also be individually shortened. To do this, unfasten the hook-and-loop fastener at the end of the strap and shorten it accordingly (Fig. 12).
- Lastly, place the shoulder strap over your neck. The neck pad can be moved to suit the individual.
- If you fasten the straps too tightly, this can cause lymphostasis. The straps should therefore be opened frequently.

---

### Care instructions

Hook and loop fasteners should be closed for washing. Fats, oils, lotions, ointments and soap residues can affect the material and cause skin irritations.

- Wash the paddings and straps by hand, preferably using medi clean washing agent.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

With normal use, the hinge mechanism should not require any maintenance, but it can be lubricated with PTFE spray.

---

### Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

**Material composition**

Aluminum, polyamide, elastane, polyurethane

**Liability**

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

**Disposal**

The product can be disposed of in the domestic waste.

**Your medi team**

Wishes you a speedy recovery!

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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