**medi GmbH & Co. KG**  
 Medicusstraße 1  
 95448 Bayreuth  
 Germany  
 T +49 921 912-0  
 F +49 921 912-780  
 ortho@medi.de  
 www.medi.de

medi Australia Pty Ltd  
 83 Fennell Street  
 North Parramatta NSW 2151  
 Australia  
 T +61-2 9890 8696  
 F +61-2 9890 8439  
 sales@mediaustralia.com.au  
 www.mediaustralia.com.au

medi Austria GmbH  
 Adamgasse 16/7  
 6020 Innsbruck  
 Austria  
 T +43 512 57 95 15  
 F +43 512 57 95 15 45  
 vertrieb@medi-austria.at  
 www.medi-austria.at

medi Belgium NV  
 Staatsbaan 77/0099  
 3945 Ham  
 Belgium  
 T: + 32-11 24 25 60  
 F: +32-11 24 25 64  
 info@medibelgium.be  
 www.medibelgium.be

medi Canada Inc / médi Canada Inc  
 104-1375 Lionel-Boulet,  
 Varennes, Québec,  
 QC Canada J3X 1P7  
 T +1 450-583-3317 / +1 800 361 3153  
 F +1 888-583-6827  
 service@medicanada.ca  
 www.medicanada.ca

MAXIS a.s.,  
 medi group company  
 Slezská 2127/13  
 120 00 Prague 2  
 Czech Republic  
 T: +420 571 633 510  
 F: +420 571 616 271  
 info@maxis-medica.com  
 www.maxis-medica.com

medi Danmark ApS  
 Vejlegardsvej 59  
 2665 Vallensbaek Strand  
 Denmark  
 T +45-70 25 56 10  
 kundeservice@sw.dk  
 www.medidanmark.dk

medi Bayreuth Espana SL  
 C/Canigo 2-6 bajos  
 Hospitalet de Llobregat  
 08901 Barcelona  
 Spain  
 T +34-932 60 04 00  
 F +34-932 60 23 14  
 medi@mediespana.com  
 www.mediespana.com

medi France  
 Z.I. Charles de Gaulle  
 25, rue Henri Farman  
 93297 Tremblay en France Cedex  
 France  
 T +33-1 48 61 76 10  
 F +33-1 49 63 33 05  
 infos@medi-france.com  
 www.medi-france.com

medi Hungary Kft.  
 Bokor u. 21.  
 1037 Budapest  
 Hungary  
 T +36 1 371-0090  
 F +36 1 371-0091  
 info@medi.hu  
 www.medi.hu

medi Japan K.K.  
 5-1 Nihonbashi Kabutocho Chuoku  
 103-0026 Tokyo  
 Japan  
 T: +81 3 6778 2590  
 F: +81 3 5847 7901  
 info@medi-japan.jp  
 www.medi-japan.co.jp

medi Nederland BV  
 Heusing 5  
 4817 ZB Breda  
 The Netherlands  
 T +31-76 57 22 555  
 F +31-76 57 22 565  
 info@medi.nl  
 www.medi.nl

medi Polska Sp. z o.o.  
 ul. Łabędzka 22  
 44-121 Gliwice  
 Poland  
 T: +48-32 230 60 21  
 F: +48-32 202 87 56  
 info@medi-polska.pl  
 www.medi-polska.pl

medi Bayreuth Unipessoal, Lda  
 Rua do Centro Cultural, no. 43  
 1700-106 Lisbon  
 Portugal  
 T +351-21 843 71 60  
 F +351-21 847 08 33  
 medi.portugal@medibayreuth.pt  
 www.medi.pt

medi RUS LLC  
 Business Center NEO GEO  
 Butlerova Street 17  
 117342 Moscow  
 Russia  
 T +7-495 374 04 56  
 F +7-495 374 04 56  
 info@medirus.ru  
 www.medirus.ru

UK Responsible Person  
 medi UK Ltd.  
 Plough Lane  
 Hereford HR4 0EL  
 Great Britain  
 T +44-1432 37 35 00  
 F +44-1432 37 35 10  
 enquiries@mediuk.co.uk  
 www.mediuk.co.uk

## medi Rehab one

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#### Important notes

This medical device **MD** is made only for multiple uses on a single patient . If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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# medi Rehab one

## Intended purpose

medi Rehab one is a round-knit medical compression  stocking for compressive care of the lower limbs, mainly for prevention and treatment of post-operative and post-traumatic oedema and general swelling.

## Indications

### Oedema

- Post-traumatic oedema
- Post-operative oedema

### Contraindications

- Advanced peripheral arterial occlusive disease (if one of these parameters is present: ABPI < 0.5, ankle arterial pressure < 60 mmHg, toe pressure < 30 mmHg or TcPO<sub>2</sub> < 20 mmHg on dorsum of foot).
- Decompensated heart failure (NYHA III + IV)
- Septic phlebitis
- Phlegmasia cerulea dolens

### Please note:

In the following instances, the therapeutic decision should be made by a doctor, who weighs up the benefits and risks and selects the most suitable compressive equipment:

- Severe weeping dermatosis
- Patient unable to tolerate compressive stocking material
- Severe paraesthesia in the limb
- Advanced peripheral neuropathy (e.g. diabetes mellitus)
- Rheumatoid arthritis

## Risks and side effects

Medical compression stockings may cause skin necrosis and pressure damage to the peripheral nerves, especially if it is handled incorrectly. If your skin is sensitive, medical compression stockings may cause itching, peeling or signs of inflammation. Therefore, adequate skincare applied under the compression stocking is a must.

## Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

## Note on wearing

Unless otherwise prescribed, you should wear medi Rehab one from morning to evening. If this duration is significantly exceeded in exceptional cases, the top band should be repositioned several times within this period.

Daily wearing and washing of your compression stockings can cause the compression and elasticity of your stockings, which are required to make them medically effective, to deteriorate over time. The stockings should

therefore only be used for a maximum of six months.

### Important note

- If pain or increased irritation of the skin occurs while you are wearing the stocking, remove it immediately and consult your doctor or the medical specialists.
- Appropriate skincare must be maintained during compression therapy. This is essential. Creams or ointments that are particularly greasy and also soap residue can cause skin irritation and material wear and can thus impair the efficacy of the compression stocking. Therefore, medi also offers a skincare range (medi day, medi night, medi soft foam lotion) that is specially formulated for compression stockings. Please visit a medical specialist retailer to discuss which product is right for you.
- Sharp fingernails or toenails, faults in your shoe's lining, jewellery or rough or calloused heels may cause damage to the knit fabric.
- The sewn-in fabric label is an important component in ensuring that the identification and traceability of the product is guaranteed. Please do not remove it from the stocking by detaching it or cutting it out.

### Fitting instructions

- Put the compression stocking on immediately after you get up. If you have a shower in the morning and then apply skin care products, then you should wait a while until your skin is completely dry. This will make it far easier to put on the compression stocking.

- Place your foot in the medi donning aid included with this product. (Fig. 1)
- Reach into the medi Rehab one, take the heel and turn it inside out. (Figs. 2 + 3)
- Now pull the stocking over the donning aid and over your heel. (Fig. 4)
- Grip the top layer of the stocking as shown in Fig. 5 and pull it onto your lower leg. (Fig. 6)
- Repeat the process step by step until the stocking is fully on. The green border should come to just under the knee. (Figs. 7 + 8)
- To remove the donning aid, pull the end of the stocking towards your heel and back over the heel. Ensure that you do not create any creases while doing this. Creases could result in painful constriction of your leg. (Fig. 9)
- The medi donning aid can now easily be pulled out. (Fig. 10)
- Finally, pull the forefoot part of the stocking back over the heel and foot. Pull it so that it is evenly distributed up to the root of the toes. Ensure that the heel is sitting correctly. (Fig. 11)
- Check the entire stocking to ensure that it is seated correctly.

### Please note the following

The stocking must lie correctly and should not be twisted (the ribbed pattern provides orientation for this) or stretched too long.

---

### Care instructions

Wash your stocking every day after you have worn it. Ideally, turn the stocking inside out and wash it separately or with items of clothing that are the same colour. We recommend using a laundry bag. Soap residue, creams or

ointments can lead to irritation of the skin and deterioration of the material. Please only use medi spot ex to remove stains.

- Wash the product by hand, preferably with medi clean soap, or in a delicates wash cycle at 40°C with mild detergent and without fabric conditioner.
- Do not bleach.
- Dry at a low temperature.
- Do not iron.
- Do not dry-clean.

- After washing, you can roll the compression stocking up in a towel and squeeze it to dry it. Do not leave the damp stocking in the wet towel, instead let it dry in the air.
- Do not place it in sunlight or on a radiator/heating appliance.

### Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

### Material composition

Polyamide, elastane, polyester, Merino wool

### Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

### Disposal

The product can be disposed of in the domestic waste.

### Your medi team

Wishes you a speedy recovery!

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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