medi GmbH & Co. KG  
Medicusstraße 1  
95448 Bayreuth / Germany  
T +49 921 912-0  
F +49 921 912-780  
ortho@medi.de  
www.medi.de

medi Australia Pty Ltd  
83 Fennell Street  
North Parramatta NSW 2151  
Australia  
T +61-2 9890 8696  
F +61-2 9890 8439  
sales@mediaustralia.com.au  
www.mediaustralia.com.au

medi Austria GmbH  
Adamgasse 16/7  
6020 Innsbruck  
Austria  
T +43 512 57 95 15  
F +43 512 57 95 15 45  
vertrieb@medi-austria.at  
www.medi-austria.at

medi Belgium NV  
Staatsbaan 77/0099  
3945 Ham  
Belgium  
T: +32-11 24 25 60  
F: +32-11 24 25 64  
info@medibelgium.be  
www.medibelgium.be

medi Canada Inc / médi Canada Inc  
104-1375 Lionel-Boulet,  
Varenes, Québec,  
QC Canada J3X 1P7  
T +1 450-583-3317 / +1 800 361 3153  
F +1 888-583-6827  
service@medicanada.ca  
www.medicanada.ca

MAXIS a.s.,  
medi group company  
Slezská 2127/13  
120 00 Prague 2  
Czech Republic  
T: +420 571 633 510  
F: +420 571 616 271  
info@maxis-medica.com  
www.maxis-medica.com

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#### Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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# medi Ankle RTS

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medi Ankle RTS is an adjustable ankle brace for stabilisation  in one plane.

---

## Indications

For all indications which require immobilisation of the upper or lower ankle in one plane, such as:

- After sprains (e.g. supination injuries)
- Chronic ligament instability in the upper and lower ankle
- As a preventative measure after supination injuries, particularly for sports
- Mild injury to the syndesmotic ligaments

## Contraindications

None known at present.

## Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paresthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in conjunction with sweating) or due to the composition of the material.

## Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

## Note on wearing

Supports are particularly beneficial during physical exercise. The supports are designed to fit comfortably and may be worn all day. It should be removed during prolonged periods of inactivity (e.g. sitting for a long time or a long car drive, sleeping).

The product should not be worn on bare skin, but instead only over a sock

or compression stocking. For greater stability during activities, we recommend wearing the brace in a shoe.

### Putting on / taking off

Please note that initial fitting should be supervised by a qualified healthcare professional.

- Undo all the straps of the brace so that you can place your foot into it easily (Fig. 1+2). If needed, you can pull the narrower strap out of the eyelet. To do this, remove the Y hook-and-loop fastening (Fig. 3). We particularly recommend doing this in the first few days, during which painful stretching of the ankle should be avoided. This allows you to place your foot into the brace more easily. This step is not necessary in the later course of treatment.
- Now put your foot into the brace and place your heel in the intended heel opening. Check that your heel is as far back in the brace as possible (Fig. 4). For the next steps, ensure that your foot is in a neutral position (at 90° to your leg) (Fig. 5).
- Grip the longer, wide strap and pull it over the arch of your foot (instep) from the inside of your foot to the outside. Then fasten it to the black fuzzy patch on the outside (Fig. 6). Continue by wrapping the strap around the back of your calf. While wrapping, ensure that the strap does not cover the grey fuzzy patch on the inside of the brace (Fig. 7). Fasten the end of the strap to the intended grey fuzzy patch on the outside (Fig. 8). The strap should neither be too loose nor pulled too tight and should make you feel secure.
- Pull the narrow strap until it is snug, redirect it and fasten it to the intended light grey fuzzy patch at the upper edge of the brace (Fig. 9+10). If this strap is too long, it can be shortened. To do this, remove the Y hook from the end of the strap and shorten the strap using scissors (Fig. 11). This adjustment can be repeated once the swelling has gone down.
- Check that the brace fits snugly and that you feel stable in it. If necessary, correct the fit.

---

### Care instructions

Hook and loop fasteners should be closed for washing. We recommend using a laundry bag. Fats, oils, lotions, ointments and soap residues can affect the material and cause skin irritations.

- Wash the product by hand, preferably using a medi clean detergent, or in delicate cycle at 30°C using a mild detergent without fabric conditioners.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

**Storage instructions**

Keep the product in a dry place and do not expose to direct sunlight.

---

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PA, PU, EVA, polyester

**Liability**

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

**Disposal**

Dispose of properly after use.

**Your medi team**

Wishes you a speedy recovery!

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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