**medi GmbH & Co. KG**  
 Medicusstraße 1  
 95448 Bayreuth  
 Germany  
 F +49 921 912-0  
 T +49 921 912-57  
 info@medi.de  
 www.medi.de

medi Austria GmbH  
 Adamgasse 16/7  
 6020 Innsbruck  
 Austria  
 P +43 512 57 95 15  
 F +43 512 57 95 15 45  
 vertrieb@medi-austria.at  
 www.medi-austria.at

medi Belgium NV  
 Staatsbaan 77/0099  
 3945 Ham  
 Belgium  
 T +32 11 24 25 60  
 F +32 11 24 25 64  
 info@medibelgium.be  
 www.medibelgium.be

medi Canada Inc / médi Canada Inc  
 104-1375 Lionel-Boulet,  
 Varennes, Québec,  
 QC Canada J3X 1P7  
 T +1 450-583-3317 / +1 800 361 3153  
 F +1 888-583-6827  
 service@medicanada.ca  
 www.medicanada.ca

MAXIS a.s.,  
 medi group company  
 Slezská 2127/13  
 120 00 Prague 2  
 Czech Republic  
 T: +420 571 633 510  
 F: +420 571 616 271  
 info@maxis-medica.com  
 www.maxis-medica.com

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# mediven®

### Intended use

Round- or flat-knitted medical compression garment used for compression of the upper extremities, mainly for the treatment of disorders of the venous or lymphatic systems.

### Features

Medical compression garments work using a graduated pressure pattern (decreasing from the end of your limb towards your heart) to exert compression on your limbs.

### Indications

#### Chronic venous conditions

- Venous malformations

#### Thromboembolic venous diseases

- Superficial vein thrombosis
- Venous thrombosis of the upper extremity
- A post-thrombotic condition
- Post-thrombotic syndrome
- Thrombotic prophylaxis in mobile patients

#### Oedema

- Lymphoedema
- Post-traumatic oedema
- Post-operative oedema
- Post-operative reperfusion oedema
- Cyclic idiopathic oedema
- Lipoedema
- Conditions involving stasis due to immobility

(paresis and partial paresis of the extremity)

- Oedema caused by taking medication which cannot be altered (no viable alternative)

### Other indications

- Conditions following burn injuries
- Scar treatment

As a rule, medi recommends round-knit compression stockings for venous conditions and flat-knit compression stockings for disorders of the lymphatic system.

However, given certain conditions, a doctor may also decide that a flat-knit compressive aid is an appropriate therapy for a venous condition, for example (e.g. if the limb is very large at one point and narrow at another or if tissue folds are very deep). Patient-dependent factors in particular, such as the patient's weight, the type and severity of the oedema and the nature of their connective tissue all play a role in this decision.

### Contraindications

Medical compression garments may not be worn if the patient has the following conditions:

- Advanced peripheral arterial occlusive disease (if one of these parameters is present: ABPI < 0.5, ankle arterial pressure < 60 mmHg, toe pressure < 30mmHg or TcPO<sub>2</sub> < 20 mmHg on dorsum of foot). When non-elastic materials are implemented, the patient can still try out compression stocking if they have ankle arterial pressure between 50 and 60 mmHg and they are under close clinical supervision.
- Decompensated heart failure (NYHA III + IV)
- Septic thrombophlebitis
- Phlegmasia cerulea dolens

In the following instances, the therapeutic decision should be made by weighing up the benefits and risks and by selecting the most suitable compressive equipment:

- Severe weeping dermatosis
- Intolerance to compressive material
- Severe paraesthesia in the limbs
- Advanced peripheral neuropathy (e.g. in diabetes mellitus)
- Primary chronic polyarthritis

Medical compression garments may incur the following, especially if handled incorrectly:

- Skin necrosis and
- Pressure damage to peripheral nerves.

If your skin is sensitive, the compressive equipment may cause itchiness, peeling or signs of in-

flammation. Adequate skincare under the compression clothing is vital. Therefore, please take particular care to observe our important notes and the instructions for putting the compression garments on.

If the following symptoms are experienced then the compression garments should be removed immediately and clinical advice sought: If the toes are blue or white in colour, if you experience paraesthesia or numbness, increasing pain, shortness of breath and break out in sweats or acute limitation in terms of your movement

#### **Intended users and patient target groups**

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

#### **Length of time that compression clothing should be worn and used for**

Unless otherwise specified by the doctor, wear your mediven compression stocking from morning to night for optimal therapeutic results. A precisely graduated compression pattern (high pressure at bottom, decreasing towards the top) is required to ensure that compression garments function optimally. Daily wearing and washing of the garments can cause the compression and elasticity of your garments, which are required to make them medically effective, to deteriorate over time. The stockings should therefore only be used for a maximum of six months. If the medical compression garments has been used for six months and a repeat prescription is issued, the patient's dimensions must be checked again by the medical specialist retailer as they could have changed due to the patient's condition or their individual circumstances. Observe our instructions for putting on the armsleeve and instructions for care to ensure that the compressive effect is maintained as well as possible.

#### **Important note**

The sewn-in fabric label must not be removed, as

this renders the claim to a guarantee or exchange null and void.

#### **Instructions for putting the stocking on**

Explanatory videos with tips and tricks for putting on the compression garments correctly can be found here:

Ideally, you should put the compression garments on immediately after you wake up. If you have a shower in the morning and then apply skin care products, then you should wait a while until your skin is completely dry. This will make it far easier to put on the compression garments. Please also ensure that you have your medi fabric gloves or rubber gloves ready. These are available from specialist medical retailers. You will use them later. They make it easier to distribute the fabric over the arm, or you can use them to put the sleeve on, right from the beginning. The gloves proffer excellent purchase and make it easier to put on the sleeve.

**Tip:** Ask about our handy donning aid, the medi Arm Butler.

When you are putting the sleeve on, ensure that your fingernails have been filed to a smooth, rounded shape so that the knit fabric does not become damaged. Only put your jewellery or watch on once you have put the compression garments on. Instructions for putting the compressive garments on using the medi Arm Butler

#### **Instructions for PUTTING the compressive garments on using the medi Arm Butler**

**Tip:** Place a non-slip surface under the medi Butler (e.g. a sponge cloth or a rubberised sheet).

**A.** Place the medi Arm Butler on a table or similar surface.

**B.** Feed the compression stocking onto the semi-circular cylinder. If the compression stocking has a seam, ensure that the seam is not central and pointing directly to you but is instead slightly offset.

- Compression sleeves for the right arm - slightly to the outer right (at approx. 4 o'clock)
- Compression sleeves for the left arm - slightly to the outer left (at approx. 8 o'clock)

**C.** Now lie the medi Arm Butler flat on the table. To stabilise it, place the front bracket on the front edge of the table.

**D.** Put your hand through the opening. Ensure that you are holding down the Butler firmly

and carefully by its cuff with your other hand.

E. Insert your arm into the compression sleeve, keeping the angle the same and without turning it.

F. Check that the compression sleeve is seated correctly again. If you wish, fold up the top band or fastening.

### Instructions for care

Wash your stocking every day after you have worn it. We recommend that you use a net washing bag. Fabric conditioners, fats, oils, lotions, balms and soap residue can damage the material and lead to skin irritation. Wash it separately or with clothing that is similar in colour.

 Wash the product by hand, preferably with medi clean soap or in a delicate wash cycle at 40°C with mild detergent and without fabric conditioner.

 Do not bleach.

 Air dry or dry in a tumble dryer using the delicate cycle.

 Do not iron.

 Do not dry-clean.

### Liability

The manufacturer's liability will become void if the product is not used as intended (including own repairs or those of external service providers). In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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### **Important notes**

This medical device **MD** is made only for multiple uses on a single patient . If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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